


May 10, 2019

U Mass - Hyperinsulinemic-euglycemic clamp

DOI

dx.doi.org/10.17504/protocols.io.xzefp3e



Jason Kim¹

¹University of Massachusetts

Mouse Metabolic Phenotyping Centers

Tech. support email: info@mmpc.org



Lili Liang

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.xzefp3e

External link: <https://mmpc.org/shared/document.aspx?id=136&docType=Protocol>

Protocol Citation: Jason Kim 2019. U Mass - Hyperinsulinemic-euglycemic clamp. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xzefp3e>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 11, 2019

Last Modified: May 10, 2019

Protocol Integer ID: 20230

Keywords: Hyperinsulinemic-euglycemic clamp, type 2 diabetes,

Abstract

Summary:

Hyperinsulinemic-euglycemic clamp is the gold-standard method to assess insulin sensitivity. The hyperinsulinemic-euglycemic clamp is widely used in clinics and laboratories to measure insulin action on glucose utilization in humans and animals for clinical and basic science research. Incorporation of radioactive-labeled glucose during hyperinsulinemic-euglycemic clamps makes it possible to measure glucose metabolism in individual organs in awake mice. Impaired insulin sensitivity (insulin resistance) is a major characteristic of obesity and an early requisite event in the development of type 2 diabetes.

Materials

MATERIALS

⊗ HelixMark Standard Silicone Tubing **Helix Medical, Inc. Catalog #0.012" ID / 0.025" OD**

⊗ [3-3H] D-glucose **Perkin Elmer Catalog #NET331C005MC**

⊗ 2-[1-14C] Deoxy-D-glucos **Perkin Elmer Catalog #NEC495001MC**

⊗ Pentobarbital **Oak Pharmaceuticals, Inc. Catalog #NDC76478-501-50**

⊗ Microdialysis pumps **CMA/Microdialysis Catalog #CMA 402**

⊗ Analox GM7 Micro-stat Rapid Multi-assay Analyser **Analox Catalog #GM7**

⊗ Insulin **Novolin Catalog #Regular human insulin, U-100**

⊗ 20 % Dextrose injection USP **Pfizer (Hospira) Catalog #NDC0409-7935-19**

⊗ 0.9 % Sodium Chloride Injection USP **B.Braun Medical Inc Catalog #NDC0264-4001-55**

⊗ 1 ml tuberculin syringes **Becton Dickinson (BD) Catalog #REF 309659**

⊗ Microhematocrit capillary tubes **Fisher Scientific Catalog #22-362-566**

⊗ Heparin-coated blue polyethylene open-top tubes **Beckman Coulter Catalog #652825**

⊗ Microcentrifuge tubes (1.5 ml) **Denville Scientific Inc. Catalog #C-2170**

Note:

Hospira, [RRID:SCR_003985](#)

Fisher Scientific, [RRID:SCR_008452](#)

BD Biosciences, [RRID:SCR_013311](#)

Beckman Coulter, [RRID:SCR_008940](#)



- 1 Survival surgery is performed to establish a chronic indwelling catheter at 5~6 days prior to experiment for intravenous infusion. (refer to M1023: Surgery-jugular vein cannulation)
- 2 Mice are fasted overnight (~15 hours) or for 5 hours prior to the start of experiment.
- 3 Place a mouse in a rat-size restrainer with its tail tape-tethered at one end.
- 4 Expose and flush the intravenous catheter using saline solution. Then, connect the catheter to the CMA Microdialysis infusion pump.
- 5 During the 2-hour acclimation period, infuse D-[3-³H] glucose at 0.05 μ Ci/min to measure the basal rate of whole body glucose turnover.
- 6 Collect a plasma sample (30 μ l) at the end for the measurement of plasma glucose, insulin, and [³H] glucose concentrations (basal parameters).
- 7 Following the basal period, start a 2-hour hyperinsulinemic-euglycemic clamp with a primed (150 mU/kg body weight) and continuous infusion of human insulin at 2.5 mU/kg/min to raise plasma insulin levels.
- 8 Infuse 20% dextrose at variable rates to maintain plasma glucose at basal concentrations (euglycemia) throughout the 2-hour clamp.
- 9 Insulin-stimulated whole body glucose turnover rates are estimated with a continuous infusion of [3-³H] glucose at 0.1 μ Ci/min throughout the clamp.
- 10 Collect plasma samples (10 μ l) at 20, 40, 60, 70, 90, 100, 110, and 120 min to measure plasma glucose concentrations.
- 11 Adjust glucose infusion rates based on the instantaneous glucose levels to maintain euglycemia.
- 12 To estimate insulin-stimulated glucose uptake in individual organs, administer a bolus injection of 10 μ Ci of 2-deoxy-D-[1-¹⁴C] glucose (2-[¹⁴C]DG) at 75 minutes after the start of the clamp.
- 13 Collect plasma samples (10 μ l) at 80, 85, 90, 100, 110, and 120 min for the measurement of plasma [³H] glucose, ³H₂O, and 2-[¹⁴C] DG concentrations. (10 μ l plasma samples are

- suspended in 20 μl distilled water [dH_2O] to make 30 μl sample solutions.)
- 14 Collect additional plasma sample (10 μl) at the end of the clamp (at 120 min) to measure plasma insulin concentrations (clamp parameter).
 - 15 At the end of hyperinsulinemic-euglycemic clamp, anesthetize mice using pentobarbital and quickly dissect and collect tissues including skeletal muscles (gastrocnemius and quadriceps) from both hindlimbs, white and brown adipose tissues, liver, and heart.
 - 16 Rapidly freeze-clamp the tissues in liquid N_2 , and store tissue samples in -80°C freezer for biochemical analysis.
 - 17 Biochemical assay is conducted using plasma samples to measure [$3\text{-}^3\text{H}$] D-glucose, $^3\text{H}_2\text{O}$, and 2- ^{14}C] DG concentrations.
 - a) Transfer 15 μl of plasma sample solutions into microcentrifuge tubes with sample time clearly labeled.
 - b) Add 25 μl BaOH and vortex samples.
 - c) Add 25 μl $\text{Zn}(\text{SO})_2$ and vortex samples.
 - d) Centrifuge samples for 5 min at 12,000g (~14,000 rpm).
 - e) Prepare 2 sets of scintillation vials labeled Dry and Non-Dry for each sample.

Non-Dry samples

- i. Prepare 60 μl of dH_2O in NON-DRY labeled scintillation vials for each sample.
- ii. Transfer 20 μl of supernatant from step (d) into respective scintillation vials and vortex samples.
- iii. Add 3 ml of Ultima scintillation cocktail and vortex thoroughly.
- iv. Measure radioactive labeling using Beckman Coulter Scintillation Counter.

Dry samples

- i. Transfer 20 μl of supernatant from step (d) into respective scintillation vials and place into vacuum oven set at room temperature for overnight drying.
 - ii. Following overnight drying, add 80 μl dH_2O and vortex thoroughly.
 - iii. Add 3 ml of Ultima scintillation cocktail and vortex samples.
 - iv. Measure radioactive labeling using Beckman Coulter Scintillation Counter.
- 18 Plasma concentrations of $^3\text{H}_2\text{O}$ will be calculated as the difference in ^3H counts between Dry and Non-Dry samples and will be used to calculate the rate of whole body glycolysis.
 - 19 For biochemical assay to measure glucose uptake in individual organs, refer to M1003: Organ-specific glucose uptake experiment.

- 20 Basal rate of hepatic glucose production (HGP) or glucose turnover can be determined as the ratio of the basal [³H] glucose infusion rate (dpm/min) to the specific activity of plasma glucose (dpm/μmol) at the end of the basal period (0 min sample before the start of hyperinsulinemic-euglycemic clamp).
- 21 Insulin-stimulated whole body glucose turnover is determined as the ratio of the clamp [³H] glucose infusion rate (dpm/min) to the specific activity of plasma glucose (dpm/μmol) during the final 30 min of the clamp (90~120 min of clamp).
- 22 Insulin-stimulated HGP (during the clamp) is determined by subtracting the glucose infusion rate from the whole body glucose turnover rate. The difference between insulin-stimulated and basal rates of HGP reflects hepatic insulin action (insulin-mediated suppression of HGP).
- 23 Whole body glycolysis is calculated from the rate of increase in plasma ³H₂O concentrations, determined by linear regression of the measurements at 80, 85, 90, 100, 110, and 120 min of clamp.
- 24 Whole body glycogen plus lipid synthesis are estimated by subtracting whole body glycolysis from whole body glucose turnover.