

May 10, 2019

U Mass - Urea/BUN

DOI

dx.doi.org/10.17504/protocols.io.x5qfq5w



Jason Kim¹

¹University of Massachusetts

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x5qfq5w>

External link: <https://mmpc.org/shared/document.aspx?id=187&docType=Protocol>

Protocol Citation: Jason Kim 2019. U Mass - Urea/BUN. protocols.io <https://dx.doi.org/10.17504/protocols.io.x5qfq5w>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 14, 2019

Last Modified: May 10, 2019

Protocol Integer ID: 20368

Keywords: Serum Urea/BUN levels, kidney failure, renal complications of diabetes, bun level, bun summary, serum urea, using roche cobas clinical chemistry analyzer, roche cobas clinical chemistry analyzer, alterations in systemic protein, renal complication, nitrogen metabolism, systemic protein, urea, spectrophotometric measurement, diabetes, serum

Abstract

Summary:

This experiment involves a spectrophotometric measurement using Roche Cobas Clinical Chemistry Analyzer. Serum Urea/BUN levels are affected by alterations in systemic protein and nitrogen metabolism. Serum Urea/BUN levels are altered in kidney failure and renal complications of diabetes.

Materials

MATERIALS

- ⊗ Urea/BUN Roche Catalog #04657616 190
- ⊗ Calibrator f.a.s. Roche Catalog #10759350 360
- ⊗ Precinorm U plus Roche Catalog # 12149435 160
- ⊗ Precipath U plus Roche Catalog #12149443 160
- ⊗ NaCl Diluent 9 % Roche Catalog #04774230 190
- ⊗ Chimneys Roche Catalog #11930630 001
- ⊗ Cleaner Roche Catalog #04774248 190
- ⊗ Micro Sample cups Roche Catalog #11406680 001
- ⊗ NERL High Quality Water Fisher Scientific Catalog #9805

Note:

Roche, [RRID:SCR_001326](#)

Fisher Scientific, [RRID:SCR_008452](#)

Troubleshooting



Before start

Notes:

- ✓ Try to use freshly prepared serum and plasma samples for this assay.
- ✓ No dilution or treatment of the sample is required, but plasma samples should be centrifuged to remove any fibrin/fibrinogen clumps.
- ✓ Samples should be stored at 2-8°C for 24 hours prior to analysis. For longer periods, store samples at -70°C, and avoid repeated freeze/thaw cycles.
- ✓ A 50 µl dead volume is required in addition to sample volume for multi-protein analysis (typically 1-5 µl).



- 1 Perform daily quality control assessment of instrumentation before analysis.
- 2 Load each sample into a specialized micro-sample cup for the clinical chemistry analyzer.
- 3 Select Urea/BUN test on display and run the analysis.
- 4 Collect and analyze the data.