

Jun 30, 2025

Ubiquitin enrichment assay

 [Nature Communications](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly4z7rlx9/v1>

Bishal Basak^{1,2}, Erika Holzbaur^{1,2}

¹Department of Physiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA, 19104;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

University of Pennsylva...



Bishal Basak

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly4z7rlx9/v1>

Protocol Citation: Bishal Basak, Erika Holzbaur 2025. Ubiquitin enrichment assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2ly4z7rlx9/v1>

Manuscript citation:

<https://www.biorxiv.org/content/10.1101/2024.06.12.598718v1.full>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 30, 2025

Last Modified: June 30, 2025

Protocol Integer ID: 221323

Keywords: ubiquitin enrichment, neuronal lysate, assay

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000350

Abstract

Ubiquitin enrichment assay from neuronal lysates

Protocol materials

 Signal-Seeker Ubiquitination Detection Kit **Cytoskeleton Inc. Catalog #BK161**



Plate neurons

- 1 Plate 3- 5 million neurons in a 10 cm dish at days in vitro (DIV) 0

Assay steps

- 2 When neurons are mature (DIV 6 or more), treat neurons depending on experimental needs
- 3 Post treatment, aspirate out media, wash neurons once with 1X fresh PBS
- 4 Lyse neurons in lysis buffer containing 50 mM HEPES (pH=7.4), 1 mM EDTA, 1 mM MgCl_2 , 25 mM NaCl, 0.5% Triton X-100, 20 $\mu\text{g}/\text{mL}$ Leupeptin, 2 mM DTT, 20 $\mu\text{g}/\text{mL}$ TAME, 2 $\mu\text{g}/\text{mL}$ Pepstatin A, 1 mM PMSF. Add Deubiquitination inhibitors provided with Ubiquitin Detection Kit (Cytoskeleton, BK161) to the lysis buffer at 1X concentration.

 Signal-Seeker Ubiquitination Detection Kit **Cytoskeleton Inc. Catalog #BK161**
- 5 Post lysis of the neurons, spin lysates at 17000x g for 10 mins, and discard the pellet.
- 6 Store 5% of the supernatant as input to be used for western blotting later.
- 7 Wash Ubiquitination affinity beads and Ubiquitination IP control beads (provided in the kit) with 1XPBST twice and 1XPBS once, and then incubate with equal volumes of the lysate for 4 hrs at 4⁰C.

Note

The assay was performed using Signal Seeker Ubiquitination Detection Kit, as mentioned in the manufacturer's protocol with some modifications.

 Signal-Seeker Ubiquitination Detection Kit **Cytoskeleton Inc. Catalog #BK161**

- 8 Post incubation, centrifuge beads at 5000x g at 4⁰C for 1 min.



- 9 Wash beads thrice with BlastR-2™ wash buffer for 5 minutes each in a rotating platform at 4°C and centrifuge at 5000Xg.
- 10 After final wash and centrifugation, remove buffer supernatant completely.
- 11 Incubate beads with Bead Elution Buffer for 5 mins at RT and then centrifuge in a spin column at 10000x g for 1 min to collect the immunoprecipitates.
- 12 Add 2-mercaptoethanol to each sample and boil samples at 95°C for 10 minutes to be used for western blotting.