

Jan 22, 2024

Ubiquitin immunoprecipitation using an anti-ubiquitin nanobody

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp36k1vzp/v1>

Cole S Sitron¹, Victoria A Trinkaus¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Cole Sitron

Max Planck Institute of Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp36k1vzp/v1>

Protocol Citation: Cole S Sitron, Victoria A Trinkaus, F Ulrich Hartl 2024. Ubiquitin immunoprecipitation using an anti-ubiquitin nanobody. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp36k1vzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 16, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 93785

Keywords: ASAPCRN, ubiquitin immunoprecipitation, ubiquitinated protein, ubiquitination, immunoprecipitation, protein, protein of interest, cell lysate, agarose bead

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol describes a method to detect ubiquitination on a protein of interest. This technique relies on immunoprecipitation (IP) of ubiquitinated proteins from a cell lysate through the use of an anti-ubiquitin nanobody coupled to agarose beads. The eluate can be run on SDS-PAGE to determine whether the protein of interest was recovered in the IP and, therefore, ubiquitinated.

Attachments



[io_Ubiquitin immunop...](#)

17KB

Materials

Materials

Cell culture

2-4 50-80% confluent wells of cells growing in a 6-well plate.

Equipment

Equipment

Bioruptor® Plus sonication device

NAME

Bioruptor®

BRAND

B01020001

SKU

<https://www.diagenode.com/en/p/bioruptor-plus-sonication-device>^{LINK}

- CLARIOstar Plus Plate Reader (BMG Labtech) (or equivalent)
- Microcentrifuge

Specialized Reagents

-  Ubiquitin Selector **NanoTag Biotechnologies Catalog #N2510**
-  MiniSpin Columns (50 units) **NanoTag Biotechnologies Catalog #A1001-L**
-  TrypLE® Express Enzyme (1X), phenol red **Thermo Fisher Catalog #12605036**
-  Pierce® Rapid Gold BCA Protein Assay Kit **Thermo Fisher Catalog #A53225**

Buffers

- Triton Lysis Buffer:

	A	B
	Tris pH 8	20 mM



A	B
NaCl	150 mM
Triton X-100	1%
cOmplete, Mini EDTA-free Protease Inhibitor Cocktail (Roche cat. no. 11873580001)	1X
Freshly-added 20 mM N-ethylmaleimide (2.5 mg/ml; Sigma Aldrich E3876-25G)	2.5 mg/ml
50 uM PR-619 (from 100 mM DMSO stock; Sigma-Aldrich 662141-25MG)	50 uM
Benzonase (Max Planck Institute of Biochemistry Core Facility)	7.5 U/ml

⊗ cOmplete™ EDTA-free Protease Inhibitor Cocktail **Roche Catalog #11873580001**

⊗ N-Ethylmaleimide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E3876**

⊗ DUB Inhibitor V, PR-619 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #662141**

▪ Triton Wash Buffer:

A	B
Tris pH 8	20 mM
NaCl	150 mM
EDTA	1 mM
Triton X-100	1%
SDS	0.1%
PR-619	50 uM

- 4X ⊗ NuPAGE™ LDS Sample Buffer (4X) **Thermo Fisher Scientific Catalog #NP0007** + 5% beta-mercaptoethanol (⊗ 2-Mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250** -100ml)
- 10% FBS (⊗ Fetal Bovine Serum **Gibco - Thermo Fisher Scientific Catalog #10270106**) in 1X PBS (diluted from ⊗ PBS (10X), pH 7.4 **Thermo Fisher Scientific Catalog #70011051**)
- 1X PBS



Cell collection and lysis

- 1 Remove the medium from the wells and trypsinize with  500 μ L TrypLE Express. 
- 2 Quench the TrypLE Express with  500 μ L 10% FBS and move the cells into a centrifuge tube. 
- 3 Pellet cells at  1500 x g for  00:03:00 at  4 $^{\circ}$ C .  3m
- 4 Wash cells with  1 mL PBS and pellet again. 
- 5 Resuspend cells in  150 μ L Triton Lysis Buffer and incubate  On ice for  00:05:00 .   5m
- 6 Sonicate in the BioRuptor for 5 cycles of  00:00:30 on and  00:00:30 off at  4 $^{\circ}$ C .  1m
- 7 Clarify the lysate by centrifugation  18000 x g for  00:10:00  4 $^{\circ}$ C .  10m
- 8 During the centrifugation, prepare BSA standards according to the Rapid Gold BCA Protein Assay Kit manual and begin equilibrating  50 μ L of Ubiquitin pan-selector resin (pipette with a cut p200 tip) in  500 μ L of lysis buffer. 
- 9 Collect supernatant and place into a new tube  On ice . 
- 10 Prepare a small 1:10 dilution of samples and perform Rapid Gold BCA Protein Assay according to manufacturer's instructions. 

Immunoprecipitation (IP)

1h 7m



11 Dilute samples into Triton Lysis Buffer in two tubes:  1 mg in  500 μ L (for IP) and  100 μ g in  25 μ L (for input).



12 Pellet the beads at  1000 x g for  00:02:00 at  4 $^{\circ}$ C .

2m



13 Pull off the supernatant, add the  1 mg of lysate to the beads, and incubate with rotation at  4 $^{\circ}$ C for  01:00:00 .

1h



14 During the incubation denature the input sample with  25 μ L 4X NuPAGE LDS Sample Buffer at  95 $^{\circ}$ C  00:05:00 .

5m



15 Pellet the beads at  1000 x g for  00:02:00 at  4 $^{\circ}$ C .

2m



16 Wash the beads with either  1 mL Triton Wash Buffer.



17 Pellet the beads and repeat for a total of 3 washes.



18 Remove the supernatant and add  100 μ L of 2X NuPAGE LDS Sample Buffer.



19 Elute by boiling at  95 $^{\circ}$ C for  00:05:00 .

5m

20 Transfer the beads to a MiniSpin column placed in a tube.



21 Centrifuge  3000 x g for  00:03:00 .

3m



22 The eluate has now been collected in the tube and is ready for SDS-PAGE analysis.



