

May 14, 2019

UC Davis - Glucose Protocol

DOI

dx.doi.org/10.17504/protocols.io.yj9fur6



Peter Havel¹

¹University of California, Davis

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yj9fur6>

External link: <https://mmpc.org/shared/document.aspx?id=88&docType=Protocol>

Protocol Citation: Peter Havel 2019. UC Davis - Glucose Protocol. protocols.io
<https://dx.doi.org/10.17504/protocols.io.yj9fur6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 26, 2019

Last Modified: May 14, 2019

Protocol Integer ID: 20833

Keywords: glucose, glucose oxidase, glucose oxidase to gluconic acid, glucose protocol summary, peroxidase with hba, peroxidase, glucose concentration, red quinoneimine dye, hydrogen peroxide reacts in the presence, hydrogen peroxide react, hydrogen peroxide, gluconic acid, proportional to the glucose concentration, oxidase, color, intensity of the color

Abstract

Summary:

Glucose is oxidized by glucose oxidase to gluconic acid and hydrogen peroxide. The hydrogen peroxide reacts in the presence of peroxidase with HBA and 4-aminoantipyrine forming a red quinoneimine dye. The intensity of the color formed is proportional to the glucose concentration and can be measured photometrically between 460 and 560 nm.

Materials


MATERIALS

 Calibrator **Fisher Diagnostics Catalog #TR1591-030**

 Reagents **Fisher Diagnostics Catalog #TR15103**

 PBS

 Microplate

 Platerreader

Reagent Preparation:

PBS – ready to use

Reagent – reconstitute with distilled water to make a 2X solution

Troubleshooting



- 1 Reconstitute powdered reagent with only 25 ml of distilled water to make a 2X solution.
- 2 Add 3 μ l of calibrator and sample to each well.
- 3 Add 150 μ l of PBS to each well. Read at 540 nm.

IMPORTANT: Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.

- 4 Add 150 μ l of 2X reagent to each well. Incubate at 37°C for 10 minutes. Read at 540 nm.
- 5 Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as $(\text{Sample Absorbance} \div \text{Calibrator Absorbance}) \times \text{Calibrator Concentration}$.