

May 14, 2019

UC Davis - Glucose Protocol

DOI

dx.doi.org/10.17504/protocols.io.yj9fur6



Peter Havel¹

¹University of California, Davis

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yj9fur6

External link: <https://mmpc.org/shared/document.aspx?id=88&docType=Protocol>

Protocol Citation: Peter Havel 2019. UC Davis - Glucose Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yj9fur6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 26, 2019

Last Modified: May 14, 2019

Protocol Integer ID: 20833

Keywords: glucose

Abstract

Summary:

Glucose is oxidized by glucose oxidase to gluconic acid and hydrogen peroxide. The hydrogen peroxide reacts in the presence of peroxidase with HBA and 4-aminoantipyrine forming a red quinoneimine dye. The intensity of the color formed is proportional to the glucose concentration and can be measured photometrically between 460 and 560 nm.

Materials


MATERIALS

 Calibrator **Fisher Diagnostics Catalog #TR1591-030**

 Reagents **Fisher Diagnostics Catalog #TR15103**

 PBS

 Microplate

 Platereader

Reagent Preparation:

PBS – ready to use

Reagent – reconstitute with distilled water to make a 2X solution



- 1 Reconstitute powdered reagent with only 25 ml of distilled water to make a 2X solution.
- 2 Add 3 µl of calibrator and sample to each well.
- 3 Add 150 µl of PBS to each well. Read at 540 nm.

IMPORTANT: Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.

- 4 Add 150 µl of 2X reagent to each well. Incubate at 37°C for 10 minutes. Read at 540 nm.
- 5 Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as $(\text{Sample Absorbance} \div \text{Calibrator Absorbance}) \times \text{Calibrator Concentration}$.