

May 14, 2019

## UC Davis - Insulin signaling pathway

DOI

<https://dx.doi.org/10.17504/protocols.io.yp2fvqe>



Fawaz G. Haj<sup>1</sup>

<sup>1</sup>University of California, Davis

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)



Lili Liang

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.yp2fvqe>

**External link:** <https://mmpc.org/shared/document.aspx?id=101&docType=Protocol>

**Protocol Citation:** Fawaz G. Haj 2019. UC Davis - Insulin signaling pathway. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.yp2fvqe>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 28, 2019

**Last Modified:** May 14, 2019

**Protocol Integer ID:** 20954

**Keywords:** Insulin signaling pathway, insulin receptor (IR), insulin substrate (IRS1/2), insulin receptor, defects in the insulin, signaling pathway, insulin, signaling pathway summary, pathway summary, map kinase, uc davi, activation state, evaluation of the activation state

## Abstract

### Summary:

This test is designated to determine defects in the insulin signaling pathway, through evaluation of the activation state of the insulin receptor (IR) and its substrate (IRS1/2), as well as downstream target, mainly Akt and MAP kinases.



## Materials

### MATERIALS

- ✕ Cell Lysis Buffer (10X) **Cell Signaling Technology Catalog #9803**
- ✕ 4-20% Tris-Glycine Gels **Invitrogen - Thermo Fisher Catalog #EC60285BOX**
- ✕ Tris-Glycine SDS Sample Buffer **Invitrogen - Thermo Fisher Catalog #LC2676**
- ✕ Tris-Glycine SDS Running Buffer **Invitrogen - Thermo Fisher Catalog #LC26755**
- ✕ Tris-Glycine Transfer Buffer **Invitrogen - Thermo Fisher Catalog #NP00061**
- ✕ Methanol **Fisher Scientific Catalog #A412P-4**
- ✕ PVDF 0.2 µm pore size **Invitrogen - Thermo Fisher Catalog #LC2002**
- ✕ WesternBreeze® Chemiluminescent Kit-Anti- Mouse **Invitrogen - Thermo Fisher Catalog #WB7104**
- ✕ WesternBreeze® Chemiluminescent Kit-Anti- Rabbit **Invitrogen - Thermo Fisher Catalog #WB7106**
- ✕ XCell SureLock® Mini-Cell and XCell IITM Blot Module Kit **Invitrogen - Thermo Fisher Catalog #EI0002**
- ✕ p-IR **Merck Millipore (EMD Millipore) Catalog #07-841**
- ✕ Total IR **Merck Millipore (EMD Millipore) Catalog #46-687**
- ✕ Insulin Receptor Substrate Antibody Sampler Kit **Cell Signaling Technology Catalog #3015**
- ✕ p-Akt Antibody **Cell Signaling Technology Catalog #4060**
- ✕ Total Akt Antibody **Cell Signaling Technology Catalog #9272**
- ✕ p44/42 MAPK (Erk1/2) (Thr202/Tyr204) Antibody Duet **Cell Signaling Technology Catalog #8201**
- ✕ Thermo Scientific Pierce\* BCA Protein Assay Kits **Thermo Scientific Catalog #23225**
- ✕ Cuvette 1.5ml **Fisher Scientific Catalog #14-955-127**

### Note:

**Cell Signaling Technology Pathway Database, [RRID:SCR\\_002071](#)**

**Thermo Fisher Scientific, [RRID:SCR\\_008452](#)**

**P-IR # 07-841, Cite this, (Millipore Cat# 07-841, [RRID:AB\\_568831](#))**

**Insulin Receptor Substrate Antibody Sampler Kit #3015, Cite this, (Cell Signaling Technology Cat# 3015, [RRID:AB\\_1196650](#))**

**Total Akt Antibody #9272, Cite this (Cell Signaling Technology Cat# 9272, [RRID:AB\\_329827](#))**

**p-Akt Antibody #4060, Cite this (Cell Signaling Technology Cat# 4060, [RRID:AB\\_2315049](#))**



**p44/42 MAPK (Erk1/2) (Thr202/Tyr204) Antibody Duet #8201**, Cite this, **(Cell Signaling Technology Cat# 8201, RRID:AB\_10695902)**

- 1 Fast mice for 16 hours by taking away food the day before (3:00pm)
- 2 The following day, give the mouse an intraperitoneal injection of insulin (10 U/kg) with a 27 G needle.
- 3 Use cervical dislocation to euthanize mice.
- 4 Collect maximum blood from portal vein and isolate plasma according to standard protocols or as desired by the P.I.
- 5 Quickly collect tissues (Liver, Muscle, Adipose and Pancreas). Each tissue should be divided into three portions, one portion should be snap frozen in liquid nitrogen, one portion should be kept into RNA later solution and the third one should be fixed into the appropriate fixative solution. Please note that the whole procedure of tissue collection should be done within 3 minutes maximum.
- 6 For western blotting, tissues will be lysed into the appropriate lysis buffer.
- 7 Overall and site-specific tyrosyl phosphorylation of individual components in insulin signaling such as IR and insulin receptor substrates (IRSs) will be determined by immunoprecipitation (IP) then probed with phospho-tyrosyl antibodies according to standard protocols. Tyrosine, Serine and/or threonine phosphorylation of other component of the insulin signaling pathway such as Akt and MAP kinases will be determined also according to the standard Western blotting protocols.
- 8 **Note:**  
Evaluation of the activation state of other component in the insulin signaling pathway is also possible upon special request. Extra charges may apply.