



May 15, 2019

## UC Davis - Intraperitoneal Insulin Tolerance Test

DOI

[dx.doi.org/10.17504/protocols.io.yp5fvq6](https://dx.doi.org/10.17504/protocols.io.yp5fvq6)



Fawaz G. Haj<sup>1</sup>

<sup>1</sup>University of California, Davis

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)



Lili Liang

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yp5fvq6>

External link: <https://mmpc.org/shared/document.aspx?id=84&docType=Protocol>

**Protocol Citation:** Fawaz G. Haj 2019. UC Davis - Intraperitoneal Insulin Tolerance Test. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.yp5fvq6>

**Manuscript citation:**

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 01, 2019

**Last Modified:** May 15, 2019

**Protocol Integer ID:** 20957

**Keywords:** Intraperitoneal Insulin Tolerance Test, glucose, intraperitoneal insulin tolerance test, intraperitoneal insulin tolerance test summary, sensitivity of insulin, insulin injection, insulin, ipitt, measurement of glucose, responsive tissues in the rodent, responsive tissue, test, sensitivity

## Abstract

**Summary:**

An intraperitoneal insulin tolerance test or ipITT is designed to determine the sensitivity of insulin-responsive tissues in the rodent. This is determined by measurement of glucose remaining in the circulation over time after a bolus ip insulin injection.

## Materials

**MATERIALS**

☒ Humalin® R Eli Lilly Catalog #R-100

☒ Insulin Syringes Fisher Scientific Catalog #14-826-79

☒ Saline Solution Fisher Scientific Catalog #L97753

☒ Easy Check Glucose test strips JRS Medical Catalog #00-101(new SKU 88982400)

☒ Easy Check Glucose monitor JRS Medical Catalog #Y4209(new SKU 88972401)

**Reagent Preparation:**

Dilute the stock solution (100 U/ml) with saline to 0.5 U/mL (1/200 dilution) by adding 5µl stock (100 U/mL) to 995 µl 0.9% (w/v) sterile saline

## Troubleshooting

- 1 Fast mice for 4 h only by taking away food early in the morning (7:00am).
- 2 Calibrate the glucose meter according to the manufacturer's instructions.
- 3 Deprive mice from water then remove approximately 5µl of blood (one drop) from the tail via a tail tip cut and transfer directly onto a glucose indicator strip.
- 4 Measure blood glucose immediately in a glucometer.
- 5 Give the mouse an intraperitoneal injection of insulin (0.5 U/kg) with a 27 G needle.
- 6 Continue to take blood samples from the initial tail cut before the insulin injection and at 15, 30, 45, 60 and 120 min.  
Between each of these time points, return the mouse to its cage and monitor it every minute.

7 **NOTE:**

1-The mouse is given an intraperitoneal injection with a 27G needle of insulin. Before performing the experiment, we will have to determine if the mouse strain is insulin resistant or glucose tolerant, so as to avoid giving the wrong dose of insulin.

2-For insulin resistant mice, a higher dose of insulin could be used as they do not respond as well as insulin sensitive mice to insulin. However, we will make sure that the dose of insulin used will not cause hypoglycemia.

3-At the end of the experiment, wipe tail with 70% alcohol and allow drying. Ensure that blood loss from the tail stopped before placing the animal back to its cage.