

May 15, 2019

UC Davis - LDL Protocol

DOI

dx.doi.org/10.17504/protocols.io.yrcfv2w



Peter Havel¹

¹University of California, Davis

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yrcfv2w

External link: <https://mmpc.org/shared/document.aspx?id=94&docType=Protocol>

Protocol Citation: Peter Havel 2019. UC Davis - LDL Protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yrcfv2w>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 01, 2019

Last Modified: May 15, 2019

Protocol Integer ID: 20996

Keywords: cholesterol, triglyceride, LDL Protocol



Abstract

Summary:

LDL and VLDL are separated from HDL using a precipitation reagent. Then the HDL fraction is measured for either TC or TG using the same reagents for total cholesterol or triglyceride.

Materials

MATERIALS

- ⊗ Calibrator **Fisher Diagnostics Catalog #TR43002**
- ⊗ TC Reagents **Fisher Diagnostics Catalog #TR13421**
- ⊗ TG Reagents **Fisher Diagnostics Catalog #TR22421**
- ⊗ 2X LDL/VLDL Precipitation Buffer **Abcam Catalog #ab105138**
- ⊗ PBS
- ⊗ Microplate
- ⊗ Platereader

Note:

Fisher Scientific, RRID:SCR_008452

Abcam, RRID:SCR_012931

- 1 Add 25µl 2X precipitation buffer to 25µl of sample using a positive displacement pipet.
- 2 Vortex and let sit at RT for 10 minutes.
- 3 Centrifuge at 2000×g for 10 minutes at 4°C.
- 4 Pipet supernatant into new tube, this is the HDL fraction.

- 5 Add 5 µl of calibrator and sample to each well.

IMPORTANT: Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.

- 6 Add 300 µl of TC or TG reagent to each well. Incubate at 37°C for 5 minutes. Read at 540 nm.

IMPORTANT: If samples are hemolyzed, pipet a blank well with 5µl sample and 300µl PBS

- 7 Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as (Sample Absorbance ÷ Calibrator Absorbance) × Calibrator Concentration.
- 8 HDL samples are diluted ½ so multiply these by 2 to get the final value. Subtract this from the total triglyceride or cholesterol value to get the LDL/VLDL value.