




May 15, 2019

UC Davis - Lipoprotein Binding Protein (LBP)-Endotoxemia Assay

DOI

dx.doi.org/10.17504/protocols.io.yrhfv36



Trina Knotts¹

¹University of California, Davis

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yrhfv36>

External link: <https://mmpc.org/shared/document.aspx?id=299&docType=Protocol>

Protocol Citation: Trina Knotts 2019. UC Davis - Lipoprotein Binding Protein (LBP)-Endotoxemia Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yrhfv36>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 01, 2019

Last Modified: May 15, 2019

Protocol Integer ID: 21001

Keywords: lipopolysaccharide binding protein (LBP), lipoprotein binding protein, lipoprotein, endotoxemia, measure of endotoxemia, bacterial lp, lbp, surrogate for bacterial lp, binding protein, protein

Abstract

Summary:

Plasma samples will be assayed for lipopolysaccharide binding protein (LBP) as surrogate for bacterial LPS/measure of endotoxemia via ELISA.

Materials

MATERIALS

Mouse LBP ELISA kit **Biometec Catalog ##43**

Shaker

Microplate spectrophotometer

MB grade water (diluent)

Reagent Preparation: *(According to the manufacturer.)*

Wash Buffer (PBS/ Tween 0.05%):

Dissolve 1 Tablet Phosphate buffered saline (PBS, vial 5) in 200ml distilled water -add 100 μ l Tween 20 (vial 7). (Prepared wash buffer is stable for 4 weeks at refrigerator).

Phosphate Buffered Saline (PBS):

Dilute 1 Tablet of vial 5 in 200 ml distilled water

Dilution Buffer:


Add content of the vial 6 to 50ml PBS (Buffer C). Prepare just before use. Store remaining dilution buffer after reconstitution at -20°C

Reference serum dilution:

Add 10 μ l distilled water to the vial 4. This contains $12.14 \pm 3.5\mu\text{g/ml}$ LBP (! new Reference for Lot #141016). For assay dilute 1:800 (10 μ l serum +7990 μ l dilution buffer and use 100 μ l/well.

LBP standards :

Firstly, pipette 30 μ l distilled water to the vial 3 for reconstitution and secondly add 270 μ l dilution buffer (C) to this vial and mix carefully, thirdly pipette 50 μ l from this vial to a new vial containing 450 μ l dilution buffer (C) and mix carefully. Finally this last vial contains 500 μ l standard dilution and containing 50ng/ml LBP = vial a. For standard curve prepare vial b-f and use vial a –f Prepare just before use. Store the standard at -20°C .





Troubleshooting

Before start

IMPORTANT: Check kit datasheet for lot-specific instructions that may modify general protocol.

- 1 Prepare kit reagents.
- 2 Dilute mouse plasma or serum samples 1:800.
- 3 Add 100 μ l of standards (50, 25, 12.5, 6.25, 3.12, 1.56 ng/ml) or diluted samples in duplicate into the corresponding wells of the precoated modules and incubate for one hour at room temperature and shaking (300rpm).
- 4 Wash 3X with Wash Buffer.
- 5 Add 100 μ l detecting antibody to each well and incubate at room temperature for 1 hour at shaker.
- 6 Wash 3X with Wash Buffer.
- 7 Add 100 μ l Substrate solution to each well. Incubate 12-14 min **in the dark** at room temperature **without** shaking. Cover with foil during incubation or place in drawer.
- 8 Add 100 μ l stopping solution to each well. Tap gently to mix.
- 9 Read absorbance at 450 nm (reference wave length 620 nm)
- 10 Calculate the LBP concentration by first plotting the OD means of standards (y-axis) and the LBP concentration (x-axis). Calculate the LBP concentration from the mean OD of samples from the standard curve and multiply with dilution factor.