

May 15, 2019

UC Davis - Metabolomics: Sample preparation for Lipidomics

DOI

dx.doi.org/10.17504/protocols.io.ytpfwmn



Oliver Fiehn¹

¹University of California, Davis

Mouse Metabolic Pheno...

Metabolomics Protocols ...



Lili Liang

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.ytpfwmn

External link: <https://mmpc.org/shared/document.aspx?id=122&docType=Protocol>

Protocol Citation: Oliver Fiehn 2019. UC Davis - Metabolomics: Sample preparation for Lipidomics. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ytpfwmn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 05, 2019

Last Modified: May 15, 2019

Protocol Integer ID: 21071

Keywords: Lipidomics

Abstract

Summary:

This SOP describes sample extraction and sample preparation for lipid profiling by liquid chromatography / quadrupole time of flight mass spectrometry (LC-QTOF) or nanoelectrospray ion trap-FTICR MS.

Materials

MATERIALS

-  Centrifuge **Eppendorf Catalog #5415 D**
-  Calibrated pipettes 1-200 ul and 100-1000ul **Catalog #1-200 ul and 100-1000ul**
-  Eppendorf tubes 1.5 mL uncolored **Eppendorf Catalog #022363204**
-  ThermoElectron Neslab RTE 740 cooling bath at – 20°C **Catalog #RTE 740 cooling bath**
-  MiniV ortexer **VWR International (Avantor) Catalog #58816-121**
-  Orbital Mixing Chilling/Heating Plate **Torrey Pines Scientific Instruments**
-  Speed vacuum concentration system **Labconco Centrivap cold trap**
-  Eppendorf tips for organic solvents such as acetonitrile methanol and MTBE
-  Glass Amber Vials **National Scientific Catalog #C4000-2W**
-  Glass Inserts **Supelco Catalog #27400-U**
-  Blue Tops for Vials **Agilent Technologies Catalog #5182-0717**
-  Crushed ice
-  Nitrogen line with pipette tip
-  Pure water
-  MTBE: Sigma Chromasolv 99.8% for HPLC 100mL (smallest available) (34875-100mL) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #34875-100ML**
-  Methanol: J.T. Baker LC/MS Grade (9830-03) **JT Baker Catalog #9830-03**
-  CUDA (12-[[[(cyclohexylamino)carbonyl]amino]- dodecanoic acid) **Cayman Chemical Company Catalog #10007923**

Note:

Sigma-Aldrich RRID:SCR_008988

Cayman Chemical RRID:SCR_008945

1 **Starting material:**

Plasma/serum: 30 µl sample volume or aliquot

2 **Sample Preparation:**

Switch on bath to pre-cool at -20°C ($\pm 2^{\circ}\text{C}$ validity temperature range)

Extraction solvents

- ♦ Purge both MeOH and MTBE for 5 min with nitrogen.
- ♦ Store solvents in the -20°C freezer to pre-chill

Homogenization and extraction

- ♦ Thaw plasma on ice, and gently rotate or invert the blood samples for about 10s to obtain a homogenized sample.
- ♦ Take out 60 µL and add 220 µL cold MeOH. Add 5 µL of QC mix as internal standard (see SOP "QC mix for LC-MS lipid analysis").
- ♦ Vortex each sample for 10s, keeping the rest on ice
- ♦ Add 750 µL MTBE
- ♦ Vortex for 10s
- ♦ Shake for 6min at 4°C
- ♦ Add 187.5 µL distilled water
- ♦ Vortex for 20s
- ♦ Centrifuge for 2 min @ 14000 rcf
- ♦ Remove supernatant, splitting into two aliquots of 300 µL, keeping one at -20°C for backup
- ♦ Dry samples to complete dryness in the speed vacuum concentration system

Preparation for analysis

- ♦ Re-suspend dry samples in 70 µL MeOH containing CUDA (10 µM), degassed using the above method.
- ♦ Transfer 30 µL to two separate amber glass vial with micro-insert. Cap vials with Agilent blue top.
- ♦ Use independent vials for positive and negative mode acquisitions.


3 **Quality assurance**

- ♦ For each sequence of sample extractions, perform one blank negative control extraction by applying the total procedure (i.e. all materials and plastic ware) without biological sample.
- ♦ Use one commercial plasma/serum pool sample per 10 authentic subject samples as control. If no suitable commercial blood sample is available, prepare a large pool sample during the thawing/mixing step by aliquoting 100 µl per 1 ml plasma sample, and aliquot such pool sample for 1 pool extract per 10 authentic subject samples.

- ◆ Prepare at least six NIST plasma extracts in the same manner as positive controls

4

Final Protocol



IMPORTANT: To prevent contamination disposable material is used. To prevent inhalation of toxic ether vapor, use fume hood during lipid extraction.

DISPOSAL OF WASTE: Collect all chemicals in appropriate bottles and follow the disposal rules. Collect residual plasma / serum samples in specifically designed red 'biohazard' waste bags.