




May 15, 2019

UC Davis - Superoxide Dismutase

DOI

<https://dx.doi.org/10.17504/protocols.io.ywdfxa6>



Peter Havel¹

¹University of California, Davis

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ywdfxa6>

External link: <https://mmpc.org/shared/document.aspx?id=129&docType=Protocol>

Protocol Citation: Peter Havel 2019. UC Davis - Superoxide Dismutase. [protocols.io](https://dx.doi.org/10.17504/protocols.io.ywdfxa6)
<https://dx.doi.org/10.17504/protocols.io.ywdfxa6>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 06, 2019

Last Modified: May 15, 2019

Protocol Integer ID: 21157

Keywords: Superoxide Dismutase, superoxide dismutase summary, superoxide dismutase, significant amounts of superoxide dismutase, dismutation of superoxide radical, superoxide radical, concentration of superoxide, superoxide, oxidative stress, antioxidant capability, xanthine oxidase, other neurodegenerative disorders such as alzheimer, quantification of sod activity, other neurodegenerative disorder, cayman chemical sod assay kit, familial amyotrophic lateral sclerosis, sod activity from plasma, enzyme in cell, mutations in sod account, alzheimer, enzyme, types of sod, reliable assay, sod activity, controlled sod standard, parkinson

Abstract

Summary:

Significant amounts of superoxide dismutase (SOD) in cellular and extracellular environments are crucial for the prevention of diseases linked to oxidative stress. Mutations in SOD account for approximately 20% of familial amyotrophic lateral sclerosis (ALS) cases. SOD also appears to be important in the prevention of other neurodegenerative disorders such as Alzheimer's, Parkinson's, and Huntington's Diseases. The reaction catalyzed by SOD is extremely fast, having a turnover of $2 \times 10^6 \text{ M}^{-1}\text{sec}^{-1}$ and the presence of sufficient amounts of the enzyme in cells and tissues typically keeps the concentration of superoxide very low. Quantification of SOD activity is therefore essential in order to fully characterize the antioxidant capabilities of a biological system. The Cayman Chemical SOD Assay kit is a fast and reliable assay for the measurement of SOD activity from plasma, serum, tissue homogenates, and cell lysates. SOD activity is assessed by measuring the dismutation of superoxide radicals generated by xanthine oxidase and hypoxanthine in a convenient 96 well format. A key feature of the kit is the inclusion of a quality-controlled SOD standard. The standard curve generated using this enzyme provides a means to accurately quantify the activity of all three types of SOD (Cu/Zn-, Mn-, and Fe-SOD). Each kit contains sufficient reagents to assay 41 samples in duplicate and includes assay buffer, sample buffer, radical detector, SOD (standard), xanthine oxidase, a 96 well plate, and complete instructions.



Materials

MATERIALS

⊗ Assay Kit **Cayman Chemical Company Catalog #706002**

⊗ Assay Buffer

⊗ Sample Buffer

⊗ Radical Detector

⊗ Standard

⊗ Xanthine Oxidase

Note:

Cayman Chemical [RRID:SCR_008945](#)

- 1 **SOD Standard Wells** - add 200 μ l of the diluted Radical Dector and 10 μ l of Standard (tubes A-G) per well in the designated wells on the plate (see sample plate format, Figure 2, page 11)

	1	2	3	4	5	6	7	8	9	10	11	12
A	A	A	S2	S2	S10	S10	S18	S18	S26	S26	S34	S34
B	B	B	S3	S3	S11	S11	S19	S19	S27	S27	S35	S35
C	C	C	S4	S4	S12	S12	S20	S20	S28	S28	S36	S36
D	D	D	S5	S5	S13	S13	S21	S21	S29	S29	S37	S37
E	E	E	S6	S6	S14	S14	S22	S22	S30	S30	S38	S38
F	F	F	S7	S7	S15	S15	S23	S23	S31	S31	S39	S39
G	G	G	S8	S8	S16	S16	S24	S24	S32	S32	S40	S40
H	S1	S1	S9	S9	S17	S17	S25	S25	S33	S33	S41	S41

A-G = Standards
S1-S41 = Sample Wells

Figure 2. Sample plate format

- 2 **Sample Wells** - add 200 μ l of the diluted Radical Dector and 10 μ l of sample to the wells.

NOTE: If using ad inhibitor, add 190 μ l of the diluted Radical Dector, 10 μ l of inhibitor, and 10 μ l of sample to the wells. The amount of sample added to the well should always be 10 μ l. Samples should be diluted with Sample Buffer (dilute) or concentrated with an Amicon centrifuge concentrator with a molecular weight cut-off of 10,000 to bring the enzymatic activity to fall within the standard curve range.

- 3 Initiate the reactions by adding 20 μ l of diluted Xanthine Oxidase to all the wells you are using. Make sure to note the precise time you started and add the Xanthine Oxidase as quickly as possible.

NOTE: If assaying sample backgrounds, add 20 μ l of Sample Buffer instead of Xanthine Oxidase.

- 4 Carefully shake the 96-well plate for a few seconds to mix. Cover with the plate cover.



- 5 Incubate the plate on a shake for 20 minutes at room temperature. Read the absorbance at 440-460 nm using a plate reader.

6 Calculation

1. Calculate the average absorbance of each standard and sample. If assayed, subtract sample background absorbance from the sample.
2. Divide standard A's absorbance by itself and divide standard A's absorbance by all the other standards and samples absorbances to yield the linearized rate (LR) (i.e., LR for Std A = Abs. Std A/Abs Std A; LR for Std B = Abs Std A/Abs Std B).
3. Plot the linearized SOD standard rate (LR) (from step 2 above) as a function of final SOD Activity (U/ml) from Table 1. See Figure 3 (on page 17) for a typical standard curve.

Tube	SOD Stock (μl)	Sample Buffer (μl)	Final SOD Activity (U/ml)
A	0	1,000	0
B	20	980	0.025
C	40	960	0.05
D	80	920	0.1
E	120	880	0.15
F	160	840	0.2
G	200	800	0.25

Table 1. Superoxide Dismutase Standards




Figure 3. Superoxide Dismutase Standard Curve

4. Calculate the SOD activity of the samples using the equation obtained from the linear regression of the standard curve substituting the linearized rate (LR) for each sample. One unit is defined as the amount of enzyme needed to exhibit 50% dismutation of the superoxide radical. SOD activity is standardized using the cytochromic c and xanthine oxidase coupled assay

$$\text{SOD (U/ml)} = \left[\left(\frac{\text{sample LR} - \text{y-intercept}}{\text{slope}} \right) \times \frac{0.23 \text{ ml}}{0.01 \text{ ml}} \right] \times \text{sample dilution}$$