

May 15, 2019

## UC Davis - Triglyceride Protocol

DOI

[dx.doi.org/10.17504/protocols.io.yw3fxgn](https://dx.doi.org/10.17504/protocols.io.yw3fxgn)



Peter Havel<sup>1</sup>

<sup>1</sup>University of California, Davis

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)



Lili Liang

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.yw3fxgn](https://dx.doi.org/10.17504/protocols.io.yw3fxgn)

External link: <https://mmpc.org/shared/document.aspx?id=91&docType=Protocol>

**Protocol Citation:** Peter Havel 2019. UC Davis - Triglyceride Protocol. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.yw3fxgn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** March 07, 2019

**Last Modified:** May 15, 2019

**Protocol Integer ID:** 21179

**Keywords:** Triglyceride,

## Abstract

### Summary:

Triglycerides are enzymatically hydrolyzed by lipase to free fatty acids and glycerol. The glycerol is phosphorylated by adenosine triphosphate (ATP) with glycerol kinase (GK) to produce glycerol-3- phosphate and adenosine diphosphate. Glycerol-3-phosphate is oxidized by dihydroxyacetone phosphate (DAP) by glycerolphosphate oxidase producing hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>). In a Trinder5 type color reaction catalyzed by peroxidase, the H<sub>2</sub>O<sub>2</sub> reacts with 4-aminoantipyrine (4-AAP) and 3,5-dichloro-2- hydroxybenzene sulfonate (DHBS) to produce a red colored dye. The absorbance of this dye is proportional to the concentration of triglycerides present in the sample.

## Materials

### MATERIALS

 Calibrator **Fisher Diagnostics Catalog #TR43002**

 Reagents **Fisher Diagnostics Catalog #TR22203**

 PBS

 Microplate

 Platereader

### Reagent Preparation:

PBS – ready to use

Reagent – reconstitute with distilled water to make a 2X solution



- 1 Reconstitute powdered reagent with only 25 ml of distilled water to make a 2X solution.
- 2 Add 3 µl of calibrator and sample to each well.
- 3 Add 150 µl of PBS to each well. Read at 540 nm.

**IMPORTANT:** *Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.*

- 4 Add 150 µl of 2X reagent to each well. Incubate at 37°C for 5 minutes. Read at 540 nm.
- 5 Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as  $(\text{Sample Absorbance} \div \text{Calibrator Absorbance}) \times \text{Calibrator Concentration}$ .