

May 15, 2019

UC Davis - Urea Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.yw6fxhe>



Peter Havel¹

¹University of California, Davis

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yw6fxhe>

External link: <https://mmpc.org/shared/document.aspx?id=100&docType=Protocol>

Protocol Citation: Peter Havel 2019. UC Davis - Urea Protocol. protocols.io <https://dx.doi.org/10.17504/protocols.io.yw6fxhe>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 07, 2019

Last Modified: May 15, 2019

Protocol Integer ID: 21182

Keywords: Urea Protocol, urea concentration, urea protocol summary, enzyme methodology, enzyme, reagent, absorbance, assay, nad, protocol summary, reaction

Abstract

Summary:

The enzyme methodology employed in this reagent is based on the reaction first described by Talke and Schubert. To shorten and simplify the assay, the calculations are based on the discovery of Tiffany, et al. that urea concentration is proportional to absorbance change over a fixed time interval. The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm as NADH is converted to NAD.

Materials

MATERIALS

 Calibrator **Fisher Diagnostics Catalog #TR43002**

 Reagents **Fisher Diagnostics Catalog #TR12003**

 Microplate

 Platereader



- 1 Add 3 μ l of calibrator and sample to each well.
- 2 Add 300 μ l of reagent to each well. Incubate at 37°C for 30 seconds. Read at 340 nm.

IMPORTANT: *Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.*

- 3 Incubate at 37°C for 60 seconds. Read at 340 nm.
- 4 Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as $(\text{Sample Absorbance} \div \text{Calibrator Absorbance}) \times \text{Calibrator Concentration}$.