

Apr 23, 2024

Ultra-deep ATAC-seq for sorted neurons

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyw35rvx9/v1>

Andrea Kriz¹, Alisa Mo¹

¹Boston Children's Hospital

SMaHT



Andrea Kriz

Boston Children's Hospital

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyw35rvx9/v1>

Protocol Citation: Andrea Kriz, Alisa Mo 2024. Ultra-deep ATAC-seq for sorted neurons. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2lyw35rvx9/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 18, 2024

Last Modified: April 23, 2024

Protocol Integer ID: 98347

Keywords: seq for sorted neurons isolation, sorted neurons isolation, preparation of atac, somatic mutation detection, atac, sequencing, seq, neuron, seq library

Funders Acknowledgements:


Somatic Mosaicism Across Human Tissues (SMaHT)

Grant ID: UG3NS132138

Abstract


Isolation of neurons from frozen post-mortem human brain tissue and preparation of ATAC-seq libraries for ultra-deep sequencing and somatic mutation detection.


Protocol materials















Isolation of Nuclei from Adult Human Brain Tissue (modified from Allen Human Brain Tissue PF0291)

1 Reagent list:

- ⊗ Tagmentase (Unloaded) **Diagenode Catalog #C01070010-20**
- ⊗ Tagmentation Buffer (2X) **Diagenode Catalog #C01019043-5000**
- ⊗ Anti-NeuN Antibody, clone A60, Alexa Fluor®488 conjugated **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAB377X**
- ⊗ DAPI Solution (1 mg/mL) **Thermo Fisher Scientific Catalog #62248**
- ⊗ DNA Clean & Concentrator-5 (capped) **Zymo Research Catalog #D4014**
- ⊗ Nextera XT Index Kit v2 Set B (96 indexes, 384 samples) **Illumina, Inc. Catalog #FC-131-2002**
- ⊗ NEBNext High-Fidelity 2X PCR Master Mix - 250 rxns **New England Biolabs Catalog #M0541L**

2 **Buffer preparation:** On the day prior to sorting, prepare the following buffers:

a. Nuclei Isolation Media (NIM):

A	B	C	D	E
Formula(in milliQ water)	Ingredients:	For 1.5 mL(1 ATAC-seq brain)	For 5 mL(3 ATAC-seq brains)	For 20 mL
0.25M Sucrose	Sucrose	0.13g	0.43 g	1.72 g
25mM KCl	2M KCl	18.8 µL	62.7 uL	250.8 uL
5mM MgCl ₂	1M MgCl ₂	7.5 uL	25 uL	100 uL
10mM Tris-HCl, pH 8	1M Tris-HCl, pH 8	15 uL	50 uL	200 uL
0.1% Triton X-100	10% Triton X-100	15 uL	50 uL	200 uL
Water	Water	~1410 uL	~4.7 mL	18.8 mL

Filter through 0.22µm filter and store at 4°C. Can keep for ~1-2 weeks.

b. Blocking buffer:



	A	B	C	D	E
	Formula (in milliQ water)	Ingredients:	For 5 mL	For 10 mL	For 20 mL
	Water	Water	4.3 mL	8.6 mL	17.2 mL
	1X PBS	10X PBS	500 uL	1000 uL	2000 uL
	0.8% BSA	20% BSA	200 uL	400 uL	800 uL

Store at 4°C. Can keep for ~1-2 weeks.

c. 1X Tagmentation buffer w/ BSA for ATAC-seq:

	A	B	C
	Ingredients:	For 4 ATAC- seq reactions	For 6 ATAC- seq reactions
	2X Tagmentation buffer (Diagenode)	220 uL	330 uL
	5% digitonin	0.88 uL	1.3 uL
	10% Tween- 20	4.4 uL	6.6 uL
	PBS	145 uL	218 uL
	20% BSA	22 uL	33 uL
	Water	25.5 uL	38.3 uL

Place at 4°C until ready to use.

- 3 **Load Tn5:** Mix 1 µL of annealed oligo A/oligo Rev with 1 µL of the annealed oligo B/oligo Rev as per Diagenode instructions.
- 4 Add 2 µL of Diagenode Tagmentase

- 5 Incubate at RT for 30 min
- 6 Add 2 μ L glycerol. Place at -20°C until ready for use. Can store for a few week.
- 7 **Tissue preparation:**
Chill tweezers, scalpel, and dissecting plate (cell culture plate) on dry ice. Chill dounce homogenizers with B pestle on ice. Make sure swinging-bucket centrifuge is cooled to 4°C .
- 8 Spray down all surfaces and pipettes with 70% Ethanol and DNA Away
- 9 Make homogenization buffer (NIM + additives): Homogenization buffer: 5mL NIM + 5 μ L 1M DTT (1 mM final concentration) + Mini Protease inhibitor (1/2 tablet) + 25 μ L 100 mM spermidine (0.5 mM final concentration).
- 10 Scrape ~10mg frozen brain tissue onto cell culture plate, either on top of dry ice or in a cryostat at -20°C . Add 1.5 mL Homogenization buffer to tissue. Add tissue/buffer mix to Dounce homogenizer.
- 11 Dounce 15-20 strokes with B pestle. If the tissue is large, use more homogenization buffer (we have used up to 5 mL) and dounce with A pestle followed by B pestle.
- 12 Filter through 40 μ M filter into Eppendorf tube.
- 13 Spin in 4°C swinging-bucket centrifuge for 10min @ 900xg
- 14 While spin is going on, clean homogenizers by rinsing thoroughly with ddH₂O, and then soaking them in 20% bleach for at least 20 minutes.
- 15 Make immunostaining buffer (Blocking buffer + 1:1000 dilution NeuN-488).
- 16 Remove supernatant. Do not disturb the pellet.
- 17 Add 1mL Blocking buffer + NeuN. Resuspend pellet gently with pipette. Rotate end-to-end in the cold room for 15-20 minutes.



- 18 Spin in 4°C swinging-bucket centrifuge for 5min @ 400xg
- 19 Make DAPI solution. First, dilute stock DAPI (1mg/mL) 1:15 in blocking buffer. Then, add 1 μ L of the diluted dapi solution to 1 mL blocking buffer.
- 20 Remove supernatant from pellet.
- 21 Add 1mL Blocking buffer + Dapi and again resuspend. Strain through 40 μ M filter into new Eppendorf tube immediately before nuclei sorting.
- 22 Prepare tubes for ATAC-seq: Add 60 μ L of 1X Tagmentation buffer into each tube for ATAC-seq (2 tubes per brain)

ATAC-seq (modified from Omni-ATAC protocol (Corces et al., 2017))

- 23 Sort 10,000 NeuN+ nuclei into 60 μ L of tagmentation buffer per ATAC-seq (2 replicates per brain sample). Proceed immediately after sorting to ATAC-seq
- 24 After sorting, spin each tube in 4 deg swinging bucket centrifuge at 500xg for 5 minutes.
- 25 Remove ~75 μ L of supernatant leaving ~10 μ L at bottom. Tap tube gently to mix pellet.
- 26 Add 45 μ L of 1X Tagmentation buffer. Do not mix.
- 27 Add 1 μ L of loaded Tn5. Pipette gently 3-4 times and tap to mix.
- 28 Incubate at 37°C for 30 minutes in the Thermomixer without mixing. Gently tap tube at ~10 minutes and at ~20 minutes into the incubation to mix.
- 29 Stop tagmentation reaction by adding 250 μ L of DNA binding buffer from Zymo Clean and Concentrator 5 kit. Add to sample tube, pipette 6x or vortex until homogenous.
- 30 Follow Zymo kit instructions to purify DNA:



1. Add solution to spin-column tube. Spin at 10,000xg for 30 seconds.
2. Wash x 200 μ L wash buffer. Spin at 10,000xg for 30 seconds.
3. Repeat wash x 200 μ L wash buffer. Spin at 10,000xg for 30 seconds.
4. Transfer spin-column tube into Eppendorf tube. Add 20 μ L of water. Let sit for 1 minute. Then spin at 10,000xg for 1 minute.

31 Prepare PCR reaction on ice:

- 2.5 μ L i5 primer
- 2.5 μ L i7 primer
- 25 μ L NEBNext High-Fidelity 2X PCR Master Mix, thaw on ice and invert until homogenous
- 20 μ L purified tagmented DNA

32 PCR amplification:

1. 72°C x 5 minutes
2. 98°C x 30 seconds
3. 98°C x 10 seconds
4. 63°C x 30 seconds
5. 72°C x 1 min
6. Loop back to step [3] 8 times (9 cycles total)
7. 72°C x 1 minute
8. 4°C hold

33 Clean the library using Zymo clean and concentrator 5 kit. Add 250 μ L of DNA binding buffer to PCR reaction and follow Zymo kit directions to purify.

34 Elute in 40-50 μ L 0.1X TE (or Zymo elution buffer). Let sit for 1 minute, spin for 1 minute.

35 Run a 1:4 dilution of library on HS5000 Tapestation


36

Expected result


Library traces from frozen brain tissues are quite variable. In general, a nucleosomal ladder should be apparent with a sub-nucleosomal peak around ~170bp (keep in mind the size of the adapters is around ~100bp, so this corresponds to an insert size of ~70bp), mono-nucleosomal peak around ~300bp, di-nucleosomal peak around ~450bp, tri-nucleosomal peak around ~600bp and so on.

The strength of the nucleosome ladder relative to larger molecular weight fragments can differ greatly between samples (and even between replicates). We have had success sequencing ATAC libraries with traces seen below. Note that for brain samples, the mono-nucleosomal peak should be larger than the sub-nucleosomal peak. An overly large sub-nucleosomal peak can indicate degradation of chromatin structure


Below are three representative traces of libraries we have successfully sequenced:



ATAC library with very apparent nucleosomal ladder



ATAC library with a greater quantity of large molecular weight fragments relative to nucleosomal ladder



ATAC library with a nucleosomal ladder as well as a comparable quantity of large molecular weight fragments

37 Store libraries at -20°C.



Protocol references

Citation

Corces MR, Trevino AE, Hamilton EG, Greenside PG, Sinnott-Armstrong NA, Vesuna S, Satpathy AT, Rubin AJ, Montine KS, Wu B, Kathiria A, Cho SW, Mumbach MR, Carter AC, Kasowski M, Orloff LA, Risca VI, Kundaje A, Khavari PA, Montine TJ, Greenleaf WJ, Chang HY (2017). An improved ATAC-seq protocol reduces background and enables interrogation of frozen tissues.

<https://doi.org/10.1038/nmeth.4396>

LINK

Citation

LeinU01 BRAIN grant (2026). Isolation of Nuclei from Adult Human Brain Tissue for 10x Genomics Platform. protocols.io.

<https://protocols.io/view/isolation-of-nuclei-from-adult-human-brain-tissue-y6rfzd6>

LINK

Citation

Matevossian A, Akbarian S (2008). Neuronal nuclei isolation from human postmortem brain tissue.

<https://doi.org/pii:914.10.3791/914>

LINK



Citations

Corces MR, Trevino AE, Hamilton EG, Greenside PG, Sinnott-Armstrong NA, Vesuna S, Satpathy AT, Rubin AJ, Montine KS, Wu B, Kathiria A, Cho SW, Mumbach MR, Carter AC, Kasowski M, Orloff LA, Risca VI, Kundaje A, Khavari PA, Montine TJ, Greenleaf WJ, Chang HY. An improved ATAC-seq protocol reduces background and enables interrogation of frozen tissues.

<https://doi.org/10.1038/nmeth.4396>

LeinU01 BRAIN grant. Isolation of Nuclei from Adult Human Brain Tissue for 10x Genomics Platform

<https://protocols.io/view/isolation-of-nuclei-from-adult-human-brain-tissue-y6rfzd6>

Matevossian A, Akbarian S. Neuronal nuclei isolation from human postmortem brain tissue.

<https://doi.org/pii:914.10.3791/914>