



Jun 27, 2024

🌐 Ultra-high field, multi-echo MRI acquisition for resting-state functional connectivity

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpbbj1lzp/v1>

Elizabeth Rizor¹, Dr. Scott Grafton²

¹UC Santa Barbara; ²University of California, Santa Barbara



Elizabeth Rizor

UC Santa Barbara

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpbbj1lzp/v1>

External link: <https://socalkinesia.org/high-resolution-connectomics/>

Protocol Citation: Elizabeth Rizor, Dr. Scott Grafton 2024. Ultra-high field, multi-echo MRI acquisition for resting-state functional connectivity. [protocols.io https://dx.doi.org/10.17504/protocols.io.dm6gpbbj1lzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gpbbj1lzp/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 02, 2022

Last Modified: June 27, 2024

Protocol Integer ID: 57737

Keywords: Functional Magnetic Resonance Imaging (f..., resting state functional connectivity (r..., Parkinson's Disease, motor circuits, basal ganglia, ASAPCRN, ultra-high field, multi-echo, minute mri protocol, minute mri protocol at the university, mri acquisition sequence, state imaging, head coil, functional connectivity, planar imaging, parkinson, state functional connectivity, imaging, epi sequence for resting, closed loop motor circuits in healthy control, cmrr spin, image acquisition

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020-519

Abstract

This is an MRI acquisition sequence used to examine resting-state functional connectivity of open and closed loop motor circuits in healthy controls and those with Parkinson's Disease. Scans include a Localizer, anatomical T1 MPRAGE, CMRR spin-echo AP and PA echo-planar imaging (EPI) images for field map generation, and a modified CMRR multi-echo EPI sequence for resting-state imaging. All participants complete a 30-minute MRI protocol at the University of Southern California's Center for Image Acquisition, where they are scanned with a Siemens 7T Terra scanner and a research-only 8-channel transmit 32-channel receive Nova 8Tx/32Rx head coil (Nova Medical Inc.)

Attachments



[highres_fc_acquisiti...](#)

73KB

Materials

Equipment

- Siemens 7T MAGNETOM Terra
- Nova Medical 8-channel pTx system and 8Tx/32Rx head coil



Protocol references

Puckett, A.M., Bollmann, S., Poser, B. A., Palmer, J., Barth, M., & Cunningham, R. (2018). Using multi-echo simultaneous multi-slice (SMS) EPI to improve functional MRI of the subcortical nuclei of the basal ganglia at ultra-high field (7T). *NeuroImage*, 172, 886–895. <https://doi.org/10.1016/j.neuroimage.2017.12.005>

Shao, X., Yan, L., Ma, S. J., Wang, K., & Wang, D. J. (2021). High-Resolution Neurovascular Imaging at 7T: Arterial Spin Labeling Perfusion, 4-Dimensional MR Angiography, and Black Blood MR Imaging. *Magnetic Resonance Imaging Clinics*, 29(1), 53-65.

Wu, X., Auerbach, E. J., Vu, A. T., Moeller, S., Van de Moortele, P. F., Yacoub, E., & Uğurbil, K. (2019). Human Connectome Project-style resting-state functional MRI at 7 Tesla using radiofrequency parallel transmission. *Neuroimage*, 184, 396-408.