

Aug 20, 2019

Ultracentrifugal separation of HDL alone and calculation of non-HDL

DOI

dx.doi.org/10.17504/protocols.io.32mgqc6



Daniel Teupser¹, Jan Breslow¹

¹Rockefeller University

Diabetic Complications Consortium
Tech. support email: rmcindoe@augusta.edu



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.32mgqc6>

External link: <https://www.diacomp.org/shared/document.aspx?id=17&docType=Protocol>

Protocol Citation: Daniel Teupser, Jan Breslow 2019. Ultracentrifugal separation of HDL alone and calculation of non-HDL. protocols.io <https://dx.doi.org/10.17504/protocols.io.32mgqc6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 11, 2019

Last Modified: August 20, 2019

Protocol Integer ID: 24365

Keywords: Ultracentrifugal separation of HDL alone, calculation of non-HDL, cardiovascular, various lipid fractions from blood plasma, various lipid fraction, using ultracentrifugation, ultracentrifugation, hdl summary, lipid, hdl, blood plasma, diabetic complication

Abstract

Summary:

This protocol is used to isolate the various lipid fractions from blood plasma using ultracentrifugation. The actual measured concentrations are performed separately once the isolations are complete.

NOTE: *This protocol IS applicable for ApoE knockout mice.*

Diabetic Complication:



Cardiovascular

Materials

MATERIALS

⊗ Beckman Optima TL tabletop ultracentrifuge **Beckman Coulter**

⊗ Beckman 7×20 mm thick walled ultracentrifuge tube **Beckman Coulter Catalog #343621**

⊗ Hamilton Syringe (100 ul)

⊗ KBr Solution

⊗ Phosphate Buffered Saline

Reagent/Material	Quantity Required
Beckman Optima TL tabletop ultracentrifuge	
Beckman 7x20 mm, thick walled ultracentrifuge tube	2
Hamilton Syringe (100 ul)	1
KBr Solution	1 ml
Phosphate Buffered Saline	1 ml

Troubleshooting

Safety warnings

! WARNING.

The use of an ultracentrifuge should only be performed by qualified technicians/personnel.

- 1 Add 60 μ l of plasma to Beckman ultracentrifugation tube (7 \times 20 mm; thick walled; polyallomer; cat. # 343621).
- 2 Layer 60 μ l of PBS on top of the plasma and place tubes in a TLA100 rotor.
- 3 Spin for 3 hours Beckman Optima TL tabletop ultracentrifuge at 70,000 rpm, 4°C.
- 4 Using a 100 μ l Hamilton syringe, carefully remove the bottom 60 μ l and transfer to a new Beckman tube labeled with the sample number. Discard the upper portion of the sample (impure VLDL). Between samples rinse the Hamilton syringe with distilled water.
- 5 Add 60 μ l KBr solution (density = 1.12 g/ml) to make a final density of 1.063 g/ml) and mix 5 to 6 times up and down with the same pipette tip.
- 6 Spin for 18 h overnight in the ultracentrifuge at 70,000 rpm at 4°C as above.
- 7 Using a rinsed 100 μ l Hamilton syringe remove the bottom 60 μ l to a new Eppendorf tube labeled HDL. Discard the upper portion of the sample containing mostly LDL.
- 8 Measure cholesterol, triglycerides or phospholipids concentrations in the HDL fraction using their respective protocols.
- 9 The non-HDL is calculated by subtracting the HDL from the total.

The density of the HDL fraction is > 1.063 g/ml