

Nov 19, 2022

Untargeted lipidomics analysis for Golgi immunopurification (Golgi-IP)

 Proceedings of the National Academy of Sciences of the United States of America

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4jq6jlo5/v1>

Wentao Dong¹, Eshaan Rawat¹, Monther Abu-Remaileh¹

¹Department of Chemical Engineering, Department of Genetics, The Institute for Chemistry, Engineering & Medicine for Human Health (ChEM-H), Stanford University, Stanford, CA 94305, USA.

Monther Abu-Remaileh: monther@stanford.edu

Metabolomics Protocols & Workflows
Tech. support email: bbmisracb@gmail.com



Monther Abu-Remaileh

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4jq6jlo5/v1>

External link: <https://doi.org/10.1073/pnas.2219953120>

Protocol Citation: Wentao Dong, Eshaan Rawat, Monther Abu-Remaileh 2022. Untargeted lipidomics analysis for Golgi immunopurification (Golgi-IP). protocols.io <https://dx.doi.org/10.17504/protocols.io.3byl4jq6jlo5/v1>

Manuscript citation:

Fasimoye R, Dong W, Nirujogi RS, Rawat ES, Iguchi M, Nyame K, Phung TK, Bagnoli E, Prescott AR, Alessi DR, Abu-Remaileh M, Golgi-IP, a tool for multimodal analysis of Golgi molecular content. Proceedings of the National Academy of Sciences of the United States of America 120(20). doi: [10.1073/pnas.2219953120](https://doi.org/10.1073/pnas.2219953120)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 18, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 72965

Keywords: Golgi, immunoprecipitation, metabolomics, lipidomics, ASAPCRN, golgi metabolite profile, golgi immunopurification, untargeted lipidomics analysis, lipidomics sample, nonpolar lipid profiling, using liquid chromatography mass spectrometry, liquid chromatography mass spectrometry, golgi apparatus function, golgiip, golgi, using immunoprecipitation, untargeted lipidomics analysis for golgi immunopurification

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: SAP-000463

NIH

Grant ID: DP2-CA271386

Abstract

The Golgi apparatus functions as a central hub in the cell that processes, packages, and distributes proteins. Despite its critical cellular function, there has been challenges to quantitatively assess Golgi metabolite profiles. To overcome this hurdle, we developed a rapid harvesting and purification method using immunoprecipitation (GolgiIP). This protocol provides details for analyzing GolgiIP lipidomics samples using liquid chromatography mass spectrometry (LC-MS) for nonpolar lipid profiling.



Materials

Reagents

- Optima LC/MS water (Fisher, cat. no. W6-4)
- Optima LC/MS acetonitrile (Fisher, cat. no. A955-4)
- Optima LC/MS 2-propanol (Fisher, cat. no. A461-500)
- Ammonium formate
- Formic acid
- EASYIC™

Equipment

- ID-X Orbitrap Tribrid Mass Spectrometer (Thermo Fisher Scientific) with a heated electrospray ionization (HESI) probe
- Ascentis Express C18 150 × 2.1 mm column (Millipore Sigma 53825-U)
- 5 × 2.1 mm guard (Sigma-Aldrich 53500-U)

Safety warnings

 Please refer to Safety Data Sheets (SDS) for health and environmental hazards.

LC/MS lipidomics settings

- 1 Set an ID-X tribrid mass spectrometer (Thermo Fisher Scientific) with a heated electrospray ionization (HESI) probe, for initial nonpolar lipid profiling.

Prepare an Ascentis Express C18 150 × 2.1 mm column (Millipore Sigma 53825-U) coupled with a 5 × 2.1 mm guard (Sigma-Aldrich 53500-U), to carry out C18-based lipid separation prior to mass spectrometry. Use EASYICTM for internal calibration.

- 2 For C18-based lipid separation, use 10 millimolar (mM) ammonium formate and 0.1 % (v/v) formic acid dissolved in 60 % (v/v) LC/MS grade water and 40 % (v/v) LC/MS grade acetonitrile for Buffer A, and 10 millimolar (mM) ammonium formate and 0.1 % (v/v) formic acid dissolved in 90 % (v/v) LC/MS grade 2-propanol and 10 % (v/v) LC/MS grade acetonitrile for Buffer B.

- 3 Set the chromatographic gradient flow rate to 0.26 mL/min.

Use Orbitrap resolution 120,000 for MS1 and 30,000 for MS2, RF lens at 40%, AGC target 4×10^5 for MS1 and 5×10^4 for MS2, and maximum injection time 50 ms for MS1 and 54 ms for MS2. Set positive ion voltage to 3250 V, negative ion voltage to 3000 V, ion transfer tube temperature to 300 °C, and vaporizer temperature to 375 °C. Set sheath gas flow to 40 units, auxiliary gas flow to 10 units, and sweep gas flow to 1 unit.

Operate the mass spectrometer in full-scan mode with data-dependent tandem mass spectrometry (ddMS2) at m/z 250–1500, with cycle time of 1.5 sec, microscans of 1 unit, isolation window of m/z 1, intensity threshold of 1×10^4 , and dynamic exclusion time of 2.5 sec. For HCD fragmentation, use step-wise collision energies of 15%, 25%, and 35%.

- 4 Perform the elution with a gradient of 40 minutes: from 0–1.5 min isocratically elute at 32% B; from 1.5–4 min linearly increase to 45% B; from 4–5 min linearly increase to 52% B; from 5–8 min linearly increase to 58% B; from 8–11 min linearly increase to 66% B; from 11–14 min linearly increase to 70%; from 14–18 min linearly increase to 75%; from 18–21 min linearly increase to 97% B; from 21–35 min hold at 97% B; from 35–35.1 min linearly decrease to 32% B; and from 35.1–40 min hold at 32% min.

Untargeted lipidomics workflow

- 5 LipidSearch and Compound Discoverer (Thermo Fisher Scientific) were used for unbiased differential analysis. Lipid annotation was acquired from LipidSearch with the

precursor tolerance at 5 ppm and product tolerance at 8 ppm.

- 6 The mass list from LipidSearch is then exported and used in Compound Discoverer for improved alignment and quantitation. Mass tolerance, 10 ppm; minimum and maximum precursor mass, 0-5,000 Da; retention time limit, 0.1-30 min; Peak filter signal to noise ratio, 1.5; retention time alignment maximum shift, 1 min; minimum peak intensity, 10,000; compound detection signal to noise ratio, 3. Isotope and adduct settings were kept at default values. Gap filling and background filtering were performed by default settings. The MassList Search was customized with 5 ppm mass tolerance and 1 minute retention time tolerance. Area normalization was performed by constant median after blank exclusion.