

Oct 05, 2023

UPitt TriState SenNet TMC H&E staining

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvo3py7v4o/v1>

Marta Bueno¹

¹Division of Pulmonary and Critical Care Medicine. Department of Medicine. School of Medicine

TriState SenNet

Cellular Senescence Net...



Marta Bueno

Division of Pulmonary and Critical Care Medicine. Department...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo3py7v4o/v1>

Protocol Citation: Marta Bueno 2023. UPitt TriState SenNet TMC H&E staining. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvo3py7v4o/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 04, 2023

Last Modified: October 05, 2023

Protocol Integer ID: 88788

Keywords: histology staining, staining hematoxylin, eosin stain, hematoxylin from red, nuclear detail in cell, hematoxylin, dye, staining, bluing reagent, coloration, depth of coloration, preferred method for histopathologic assessment, histopathologic assessment, nuclei of cell, cytoplasm, color change, eosin, cell, upitt tristate sennet tmc, connective tissue fiber, tissue section, different types of connective tissue fiber, tissue, traditional blue color, red, nuclear detail, dna, nuclei, different shades of pink

Funders Acknowledgements:

TriState SenNET (Lung and Heart) Tissue Map and Atlas consortium - NIA

Grant ID: U54AG075931

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Hematoxylin and Eosin stains are the preferred method for histopathologic assessment of tissue sections.

Hematoxylin is used to illustrate nuclear detail in cells. Depth of coloration is not only related to the amount of DNA in the nuclei but also to the length of time the sample spends in hematoxylin. Eosin is the most commonly used counterstain that distinguishes between the cytoplasm and nuclei of cells. It is typically pink, with different shades of pink for different types of connective tissue fibers. (Note: Bluing reagents, such as Scott's Tap Water, are used to change the hematoxylin from red to the traditional blue color we expect. These slightly basic solutions chemically alter the dye to produce this color change.)

Guidelines

- The water steps act as good pause points if a short break (~15 minutes) is needed.
- If the final water (before the 95% Ethanol) has a slight purple color to it, it can be changed and do a short dip into the clean water.
- Check the Hematoxylin staining by looking at it quickly under a light microscope after the final water step making sure the section does not wash out. If the Hematoxylin is not dark enough, it can be placed back into the Hematoxylin for another 30 seconds, followed by the subsequent water steps. (*Note: do not leave in hematoxylin for too long*)
- Before final Xylene baths + mounting, the quality of the staining can be checked under the light microscope. This should be very brief to prevent the slides from drying out. If the Eosin is not dark enough the slides can be placed back into the first 95% Ethanol before the Eosin/Phloxine step and processed through once more.

Materials

·Coverglass

 Epredia™ Shandon-Mount™ **Fisher Scientific Catalog #1900333**

·Cotton Tipped Applicators

·Paper Towel

·Slide Folders

·Kimwipes

·Tissue-Tek Slide Staining Set

 Eosin Y Solution, Alcoholic, with Phloxine **Merck MilliporeSigma (Sigma-Aldrich) Catalog #HT110316**

 Mayer's Hematoxylin Solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MHS16-500ML**

 Epredia™ Xylene **Fisher Scientific Catalog #22-050-283**

 Sodium hydroxide solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #415413**

 Epredia™ Dehydrant Alcohol **Fisher Scientific Catalog #22-050-106**

·Water

·Whatman filter paper

·1 funnel

Safety warnings

- ❗
 - Use universal safety precautions when handling human samples and personal protective equipment (e.g., face mask with shield, gloves, lab coat or apron).
 - Xylenes and ethanol are flammable. All Xylene procedures need to be performed in a fume hood.
 - Gloves should be worn when performing staining process.



Before start

- Make sure that the level of the reagent in each staining line station is the required one.
 - Verify that all slides are completely dry before the staining starts.
-
1. Prepare the bluing solution by adding 1 mL of 28% Sodium Hydroxide to 400 mL of water.
 2. The Hematoxylin should be gravity filtered through the Whatman filter paper before use.
 3. The Eosin is in a ready-to-use state.



Staining Procedure

- 1 Prepare staining line with reagents in the staining train in the order of use.
- 2 Place slides into a xylene compatible slide rack and start the deparaffination and rehydration steps.

A	B
Xylene (I)	4 min
Xylene (II)	4 min
Xylene (III)	4 min (1*)
100% Ethanol (I)	1 min
100% Ethanol (II)	1 min
95% Ethanol	1 min
70% Ethanol	1 min
Water (I)	5 min
Water (II)	5 min
Water (III)	5 min (2*)

Note 1: Blot excess xylene from rack before transferring into ethanol (don't let the slide dry).*

Note 2: Blot excess water from rack before transferring into hematoxylin (don't let the slide dry).*

- 3 Place the slides in Hematoxylin for 1 minute. Adjust timing by 30 second intervals to lighten or darken the stain.
- 4 Wash slides off with water until the purple color is no longer coming off into the water. Leave them in the final water for 5 minutes.
- 5 Transfer to bluing reagent for a brief dip.
- 6 Wash the slides with 3 changes of water and leave in the final water for 5 minutes.



- 7 Transfer to 95% ethanol for 1 minute.
- 8 Transfer to Eosin-Y for 1 minute.
- 9 Transfer to 70% ethanol for 5 minute.
- 10 Follow the dehydration steps:

	A	B
	95% Ethanol	5 min
	100% Ethanol (I)	5 min
	100% Ethanol (II)	5 min (*)
	Xylene (I)	5 min
	Xylene (II)	5 min

*Note *: Blot excess ethanol from rack before transferring into xylene (don't let the slide dry).*

- 11 Coverslip the slides by removing them from the final xylene (1 at a time) and dabbing the edges with paper towel to absorb some excess Xylene. Do not allow the section to dry out.
Place 1-2 drops of mounting media on the tissue section and carefully place the coverslip starting on at an angle and laying over the tissue to avoid air bubbles.
Carefully use the cotton tipped applicator to push out any excess mounting media and air bubbles.
- 12 Place the slides, tissue side up, in a slide folder to dry overnight.

Note: A Kimwipe tissue can be used to wipe off the applied Xylene and cleaning any residual debris from the slide.

Imaging



- 13 Slides will then be scanned by the Center for Biologic Imaging (CBI) using a VS200 slide scanner. Images can be acquired with brightfield scans at (10X, 20X).