

Aug 08, 2024

URM TriState SenNet TMC Mouse Sacrifice and Organ Harvest

DOI

dx.doi.org/10.17504/protocols.io.rm7vzj2j8lx1/v1

Gagandeep Kaur¹, Irfan Rahman¹

¹Department of Environmental Medicine, University of Rochester Medical Center, Rochester, NY

TriState SenNet



Gagandeep Kaur

University of Rochester Medical Center

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzj2j8lx1/v1>

Protocol Citation: Gagandeep Kaur, Irfan Rahman 2024. URM TriState SenNet TMC Mouse Sacrifice and Organ Harvest. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzj2j8lx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 08, 2024

Last Modified: August 08, 2024

Protocol Integer ID: 105002

Keywords: Mouse handling, Sacrifice, Lung harvest, urmc tristate sennet tmc mouse sacrifice, tristate sennet tmc biospecimen core, cellular senescence, part of the cellular senescence network program, mouse lung tissue, organ harvest this protocol, cellular senescence network program, organ harvest, sennet, tissue, mice

Funders Acknowledgements:

TriState SenNet

Grant ID: Grant ID: U54 AG075931

Disclaimer

The authors have no conflict of interest to declare.

Abstract

This protocol outlines the method used to sacrifice and harvest the mouse lung tissue to perform downstream analyses for the TriState SenNet TMC Biospecimen Core at the University of Rochester, as part of the Cellular Senescence Network Program (SenNet).


Care was taken to reduce the pain and distress caused to the mice by performing all the procedures after anesthetising the mice.

Materials

1. Sodium Pentobarbital (Cat# NDC 25021-676-20, Sagent)
2. Insulin Syringes
3. Saline
4. 3mL Syringes
5. Eppendorf tubes
6. Scissors
7. 1X PBS

Troubleshooting



- 1 Sodium Pentobarbital (60-90 mg/kg of body weight) was administered intraperitoneally to anaesthetize the mice prior to organ harvest.
- 2 After 5 min of drug administration, the animal response was assessed by pinching the mouse's foot with forceps and observing the withdrawal response/pedal reflex.
- 3 Upon successful anesthetization, the animals abdomen and chest was cut open and mice was exsanguinated by cardiac puncture.
- 4 Thereafter, the thoracic cavity was cut open to harvest the mouse lung.
- 5 The harvested organ was stored in PBS and kept  On ice to perform future assays/assessments.