

Dec 04, 2024

USDA LTAR Common Experiment measurement: Dissolved nitrate (NO_3^-) concentration

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2qm1pl1y/v1>



Richard Lizotte¹, Oliva Pisani², Kristen S. Veum³, John L. Kovar⁴, Robert W. Malone⁴, Kevin J. Cole⁴

¹USDA Agricultural Research Service, Water Quality and Ecology Research, Oxford, MS;

²USDA Agricultural Research Service, Southeast Watershed Research Laboratory, Tifton, GA;

³USDA Agricultural Research Service, Cropping Systems and Water Quality Research, Columbia, MO;

⁴USDA Agricultural Research Service, National Laboratory for Agriculture and the Environment, Ames, IA



Lori J. Abendroth

USDA ARS Cropping Systems and Water Quality Research Unit

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2qm1pl1y/v1>

External link: <https://ltar.ars.usda.gov>



Protocol Citation: Richard Lizotte, Oliva Pisani, Kristen S. Veum, John L. Kovar, Robert W. Malone, Kevin J. Cole 2024. USDA LTAR Common Experiment measurement: Dissolved nitrate (NO₃⁻) concentration. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l2qm1pl1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 12, 2024

Last Modified: December 04, 2024

Protocol Integer ID: 97167

Keywords: Long-Term Agroecosystem Research, LTAR, USDA LTAR, Common Experiment, crops, nitrate, surface water, groundwater, drinking water standard, nitrogen, nitrate nitrogen, colorimetry, concentration nitrate, nitrate reductase, cadmium with nitrate reductase, nitrogen, ionized form of nitrogen, colorimetric cadmium reduction method, nitrite, mg no₃, nutrient, hypoxia of surface water, photosynthesis in aquatic autotroph, dissolved organic matter, essential nutrient for growth, presence of cadmium, essential nutrient, cadmium reduction, contributes to eutrophication, high levels in some surface water, colored azo dye, flow analyzer method with cadmium reduction, drinking water standard, identical colorimetric quantification reaction, cadmium method, copperized cadmium, groundwater, cadmium, photosynthesis,

Funders Acknowledgements:

United States Department of Agriculture

Grant ID: -

Abstract

Nitrate (NO_3^-) is an ionized form of nitrogen (N) that is of interest in surface and groundwater. Nitrate often occurs in trace amounts (<0.01 mg/L) in rivers and lakes but can reach high levels in some surface water and groundwater. In excessive amounts, it can elicit an illness known as methemoglobinemia in infants, and the drinking water standard is limited to 10 mg NO_3^- -N/L to prevent this illness. In addition, NO_3^- is an essential nutrient for growth and photosynthesis in aquatic autotrophs, and elevated NO_3^- loading contributes to eutrophication and hypoxia of surface waters, particularly in coastal marine ecosystems. The recommended methods for measuring NO_3^- concentration in water include the colorimetric cadmium reduction method, whereby NO_3^- is reduced quantitatively to nitrite (NO_2^-) in the presence of cadmium. Following the addition of appropriate reagents, a highly colored azo dye is produced and measured colorimetrically. An alternative enzymatic method utilizes “greener” chemistry by replacing the toxic copperized cadmium with nitrate reductase (specifically, AtNaR2) in combination with two coenzymes (nicotinamide adenine dinucleotide and nicotinamide adenine dinucleotide phosphate) to reduce NO_3^- to NO_2^- , followed by the identical colorimetric quantification reaction used in the cadmium method. This enzyme reduction method overcomes interferences from phenolics and humic substances noted in earlier enzymatic methods. The method detection levels (or limits; MDL) for standard and low $\text{NO}_3^- + \text{NO}_2^-$ concentrations with enzyme reduction are comparable to the long-term MDLs of the continuous-flow analyzer method with cadmium reduction. In addition, NO_3^- concentration differences were statistically equivalent to zero or analytically equivalent to zero (i.e., concentration differences were less than MDLs [Patton and Kryskalla, 2011]). Recent advances in field sensor technology have afforded commercial sensors that allow for in situ NO_3^- measurement using an optical sensor. These sensors can record values over short durations (seconds to minutes) and use an ultraviolet spectrophotometer (Johnson and Coletti, 2002). Although these sensors provide high-resolution data, their accuracy can degrade because of background interference from colored dissolved organic matter (DOM) and high turbidity. Thus, sensor readings must be validated at least monthly by conducting water sample collection and lab analyses.

Safety warnings

-  Follow appropriate safety, health, and environmental precautions based on the selected methods, instrumentation, and workflow. Laboratory supervisors are responsible for knowledge of these precautions and their implementation.



Sample collection and filtration

- 1 Return samples to the laboratory  On ice and filter them on collection day if possible.
- 2 For NO_3^- measurements, filter water samples through a  0.45 μm pore-size filter to minimize interference from particulates before chemical analysis. If desired, remove the microbial community using a  0.22 μm filter. Filters can be membrane or glass fiber, and may require pre-washing to prevent contamination.
- 3 Occasional checks of filter blanks by filtering deionized water with the same equipment are always prudent. The use of a blank is described in the *USDA LTAR Common Experiment measurement: Best practices for collection, handling, and analyses of water quality samples* protocol (Pisani et al., 2024b).

Sample storage and preservation

- 4 Analyze samples as soon as possible after collection.
- 5 If storage is required, maintain filtered samples at  4 °C (not frozen) for up to 48 hours (US EPA, 1993). If lengthier storage is required, prepare samples with concentrated sulfuric acid to a pH < 2, and maintain them at  4 °C .
- 6 Hold acidified samples for up to 28 days in clean (N-free) plastic bottles. For acidified samples, NO_3^- and NO_2^- cannot be determined as individual species but only as total inorganic N.

Archiving

- 7 The common practice is to store water samples for NO_3^- analysis at  4 °C (not frozen) until data certification (QA/QC verification).
- 8 Storage of preserved water samples should not exceed 28 days.

Sample analysis

- 9 Using the cadmium reduction method, determine NO_3^- in water by passing a filtered sample through a column containing granulated copper–cadmium to reduce NO_3^- to NO_2^- .
- 10 The sum of NO_2^- initially present plus the reduced NO_3^- is measured colorimetrically (APHA, 2005; USEPA, 1993).
- 11 The sum of $\text{NO}_3^- + \text{NO}_2^-$ is sometimes called NO_x^- .
- 12 To obtain separate NO_3^- and NO_2^- values, perform the procedure with the reduction step (total NO_x^-) and then without the reduction step (NO_2^- only), and calculate NO_3^- by difference.
- 13 The NO_2^- is usually, but not always, a relatively small component of the total.
- 14 To reduce NO_3^- to NO_2^- using the enzymatic reduction method, nitrate reductase (specifically, AtNaR2) is used in combination with two coenzymes (nicotinamide adenine dinucleotide and nicotinamide adenine dinucleotide phosphate). Enzymatic reduction is an alternate method to cadmium reduction.
- 15 This effort is followed by the identical colorimetric quantification reaction used in the cadmium method (ASTM, 2014; Patton and Kryskalla, 2011), as described in Step 10.

Covariate metrics to be sampled concurrently

- 16 Dissolved ammonia (NH_3), total dissolved N (TDN), total dissolved phosphorus (TDP), and total suspended solids (TSS). These metrics are described in the following protocols:
 - USDA LTAR Common Experiment measurement: Dissolved ammonia (NH_3) concentration
 - USDA LTAR Common Experiment measurement: Total nitrogen (TN) and total dissolved nitrogen (TDN) concentration
 - USDA LTAR Common Experiment measurement: Total phosphorus (TP) and total dissolved phosphorus (TDP) concentration
 - USDA LTAR Common Experiment measurement: Total suspended solids (TSS)

Calculations

- 17
- Prepare a calibration curve by plotting the blank-corrected area of each standard peak against its respective NO_3^- or NO_x^- ($\text{NO}_3^- + \text{NO}_2^-$) concentration. The laboratory blank is used for this step.
- 18
- Compute the NO_3^- or NO_x^- concentration of the sample by comparing the sample peak area to the calibration curve. Report the concentration(s) as $\text{mg NO}_x^- \text{-N/L}$.

Recommendations for data collection

- 19
- Table 1. Summary of recommendations for the collection and measurement of dissolved NO_3^- concentration.

A	B	C	D
Attribute	Preferred	Minimum	Comments
Spatial scale	Field	Plot	
Frequency	Event-driven	Event-driven	More frequent (weekly) measurements may be preferred when the flow regime varies seasonally or following precipitation events. Sampling in this protocol is event-driven to enable cross-site comparisons
Covariate metrics	NH3-N, TDN, TDP, TSS	NH3-N, TDN, TDP	
Sample preservation and storage	Filter with a 0.45 μm pore-size filter (membrane or glass fiber), refrigerate as soon as possible, and analyze within 48 hours	Filter with a 0.45 μm pore-size filter (membrane or glass fiber), freeze as soon as possible, and analyze within 28 days	Alternatively, use a 0.7 μm GF/F. Use a 0.2 μm pore-size filter to exclude bacteria from samples but at much slower flow rates
Sample analysis	Nitrate reductase with colorimetric analysis	Cadmium reduction with colorimetric analysis	Manually set up reduction columns or use them as part of an auto-analyzer system
Water quantity	Discharge or flow rate	Discharge or flow rate	Calculate $\text{NO}_3\text{-N}$ loads by linking this metric to the

	A	B	C	D
				water quantity metric "flow"

Covariate metrics = NH₃-N = ammonia-N; TDN = total dissolved N; TDP = total dissolved phosphorus; TSS = total suspended solids.

Protocol references

American Public Health Association (APHA), (2005). Standard Methods for the Examination of Water and Wastewater, 21st ed. Washington, DC: American Public Health Association, American Water Works Association, and Water Environment Federation.

American Society for Testing and Materials (ASTM) International, (2014). Standard Test Method for Nitrite-Nitrate in Water by Nitrate Reductase. Standard ASTM-D7781-14.

Dalzell, B. & Pisani, O. (2024). USDA LTAR Common Experiment measurement: Total suspended solids (TSS). protocols.io dx.doi.org/10.17504/protocols.io.261ge5pjog47/v2

Hamilton, S. K., Pisani, O., Kovar, J. L., Malone, R. W., Morrow, A. J., & Cole, K. J. (2024). USDA LTAR Common Experiment measurement: Total phosphorus (TP) and total dissolved phosphorus (TDP) concentration. protocols.io dx.doi.org/10.17504/protocols.io.8epv5r7m6g1b/v1

Johnson, K. S. & Coletti, L. J. (2002). In situ ultraviolet spectrophotometry for high resolution and long-term monitoring of nitrate, bromide and bisulfide in the ocean. Deep Sea Research Part I: Oceanographic Research Papers, Vol 49, Issue 7, Pages 1291-1305.

Malone, R. W., Morrow, A. J., Pisani, O., Kovar, J. L., Hamilton, S. K., & Cole, K. J. (2024). USDA LTAR Common Experiment measurement: Total nitrogen (TN) and total dissolved nitrogen (TDN) concentration. protocols.io dx.doi.org/10.17504/protocols.io.5jyl82rkrl2w/v1

Patton, C. J. & Kryskalla, J. R. (2011). Colorimetric determination of nitrate plus nitrite in water by enzymatic reduction, automated discrete analyzer methods: U.S. Geological Survey Techniques and Methods, Book 5, Chapter B8, Page 34.

Pisani, O., Kovar, J. L., Malone, R. W., Morrow, A. J., & Cole, K. J. (2024a). USDA LTAR Common Experiment measurement: Dissolved ammonia (NH₃) concentration. protocols.io dx.doi.org/10.17504/protocols.io.j8nlk8b6l15r/v1

Pisani, O., Lizotte, R., Veum, K. S., Kovar, J. L., Hamilton, S. K., & Malone, R. W. (2024b). USDA LTAR Common Experiment measurement: Best practices for collection, handling, and analyses of water quality samples. protocols.io dx.doi.org/10.17504/protocols.io.q26g71z68gwz/v1

US Environmental Protection Agency (US EPA), (1993). Method 353.2, Revision 2.0: Determination of Nitrate-Nitrite Nitrogen by Automated Colorimetry.



Acknowledgements

This research is a contribution from the Long-Term Agroecosystem Research (LTAR) network. LTAR is supported by the United States Department of Agriculture. The use of trade, firm, or corporation names in this publication is for the information and convenience of the reader. Such use does not constitute an official endorsement or approval by the United States Department of Agriculture or the Agricultural Research Service of any product or service to the exclusion of others that may be suitable. USDA is an equal opportunity provider and employer.