

Oct 18, 2021

UW Virology Swift SNAPv2 protocol

This protocol is a draft, published without a DOI.

Lasata Shrestha¹, Hong Xie¹, Shah A. Mohamed Bakhsh¹, Robert J. Livingston¹, Meei-Li Huang¹, Alexander L. Greninger¹, Pavitra Roychoudhury¹

¹Department of Laboratory Medicine and Pathology, University of Washington, Seattle, WA



Lasata Shrestha

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



Protocol Citation: Lasata Shrestha, Hong Xie, Shah A. Mohamed Bakhsh, Robert J. Livingston, Meei-Li Huang, Alexander L. Greninger, Pavitra Roychoudhury 2021. UW Virology Swift SNAPv2 protocol. **protocols.io** <https://protocols.io/view/uw-virology-swift-snapv2-protocol-byw4pxgw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 09, 2021

Last Modified: October 18, 2021

Protocol Integer ID: 53948



Keywords: whole genome sequencing, SARS-CoV-2, SWIFT Biosciences, PerkinElmer, automation, genome, protocol viral whole genome sequencing, multiplex amplicon method, amplicon, amplifies 345 amplicon, complete genome, multiplex amplicon panel, instrumental in outbreak investigation, outbreak investigation, vaccine, evaluation of vaccine

Abstract

Viral whole genome sequencing (WGS) has been instrumental in outbreak investigations, deployment of public health interventions, development as well as evaluation of vaccines and therapeutics. While multiple methods are commercially available for WGS, multiplex amplicon method has proven to be faster, more efficient, scalable, and more cost-effective compared to other methods. Here, we describe the automation of a multiplex amplicon panel for WGS of Severe Acute Respiratory Syndrome Coronavirus-2 (SARS-CoV-2), the causative agent of the current COVID-19 pandemic. The SWIFT Biosciences' primer set amplifies 345 amplicons designed against the SARS-CoV-2 Wuhan-Hu-1 complete genome (NC_045512.2), in a single tube to cover the ~30 kb SARS-CoV-2 genome in less than three hours.



Materials

- SuperScript IV Kit **[Kit content](#)**
- RNaseOUT™ Recombinant Ribonuclease Inhibitor **[Kit details](#)**
- SNAP SARS-CoV-2 Kit **[Kit content Page: 4](#)**
- SPRIselect beads <https://www.beckman.com/reagents/genomic/cleanup-and-size-selection/size-selection>
- 200 Proof Ethanol
- Molecular-grade water
- Qubit®, Nanodrop, or other similar input RNA quantification assay
- qPCR-, electrophoretic-, or fluorometric-based library quantification assay for Illumina® libraries
- Microcentrifuge
- Vortex
- Programmable thermocycler
- Aerosol-resistant tips and pipettes ranging from 1 to 1000 µL
- Pipette tips (e.g., 8-channel or 12-channel), 8-tube strips, an un-skirted 96 well plate, or plate puncher for pre-piercing the foil seal if using single-use UD indexing plates.
- Qubit dsDNA HS Assay Kit <https://www.thermofisher.com/order/catalog/product/Q32851#/Q32851>
- TapeStation DNA ScreenTape & Reagents <https://www.agilent.com/en/product/automated-electrophoresis/tapestation-systems/tapestation-dna-screentape-reagents/dna-screentape-analysis-228260>
- 80 uL - Barrier Sterile 96 Rack Tips <https://www.perkinelmer.com/product/80ul-art-sterile-96-rack-10-racks-111624>
- 150 uL - Barrier Sterile 96 Rack Tips <https://www.perkinelmer.com/product/150ul-96-art-tip-box-10-racks-111426>
- Polypropylene 384-well microplate, U-Bottom, case of 50 <https://www.perkinelmer.com/product/pp-microplate-384-35ul-u-bottom-50-6008890>
- Polypropylene 96-well Microplate, Deep Well U-bottom, 2 mL <https://www.perkinelmer.com/product/pp-microplate-96-2ml-v-bottom-25-6008880>
- HARDSHELL PCR PLATE-96, BLUE/ 50 <https://www.perkinelmer.com/product/hardshell-pcr-plate-96-blue-50-6008870>
- StorPlate-96V, PP, 96 well, V-bottom, (V), 450µL, 200/box <https://www.perkinelmer.com/product/storplate-96-v-450-l-200-6008299>
- PP RESERVOIR, DW, V, 12 COL, 21mL /25 <https://www.perkinelmer.com/product/pp-reservoir-dw-v-12-col-21ml-25-6008700>

Before start

The critical steps are tagged. Steps 5-108 are optional (for troubleshooting).

Single-strand cDNA synthesis Reagent Set-up

- 1 Prepare reagents for required number of samples as per the manufacturer's instructions (including 10% overage) linked below.

[https://assets.thermofisher.com/TFS-](https://assets.thermofisher.com/TFS-Assets/LSG/manuals/SSIV_Reverse_Transcriptase_UG.pdf)

[Assets/LSG/manuals/SSIV_Reverse_Transcriptase_UG.pdf](https://assets.thermofisher.com/TFS-Assets/LSG/manuals/SSIV_Reverse_Transcriptase_UG.pdf)



Single-strand cDNA (sscDNA) synthesis Reagent Sciclone Set-up

- 2 Select the correct protocol on the liquid handler (Perkin Elmer Sciclone) and load the Sciclone as instructed.



This program will perform steps corresponding to steps 3 to 20 in the manual protocol. The code can be obtained from PerkinElmer.

Note

The SuperScript RT IV protocol was automated and optimized by Perkin Elmer personnel following ThermoFisher's commercially available protocol linked below.

SuperScript RT IV




Fig. 01: Empty Sciclone deck.

Equipment

Sciclone G3 NGSx iQ Workstation

NAME

Automated Liquid Handling

TYPE

Perkin Elmer

BRAND

CLS145321

SKU

<https://www.perkinelmer.com/category/automation-liquid-handling-instruments>

LINK








Fig. 02 Final (fully set-up) Sciclone desk set-up for sscDNA synthesis.

Sciclone movement: sscDNA synthesis


- 3 Moved Plate from D4 to B4 (magnet with no spacer)





- 4 Moved Plate from D2 to D4
- 5 Aspirate  11 μL of RNA from B4 and Dispense to plate on D4
- 6 Pipette mix and shake for  00:01:00 1m
- 6.1 Plate on B4 moved to C0 concurrently
- 6.2 Plate on C0 discarded
- 7 Moved plate from D4 to onboard thermocycler
- 8 Move ODTc lid from A1 to thermocycler
- 9 Thermocycler door closes and runs  00:07:40 7m 40s
- 10 Thermocycler door opens
- 11 ODTc lid moved from Thermocycler to A1
- 12 Moved plate from Thermocycler to D4
- 12.1 Empty tip box discarded. Fresh tip box placed onto C3
- 13 Load new tips
- 14 Aspirate  7 μL of PCR master mix from A4 and dispense into plate on D4




14.1 Pipette mix and shake for  00:01:00

1m

15 Move plate from D4 to onboard thermocycler

16 Move ODT lid from A1 to thermocycler

17 Thermocycler door closes and runs for  00:53:04

53m 4s

18 Thermocycler door opens

19 Plate moved from thermocycler to D4

20 End product of  20 μL

Swift SNAP v2 Reagent Set-up

21 Prepare reagents for required number of samples as per the manufacturer's instructions (including 10% overage) linked below.



Swift Normalase™ Amplicon Panels (SNAP) SARS-CoV-2, Additional Genome Coverage, and SARS-CoV-2 S Gene Panel

Swift SNAP v2 Sciclone Set-up

22 Select the correct protocol on the liquid handler (Perkin Elmer Sciclone) and load the Sciclone as instructed.



The selected program will perform steps corresponding to steps 23 to 127. The code can be obtained from PerkinElmer.

Note

The SWIFT SNAP v2 protocol was automated and optimized by Perkin Elmer personnel following SWIFT Bioscience's commercially available protocol linked below.

SWIFT SNAP v2

Insert labeled photo of the deck without components (deck only)




Fig. 03. Empty Sciclone deck.

Insert labeled photo of the deck with labeled components and volume information




Fig. 04. Final (fully-loaded) deck set-up for SWIFT SNAP v2 protocol on the Sciclone.

Equipment	
Sciclone G3 NGSx iQ Workstation	NAME
Automated Library Preparation	TYPE
Perkin Elmer	BRAND
*****	SKU
https://www.perkinelmer.com/category/automation-liquid-handling-instruments	LINK

Sciclone Movements: Swift SNAP v2 Multiplex PCR

1h 20m



Note

Steps 23-127 include detailed Sciclone deck movements for troubleshooting.

- 24 Moved Lid from A4 to A2
- 25 Pipette multiplex master-mix (A4 to samples on D4)
- 26 Dispose off tips
- 27 Lid from A2 to A4
- 28 Mix samples on D4
- 29 Tips disposed
- 30 Shake mix of D4 (Noted: 30uL)
- 31 D4 to Thermocycler
- 32 ODTC lid placed
- 33 Thermo-cycler cover closed ⌚ 01:15:00
- 34 PCR door opened
- 35 ODTC lid removed


1h 15m



36 The PCR plate was moved to D4

37 Moved plate from B2 to A3


38 Mix beads in B2, 10 X

39 Picked up 30 uL Beads from B2 to D4  30 μ L Beads

40 Mixed D4 wells, 10 X

41 Shake mix D4 for several minutes

42 Plate moved from D4 to B4

43 Allow beads to sediment  00:05:00

5m

44 Supernatant discarded

45 Tips discarded



46 Alcohol cover B5 removed to C5

47 Fresh tips picked.

48 Picked up 80%ETOH  150 μ L 80% Ethanol

49 Dispensed to D2



- 50 Supernatant removed to Wash plate
- 51 Tips disposed
- 52 Fresh tips picked
- 53 Alcohol (150uL) dispensed to plate B4  150 μ L
- 54 Removed supernatant from B4
- 55 Tips discarded
- 56 Fresh tips picked
- 57 Alcohol transferred to B4  150 μ L 80% EtOH
- 58 Supernatant collected and disposed
- 59 Tips discarded
- 60 Alcohol lid replaced back
- 61 Sample plate moved to D4
- 62 New tips
- 63 TE buffer lid



64 TE buffer added  17.4 μ L NEED TO CHECK

65 Mixed

66 Tips discarded

67 Plate D4 shake mixed

68 Lid from A4 lifted

69 Transfer master mix A4 to A3  28.9 μ L Indexing PCR Reaction Mix

70 Tips disposed

Sciclone Movements: SWIFT SNAP v2 Indexing PCR

1h 20m

71 Pipette Index from D2 to A3 (3.7 μ L of indexes)



72 Pipette master mix plus indexes from A3 to D4

73 Shake mixed D4

74 Lid of A2 to A4

75 Fresh tips

76 D4 mixed by pipetting



77	Tips disposal		
78	D4 shake mixed		
79	D4 plate moved to thermocycler		
80	21 min PCR run	00:21:00	21m
81	Incubation 26 min	00:26:00	26m
82	Thermocycler door open		
83	OPTC lid removed		
84	Sample plate moved to D4		
85	Beads plate moved from B2 to D3		
86	Pipetted PEG NaCl from plate on B2 to plate on D4	32.5 μ L PEG NaCl (ratio: 0.65)	
87	Mixing in D4 by pipetting (10 times)		
88	Shake mixed D4, for 4 min 35 sec	00:04:35	4m 35s
89	Plate moved from D4 to C4 (on the magnet)		
90	Incubation at room temperature	00:05:00	5m




91 Supernatant from C4 discarded

92 Tips discarded

93 Alcohol lid removed

94 New tips picked

95 Picked up alcohol  150 μ L 80% EtOH


96 Transferred alcohol to C4

97 1 min wait

98 Supernatant discarded

99 Tips dropped off

100 New tips picked up



101 Picked up alcohol  150 μ L 80% EtOH

102 Dispensed to C4


103 1 min wait

104 Supernatant discarded



- 105 Tips discarded
- 106 Moved Tip container from C3 to D0
- 107 New tip container picked up from A0
- 108 New tips picked up
- 109 Alcohol picked up  150 μ L 80% EtOH
- 110 Dispensed alcohol to C4
- 111 1 min wait
- 112 Supernatant from C4 pipetted up and discarded
- 113 Tips discarded
- 114 Alcohol cover replaced
- 115 2 min wait (incubate at room temperature until residual alcohol evaporates)
- 116 Moved Plate from C4 to D4
- 117 Picked up TE
- 118 Transfer TE to D4  20 μ L TE



- 119 Mixing
- 120 Tips discarded
- 121 Mix by pipetting and shaking in plate D4
- 122 Moved plate D4 to B4
- 123 Moved empty Plate A2 to D4
- 124 Tips picked up
- 125 Transfer from B4 to D4  20 μ L Eluate
- 126 Tips discarded
- 127 Plate lid from A4 to D4

Library Quality Control I

- 128
Following the manufacturer's instructions for library quality control, quantify (Qubit or other fluorometric instruments) and determine the size (TapeStation or other electrophoretic instruments) of the library.





Equipment

4200 TapeStation System

NAME

Electrophoresis tool for DNA and RNA sample quality control.

TYPE

TapeStation Instruments

BRAND

G2991AA

SKU

<https://www.agilent.com/en/product/automated-electrophoresis/tapestation-systems/tapestation-instruments/4200-tapestation-system-228263>

LINK

Equipment

Invitrogen™ Qubit™ 3 Fluorometer

NAME

Accurately measures DNA, RNA, and protein using the highly sensitive fluorescence-based Qubit quantitation assays

TYPE

Invitrogen™ Q33216

BRAND

Q33216

SKU

<https://www.fishersci.co.uk/shop/products/qubit-3-0-fluorometer/15387293>

LINK

Normalase

- 129 Following manufacturer's instructions as linked below, normalize the libraries using SWIFT's proprietary enzymatic normalization.

Normalase Pages: 9-12



Pooled Library Quality Control II



- 130 Determine the library concentration for final library quality control before loading on the sequencer.



Quality Control Page: 12

Equipment

Invitrogen™ Qubit™ 3 Fluorometer

NAME

Accurately measures DNA, RNA, and protein using the highly sensitive fluorescence-based Qubit quantitation assays

TYPE

Invitrogen™ Q33216

BRAND

Q33216

SKU

<https://www.fishersci.co.uk/shop/products/qubit-3-0-fluorometer/15387293>

LINK

Sequencing

- 131 Following manufacturer's recommendation, load the pooled library for sequencing on the selected/preferred sequencer.



Sequencing Page: 13

Equipment

NextSeq 500 System

NAME

Sequencer

TYPE

Illumina

BRAND

SKU

<https://www.illumina.com/systems/sequencing-platforms/nextseq.html>

LINK