



Oct 02, 2024

VAMPseq Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2q77xl1y/v1>

Raining Wang¹, Melinda Wheelock¹

¹University of Washington

IGVF



Sarah Heidl

University of Washington, Starita Lab

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2q77xl1y/v1>

Collection Citation: Raining Wang, Melinda Wheelock 2024. VAMPseq Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l2q77xl1y/v1>

License: This is an open access collection distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this collection and it's working

Created: May 10, 2024

Last Modified: October 02, 2024

Collection Integer ID: 99791

Keywords: MAVE, functional assay, variant abundance by massively parallel sequencing, missense, variation, amino acid, vampseq protocol, vampseq protocol this protocol, egfp reporter system, sequencing, egfp, experiment of variant abundance, variant abundance, gene, vamp, designation plasmid, assay, seq

Funders Acknowledgements:

NHGRI (IGVF)

Grant ID: HG011969

Abstract

This protocol describes the materials and method to conduct an experiment of Variant Abundance by Massively Parallel sequencing (VAMP-seq) assay. There are two designation plasmids available for a gene of interest (GOI) to fuse with a EGFP reporter system fused either N-terminally (N-terminal; EGFP-GOI) or C-terminally (C-terminal; GOI-EGFP).

Guidelines

Pre-barcoded Vector Design:

The current (as of 2023) working system (VAMP-seq v2) is modified based on the published version (Matreyek et al, 2018). There are two major modifications made for a two-step “plug and play” system compatible with Golden Gate Assembly. First, a pair of BsmBi/ Esp3I recognition sequences were added before the T7 promoter sequence in both of the N- & C- terminus systems. This is used to add barcode sequences to the vector. Secondly, a pair of SapI recognition sequences were added to assist cloning of the GOI variant library. In the N-terminal vector, the SapI sites are located between the linker sequences and IRES. For the C-terminal vector, the SapI sites are located between the Kozake site and the linker sequence of EGFP.

Barcodes and GOI Design

Barcode sequence design:

Below is the final amplicon sequence for adding barcodes to either N- & C- terminal plasmids.

5'-
GGCTAC**CGTCTC***Caggc****JTAAGTTGGGTAACGCCAGG****ACAACCGGTTAGAGCTCGTTTATGATACTA*
*GTATCGGCTAGCAGGACAGTTCTGCAATTGCGTGAGTAGGCAAGAACCGCTAGAAG****CGTCGCTGTA***
*CAAAATAGTTNNNNNNNNNNNNNNNNNNNN**acga****TGAGACG****GTAGCC*-3'

The Italicized template sequence is a double-stranded DNA fragment. BsmBi/Esp3I sequences (in red) are added to both ends of the template by a primer set (in bold) using PCR amplification. The lowercase indicates fragment-specific overhangs for backbone assembly. In particular, a degenerative reverse primer has designed a 18nt bases (N) barcode. The lowercase indicates fragment-specific overhangs for backbone assembly. The underlined sequence is important for downstream genomic sequencing.

Gene of Interest flanking design:

VAMP-seq v2 vectors have two pre-existing SapI sites and pre-selected fragment specific overhangs ready to receive GOI sequence with complimentary flanking ends. The specific flanking sequences are as below:

In the N-terminus system, 5'-atg and 5'-CTA (reverse complement of 5'-tag") are the pre-selected fragment-specific overhangs in the backbone. The insert fragment must contain the SapI recognition sequences (in red) and overhangs (indicated as lowercase) as the following:

5'-AAGACC**GCTCTTC***Catg*NNNNNNNNNNNNNNNNNNNN*tag***GAAGAGC**GGTCTT-3'



In the C-terminus system: 5'-atg and 5'-CGT (reverse complement of 5'-acg) are the pre-selected fragment-specific overhangs in the backbone to receive a compatible insert fragment. The insert fragment must contain the SapI recognition sequences (in red) and overhangs (indicated as lowercase) as the following:

AAGACC**GCTCTTC**CatgNNNNNNNNNNNNNNNacgG**GAAGAGC**GGTCTT-3'

*A stop codon should not be placed at the end of GOI sequence when assembling to the C-terminus backbone. It is already placed at the end of the EGFP. Albeit, in the N-terminus system Albeit, the stop codon is necessary for the GOI insert. It is a part of the SapI overhang.

Both Start and Stop codons cannot be mutated.

Wild-type and destabilizing variants design:

Wild-type of Gene of Interest:

Wild type sequences for gene of interest(s) can be obtained from PCR amplification of existing template, or synthesized as double-stranded DNA fragments. In both cases, the interior sequence should be free of SapI recognition sequences and flanked by the above sequences suitable for the desired designation vectors.

Destabilizing Variants Selection:

Variant Library design:

NGS primer/region design

Files

 SEARCH

Protocol

NAME

Cloning with Golden Gate

VERSION 1

CREATED BY



Sarah Heidl
University of Washington, Starita Lab

OPEN →

Protocol

NAME

Barcoded vector cloning

VERSION 1

CREATED BY



Sarah Heidl
University of Washington, Starita Lab

OPEN →

Protocol

NAME

Library cloning

VERSION 1

CREATED BY



Sarah Heidl
University of Washington, Starita Lab

OPEN →

Protocol

NAME

Tissue culture- Purity sorting HEK293T LLP iCasp9 cells

VERSION 1

CREATED BY



Sarah Heidl
University of Washington, Starita Lab

OPEN →

Protocol

NAME

Library transfection

VERSION 1

CREATED BY



Sarah Heidl
University of Washington, Starita Lab

OPEN →

Protocol

NAME

Genomic DNA Extraction from Sorted Cells

VERSION 1

CREATED BY



Sarah Heidl
University of Washington, Starita Lab

OPEN →

Protocol

NAME

PCR amplification of the barcode region

VERSION 1

CREATED BY



Sarah Heidl
University of Washington, Starita Lab

OPEN →