

Feb 02, 2025

Version 3

Vector Digestion and Purification V.3

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge5dyog47/v3>

Carolina Lopez¹

¹Washington University



Cecilia Escudero

wustl

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge5dyog47/v3>

Protocol Citation: Carolina Lopez 2025. Vector Digestion and Purification. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ge5dyog47/v3> Version created by **Carolina Lopez**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 14, 2025



Last Modified: February 02, 2025

Protocol Integer ID: 118321

Keywords: purification protocol for plasmid digestion, plasmid digestion, purification protocol, vector digestion, purification, plasmid

Abstract

Protocol for plasmid digestion and purification

Materials

Reagents:

- Restriction Enzymes (New England Biolabs)
- 10x Cutsmart Buffer (New England Biolabs)
- Agarose
- EtBr
- NEB Monarch DNA Gel Extraction kit (NEB #T1020)



Digestion:

1h 10m

- 1 1. Double digest Vector with New England Biolabs Restriction Enzymes: 80uL reaction

1h 10m

A	B
Component	Volume (uL)
DNA plasmid	X ul for 4 ug of plasmid
10x Cutsmart Buffer	8 uL
Enzyme 1	4 uL
Enzyme 2	4 uL
NucleaseFree H2O	64-X uL

2. Incubate for 00:30:00 at 37 °C .
3. Add 16 µL of 6x loading Dye to reaction and vortex briefly to mix.
4. Make 1% agarose gel.
- Mix 1 g of Agar with 100 mL of TAE Buffer.
 - Microwave to boil agarose and let cool until you can touch bottle, but gel is not solid.
 - Add 2 µL of EtBr to agarose and pour into DNA gel mold with 10 well comb.
 - Let gel solidify.
5. Load 100 µL of reaction into well of gel
6. Run gel for 00:40:00 at 120V.
7. Visualize band with UV light and cut out section of gel with band with new razor blade and place in 1.5 mL tube.

Notes:

- Not all two enzymes can be used together in a double digestion. Think about it ahead and you can test your two enzymes here: <https://nebcloner.neb.com/#!/redigest>.
- Enzymes are different, but most can be digested in CutSmart buffer at 37°C within 30 minutes. This is a general protocol. Please check the best digestion conditions for your enzymes in this web : <https://nebcloner.neb.com/#!/redigest>. It will provide the temperature, buffer, and incubation time specific to the two enzymes you are using.

Purification bands from gel with Monarch DNA Gel Extraction kit:

- 2 1. Weigh the gel slice.

2. Add 4 volumes of Monarch Gel Dissolving Buffer to the tube with the gel slice (e.g., 400 μ l buffer per 100 mg agarose). If the gel slice is >150 mg, consider reducing the amount of Gel Dissolving Buffer to 3 or 3.5 volumes to minimize the guanidine salt present in the workflow.
3. Incubate the sample between 37–55°C (typically 50°C), inverting periodically until the gel slice is completely dissolved (generally 5–10 minutes).
4. Insert the column into collection tube and load sample onto the column. Spin for 1 minute, then discard flow-through.
5. Re-insert column into collection tube. Add 200 μ l DNA Wash Buffer and spin for 1 minute. Discarding flow-through is optional.
6. Repeat wash (Step 5).
7. Re-spin the column and tubes.
8. Transfer column to a clean 1.5 ml microfuge tube. Use care to ensure that the tip of the column has not come into contact with the flow-through. If in doubt, re-spin for 1 minute before placing into clean microfuge tube.
9. Add ≥ 6 μ l of DNA Elution Buffer to the center of the matrix. Wait for 1 minute, and spin for 1 minute to elute DNA.

Notes:

- *Step2: If the volume of the dissolved sample exceeds 800 μ l, the loading of the sample onto the column should be performed in multiple rounds to not exceed the volume constraints of the spin column.*
- *Step3: For DNA fragments > 8 kb, an additional 1.5 volumes of water should be added after the slice is dissolved to mitigate the tighter binding of larger pieces of DNA (e.g., 100 mg gel slice: 400 μ l Gel Dissolving Buffer: 150 μ l water). Failure to dissolve all the agarose will decrease the recovery yield due to incomplete extraction of the DNA and potential clogging of the column by particles of agarose.*

Protocol references

Monarch DNA Gel Extraction Kit Protocol

<https://www.neb.com/en-us/protocols/2015/11/23/monarch-dna-gel-extraction-kit-protocol-t1020?srsId=AfmBOopHTsWwc4YwVyiO2Plq7zqd6dZL-3cMFxWIUQeZdJGdAVZ68Ltz>