



Apr 05, 2024

Ventral Midbrain Genomic PCR

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg324zzv25/v1>

madalynn.erb Erb¹

¹Van Andel Research Institute

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg324zzv25/v1>

Protocol Citation: madalynn.erb Erb 2024. Ventral Midbrain Genomic PCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg324zzv25/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 12, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 97844

Keywords: ASAPCRN, ventral midbrain genomic pcr this protocol detail, ventral midbrain genomic pcr

Disclaimer

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

This protocol details ventral midbrain genomic PCR.

Materials

 Stainless Steel Brain Matrices, 1.0mm **Stoelting Catalog #51386**

 QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**

 GeneJET PCR Purification Kit **Thermo Fisher Scientific Catalog #K0702**

	A	B	C
	Reagent	Final conc.	ul/sample
	5X KAPA2G Buffer A	1.3X	6.5 µl
	25 mM MgCl ₂	2.60 mM	2.6 µl
	10 mM KAPA dNTP Mix	0.26 mM	0.65 µl
	Forward primer (10 µM) CTGCAGCTTCGAGAGGAAAG	0.5 µM	0.5 µl
	Flox reverse primer (10 µM) CACTCTGTCCTCAGGCTTTC	0.5 µM	0.5 µl
	KO reverse primer (10 µM) AGGTGGGAATCGGGCTAGAG	0.5 µM	0.5 µl
	50% Glycerol	6.50 %	3.25 µl
	5 U/ µl KAPA2G Fast Hotstart DNA Polymerase	0.5 U/ul	0.1 µl
	DNA (diluted to 25 ng/uL)	150 ng total	6.0 µl
	H ₂ O	to 25 ul	4.4 µl



Brain tissue collection

- 1 Euthanize mouse via cervical dislocation.
- 2 Isolate $\pm$ 2 mm coronal midbrain section using stainless steel brain matrix (Stoelting 51386).
 - 2.1 Remove cortex and dorsal midbrain tissue.
 - 2.2 Separate ipsilateral and contralateral ventral midbrain regions.
 - 2.3 Immediately freeze on dry ice and store tissue at -80°C until DNA extraction.

DNA extraction from frozen brain tissue

- 3 Use Qiagen DNeasy Blood and Tissue Kit (Cat 69504).
- 4 Equilibrate ventral midbrain tissue to Room temperature .

Note

Starting material amount ~ 15 mg of tissue.

- 5 Cut tissue into small pieces.
 - 5.1 Use a pipette tip to transfer tissue to a sterile 6cm dish.
 - 5.2 Add $180\text{ }\mu\text{L}$ of **Buffer ATL** to the tissue.





5.3 Dice up tissue with sterile scalpel.

5.4 Transfer tissue and buffer to 1.5mL etube.



6 Add  20 μL **Proteinase K** to each sample.



6.1 Mix thoroughly by vortexing.



6.2 Incubate at  56 °C until samples are completely lysed.



Note

Use thermomixer – check samples after  01:00:00 , may take up to  03:00:00 .

7 Add  4 μL **RNase A** ( 100 μL) to each sample.



7.1 Mix by vortexing.



7.2 Incubate at  Room temperature for  00:02:00 .

2m



8 Mix **Buffer AL** with ethanol 1:1.



8.1 Add  400 μL Buffer AL / ethanol mix to each sample.



8.2 Vortex for  00:00:15 .

15s



- 9 Transfer samples to DNeasy Mini spin columns (placed in 2mL collection tubes). 
- 9.1 Centrifuge at $\geq 6000 \times g$ for 00:01:00. 

- 9.2 Discard flow through.
- 10 Transfer column to new collection tube. 
- 10.1 Add $500 \mu\text{L}$ **Buffer AW1**. 
- 10.2 Centrifuge at $\geq 6000 \times g$ for 00:01:00. 

- 10.3 Discard flow through and collection tube.
- 11 Place the column in a new collection tube – add $500 \mu\text{L}$ **Buffer AW2**. 
- 11.1 Centrifuge $20000 \times g$ for 00:03:00 to dry the column. 

- 11.2 Carefully remove the column to avoid contamination with residual ethanol.
- 11.3 If the column touches ethanol flow through, spin again in new collection tube for 00:01:00. 

- 12 Place the column in clean 1.5mL etube – add $150 \mu\text{L}$ of elution buffer (**Buffer AE**). 
- 12.1 Incubate at Room temperature for 00:01:00. 



12.2 Centrifuge at $\geq 6000 \times g$ for 00:01:00 to elute DNA.

1m



Genomic PCR

13 DNA concentration measured using a NanoDrop One Spectrophotometer (Thermofisher Scientific).

13.1 All samples diluted to 25 μL with **Buffer AE**.



14 Use the Kapa2g Fast HotStart PCR Kit (Roche 07960930001) according to the manufacturer's instructions.



	A	B	C
	Reagent	Final conc.	ul/sample
	5X KAPA2G Buffer A	1.3X	6.5 μl
	25 mM MgCl_2	2.60 mM	2.6 μl
	10 mM KAPA dNTP Mix	0.26 mM	0.65 μl
	Forward primer (10 μM)	0.5 μM	0.5 μl
	CTGCAGCTTCGAGAGGAAAG		
	Flox reverse primer (10 μM)	0.5 μM	0.5 μl
	CACTCTGTCCTCAGGCTTTC		
	KO reverse primer (10 μM)	0.5 μM	0.5 μl
	AGGTGGGAATCGGGCTAGAG		
	50% Glycerol	6.50%	3.25 μl
	5 U/ μl KAPA2G Fast Hotstart DNA Polymerase	0.5 U/ μl	0.1 μl
	DNA (diluted to 25ng/ μL)	150ng total	6.0 μl
	H2O	to 25 μl	4.4 μl

PCR Cycle:

	A	B	C	D
	Step	Temp (°C)	Time	Note
	1	94.0 °C	5 min.	
	2	94.0 °C	30 sec.	
	3	65.0 °C	15 sec.	-0.5 °C per cycle decrease
	4	68.0 °C	1 sec.	
	5			repeat steps 2-4 for 10 cycles (touchdown)
	6	94.0 °C	30 sec.	
	7	60.0 °C	15 sec.	
	8	72.0 °C	1 sec.	
	9			repeat steps 6-8 for 20 cycles
	10	72.0 °C	5 min.	
	11	4.0 °C	hold	Hold

15 Run samples on 2% Agarose gel at 100V for 35-45 minutes.



	A	B
	Wild type allele	400 bp
	Flox allele	500 bp
	KO allele	270 bp