



Feb 02, 2024

Very low-density lipoprotein receptor (VLDLR)-C-tag purification from HEK293E cells

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkoqw1v5r/v1>

Andreas Bracher¹, Patricia Yuste-Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkoqw1v5r/v1>

Protocol Citation: Andreas Bracher, Patricia Yuste-Checa, F Ulrich Hartl 2024. Very low-density lipoprotein receptor (VLDLR)-C-tag purification from HEK293E cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkoqw1v5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 25, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94601

Keywords: ASAPCRN, tag purification from hek293e cell, tag protein from hek293e cell, hek293e cell, density lipoprotein receptor, lipoprotein receptor, tag protein, vldlr, tag purification, purification, protein

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to purify recombinant Very low-density lipoprotein receptor (VLDLR)-C-tag protein from HEK293E cells.

Attachments



VLDLR-C-tag.purifica...

317KB

Materials

Buffers

■ Binding buffer:

	A	B
	Tris-HCl pH 7.2	20 mM
	NaCl	100 mM
	CaCl ₂	0.5 mM

■ Elution buffer:

	A	B
	Tris-HCl pH 7.0	20 mM
	MgCl ₂	2 M
	CaCl ₂	2 mM

 FreeStyle[®]; 293 Expression Medium **Thermo Fisher Catalog #12338018**

 CaptureSelect[®]; tPA Affinity Matrix **Thermo Fisher Catalog #2943430005**

 Amersham NAP-25 Columns **Cytiva Catalog #17-0852-01**




VLDLR-C-tag expression

- 1 Express VLDLR-C-tag in HEK293E cells cultured in FreeStyle 293 Expression Medium for  96:00:00 
- 2 Centrifuge culture and keep conditioned medium. 
- 3 Dialyze  300 mL conditioned medium  Overnight against  10 L Binding buffer. 

CaptureSelect C-tag affinity chromatography

- 4 Load dialyzed conditioned medium onto a CaptureSelect C-tag affinity column previously equilibrated with binding buffer (column volume, CV: 4 mL slurry for 300 mL dialyzed media) by gravity flow at  4 °C .
- 5 Wash the column with 5 CV of Binding buffer. 
- 6 Elute VLDLR-C-tag protein with 6× 1 mL of Elution buffer. Collect fractions of  1 mL .
- 7 Analyze eluted fraction by SDS-PAGE and Coomassie blue staining.



- 8 Pool fractions containing VLDLR-C-tag protein.
- 9 Exchange protein buffer to Binding buffer with a NAP-25 desalting column previously equilibrated with Binding buffer.
- 10 Collect the fractions containing protein, concentrate by ultrafiltration to $>1 \text{ mg mL}^{-1}$, aliquot and flash-freeze purified VLDLR-C-tag in liquid nitrogen for storage at -70°C .

Note

Approximate yield: From 300 ml of conditioned media around 0.6 mg of pure VLDLR-C-tag were obtained. Yield can be significantly increased if the VLDLR chaperone, RAP is co-overexpressed during protein production.

