

Dec 17, 2025

Viable PBMCs isolation using density-gradient separation

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl48ejrvo5/v1>

Layla Drwesh¹, Loris Bandirali², Julia Fitzgerald¹

¹Hertie Institute for Clinical Brain Research, University of Tübingen;

²Department of Brain and Behavioural Sciences, University of Pavia, Pavia, Italy.



Julia Fitzgerald

Hertie Institute for Clinical Brain Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl48ejrvo5/v1>

Protocol Citation: Layla Drwesh, Loris Bandirali, Julia Fitzgerald 2025. Viable PBMCs isolation using density-gradient separation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl48ejrvo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: December 16, 2025

Last Modified: December 17, 2025

Protocol Integer ID: 235141

Keywords: high quality pbmcs from whole blood, gradient separation collection of pbmc, viable pbmcs isolation, alignment across ongoing parkinson, biomarker development, gradient separation collection, high quality pbmc, pbmc, ongoing parkinson, incorporating sepmate tube, stem cell technology, different lab, density gradient method, sepmate tube, blood collection

Funders Acknowledgements:

Michael J Fox Foundation

Grant ID: 026376

Abstract

Collection of PBMCs for biomarker development is an essential part of translational research. Several methods are employed by different labs, including our own. Here we describe in detail, a method (in place in our lab from 2025) to isolate high quality PBMCs from whole blood (EDTA blood) within ~45 minutes from blood collection. We use a density gradient method, incorporating SepMate tubes (Stem Cell Technologies). We would like to share this protocol for the purpose of standardization, comparison and alignment across ongoing Parkinson's disease cohort studies.

Guidelines

Working with blood. Before starting, make sure that all local safety and occupational health rules are followed. For example, vaccinations, safety instructions.

Materials

Guidelines

Except for the centrifugation steps, all sample processing was carried out under a biosafety cabinet in sterile conditions.

Store Lymphoprep, CryoStor CS10, PBS+2%FBS and Erylysis buffer at 4°C.

Before start

Whole peripheral blood samples are collected in EDTA tubes by professional healthcare workers in medical examination rooms.

Prepare 100mL of PBS+2%FBS for two blood vials (~15mL) under the hood.

Prepare 50mL of Erylysis Buffer as follows:

- 0.1mM EDTA
 - 155mM NH_4Cl
 - 10mM KHCO_3
 - Adjust the pH to 7.3-7.4 with HCl or NaOH
- Bring Lymphoprep, PBS+2%FBS solution, Erylysis buffer, CryoStor CS10 and Mr. Frosty at room temperature (18°C-25°C).

Materials

SepMate™ tubes (50mL)

for density gradient separation **STEMCELL Technologies Catalog #85450.**

Falcon tubes (50mL) **Corning**

Catalog #352070

Falcon Pipettes (25mL) **Corning**

Catalog #357525

Falcon Pipettes (10mL) **Corning**

Catalog #356551

Falcon Pipettes (5mL) **Corning**

Catalog #357543

Cryo. S (2mL) **Greiner Bio-One Catalog**

#122263

Lymphoprep™ **STEMCELL**

Catalog #07851/07861

CryoStor CS10 **STEMCELL Catalog #07930**

1X Dulbecco's Phosphate Buffer

Saline **Sigma-Aldrich Catalog #D8537**

Fetal Bovine Serum, Premium **Thermofisher**

Catalog #A5670701

Protocol materials

⊗ Lymphoprep™ 500 mL **STEMCELL Technologies Inc. Catalog #7851**

⊗ Dulbecco's Phosphate Buffer Saline **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8537**

⊗ Fetal Bovine Serum, Premium **Gibco - Thermo Fisher Scientific Catalog #A5670701**

⊗ CryoStor CS10-100ml **STEMCELL Technologies Inc. Catalog #07930**

⊗ 50ml Falcon tubes **Corning Catalog #352070**

⊗ SepMate-50 **STEMCELL Technologies Inc. Catalog #85450**

⊗ DMSO **Panreac AppliChem Catalog #A3672**

⊗ CRYO.S 2 ML PP ROUND BOTTOM INTERNALTHREAD **greiner bio-one Catalog #122263-2DG**

Safety warnings

⚠ Please wear protective clothing, PPE and follow local safety guidelines.

Ethics statement

All research must adhere to the current, updated Helsinki declaration.

Prior ethical approval must be granted.

Informed consent must be given.

Before start

Ethics approvals, consent and local safety rules in place.

See section preparation below.

Preparation

- 1 Prepare  100 mL of
 Dulbecco's Phosphate Buffer Saline **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8537**
+2%
 Fetal Bovine Serum, Premium **Gibco - Thermo Fisher Scientific Catalog #A5670701**
for two blood vials (~  15 mL total) under the hood
- 2 Prepare  50 mL stock of Erylysis Buffer as follows: EDTA  0.1 millimolar (mM) ,
NH₄Cl  155 millimolar (mM) , KHCO₃  10 millimolar (mM) . Adjust the pH to 7.3 - 7.4 with HCl or NaOH. Erylysis buffer can be stored long term at  4 °C
- 3 Bring  Lymphoprep™ 500 mL **STEMCELL Technologies Inc. Catalog #7851** ,
 Dulbecco's Phosphate Buffer Saline **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8537**
+2%
 Fetal Bovine Serum, Premium **Gibco - Thermo Fisher Scientific Catalog #A5670701**
solution, Erylysis buffer,
 CryoStor CS10-100ml **STEMCELL Technologies Inc. Catalog #07930** and
Mr. Frosty at room temperature ( 18 °C -  25 °C)
- 4 Add  15 mL
 Lymphoprep™ 500 mL **STEMCELL Technologies Inc. Catalog #7851**
to the  SepMate-50 **STEMCELL Technologies Inc. Catalog #85450** using a
 25 mL pipette by carefully pipetting it through the central hole of the
 SepMate-50 **STEMCELL Technologies Inc. Catalog #85450** insert (plastic insert inside the tube).

**Note**

If needed, use a sterile (non-filter tip)  1000 μ L pipette tip placed into the insert hole as a way to stabilize and guide the strip pipette and Lymphoprep solution through the hole.

PBMC isolation

- 5 Add  15 mL -  17 mL blood into a  50ml Falcon tubes **Corning Catalog #352070** . Add equal volume of  Dulbecco's Phosphate Buffer Saline **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8537** +2%  Fetal Bovine Serum, Premium **Gibco - Thermo Fisher Scientific Catalog #A5670701** to the  50ml Falcon tubes **Corning Catalog #352070** . Close the lid and mix gently by inverting 2 times.
- 6 Keeping the  SepMate-50 **STEMCELL Technologies Inc. Catalog #85450** vertical, add the diluted blood sample by slowly pipetting it down the side of the tube. The sample will mix with the density gradient medium above the insert. Close the lid.
- 7 Centrifuge the  SepMate-50 **STEMCELL Technologies Inc. Catalog #85450** at 1200xg for  00:10:00 at  Room temperature with the brakes on (medium deceleration).
- 8 With care, remove the  SepMate-50 **STEMCELL Technologies Inc. Catalog #85450** lid. To reduce platelet contamination, pipette off and discard about  10 mL of supernatant from the surface using a  10 mL pipette – well above the enriched mononuclear cell layer.



- 9 Pour the remaining top layer containing the enriched mononuclear cells into a new  50ml Falcon tubes **Corning Catalog #352070** by inverting the  SepMate-50 **STEMCELL Technologies Inc. Catalog #85450** to 180 degrees. To avoid pouring the red blood cells (RBCs), do not hold the  SepMate-50 **STEMCELL Technologies Inc. Catalog #85450** in the inverted position for longer than  00:00:02
- 10 Top up the enriched mononuclear cells in tube with  Dulbecco's Phosphate Buffer Saline **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8537**
+2%  Fetal Bovine Serum, Premium **Gibco - Thermo Fisher Scientific Catalog #A5670701**
to  50 mL and centrifuge at 300xg for  00:05:00 to pellet PBMCs
- 11 Carefully discard the supernatant using a vacuum aspirator without disturbing the PBMC pellet
- 12 To get rid of residual red blood cells (RBCs), resuspend the pellet in  5 mL -  7 mL of Erylysis Buffer and incubate at  Room temperature for  00:05:00
- 13 Top up the remaining enriched mononuclear cells with  Dulbecco's Phosphate Buffer Saline **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8537**
+2%  Fetal Bovine Serum, Premium **Gibco - Thermo Fisher Scientific Catalog #A5670701**
in  50ml Falcon tubes **Corning Catalog #352070** and centrifuge at 300xg for  00:05:00 to pellet PBMCs
- 14 Carefully discard the supernatant using a vacuum aspirator without disturbing the PBMC pellet



- 15 Resuspend the PBMC pellet in  3 mL -  4 mL of  CryoStor CS10-100ml **STEMCELL Technologies Inc. Catalog #07930** or  Fetal Bovine Serum, Premium **Gibco - Thermo Fisher Scientific Catalog #A5670701**
+10%  DMSO **Panreac AppliChem Catalog #A3672**
- 16 Count the cell suspension manually or using an automated cell counter . Adjust the volume of freezing media to reach the desired concentration of PBMCs (1 - 20*10⁶ cells/mL)
- 17 Prepare  0.5 mL -  1 mL aliquots in  CRYO.S 2 ML PP ROUND BOTTOM INTERNALTHREAD **greiner bio-one Catalog #122263-2DG**
and place them in Mr. Frosty. Cryopreserved PBMCs should be transferred into liquid nitrogen for long term storage after at least 24 hours at  -80 °C . If cells are being used for experimentation within a few weeks of harvesting, they can be kept at  -80 °C .

Documentation

- 18 Document the following information:
ID number and date
Time the blood was received. Time at first centrifugation step. Time the PBMCs were frozen. Time/date of transfer to liquid nitrogen.
Cell number, concentration and number of cryovials.

Acknowledgements

Claudia Schulte, Christian Deuschle and Ann-Kathrin Hauser at The Hertie Institute for Clinical Brain Research and The Neurobiobank of The Hertie Institute for Clinical Brain Research for their constant support and teamwork across projects.

The Michael J Fox Foundation

The Michael J Fox Foundation PINK1/PRKN consortium