

Nov 16, 2020

# Video protocol for sorting *Drosophila* pupae

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.bpq3mmyn](https://dx.doi.org/10.17504/protocols.io.bpq3mmyn)

Carolyn Elya<sup>1</sup>

<sup>1</sup>Harvard University



Carolyn Elya

Harvard University



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bpq3mmyn](https://dx.doi.org/10.17504/protocols.io.bpq3mmyn)

**Protocol Citation:** Carolyn Elya 2020. Video protocol for sorting *Drosophila* pupae. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bpq3mmyn>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** November 16, 2020

**Last Modified:** November 16, 2020


**Protocol Integer ID:** 44539

## Abstract

This is my favorite way to sort pupae.

Credit to Tom Alisch & Dave Zucker (@FlySorter) for discovering a new use for transparencies!!

1 Video summary of the method:



- 2 Cut out a small piece of a projector sheet (e.g. Hygloss Overhead Projector Sheets, [https://www.amazon.com/Hygloss-75910-Products-Transparency-Projectors/dp/B077762QK5/ref=sr\\_1\\_8?dchild=1&keywords=transparency&qid=1605473049&sr=8-8](https://www.amazon.com/Hygloss-75910-Products-Transparency-Projectors/dp/B077762QK5/ref=sr_1_8?dchild=1&keywords=transparency&qid=1605473049&sr=8-8)).

You'll want to make sure the piece is long enough to completely cover the inside of the vial and not so wide as to reach beyond the surface of the food. Measure diameter of the vial mouth to determine a suitable length ( $d \cdot \pi + \sim 1 \text{ cm}$ ) and space between food and cap to determine width.

- 3 Roll up the projector sheet and place into a vial any time before flies reach 3rd instar.
- 4 Pull out the transparency after pupation.
- 5 Sort pupae using a paintbrush, using water to dissolve the glue holding pupae to the projector sheet. Placing a Kimwipe under the projector sheet will help you maintain a clean workspace and prevent inter-vial contamination.
- 6 Pupae can be directly transferred to a new vial with food or placed on a new piece of projector sheet which can then be rolled up and placed into a vial.



- 7 Between genotypes, check paintbrush bristles for any adhered larvae or pupae and remove before proceeding. Rinse paintbrush thoroughly with water and 70% ethanol to decontaminate.