Oct 30, 2023

Village Nuclei Isolation With Optiprep

DOI

dx.doi.org/10.17504/protocols.io.36wgq3bmxlk5/v1

liv_spina¹, Tara McDonald¹, Nora Reed¹, Alyssa Lutservitz¹

¹Harvard Medical School

McCarroll Lab







DOI: dx.doi.org/10.17504/protocols.io.36wgq3bmxlk5/v1

Protocol Citation: liv_spina, Tara McDonald, Nora Reed, Alyssa Lutservitz 2023. Village Nuclei Isolation With Optiprep. protocols.io https://dx.doi.org/10.17504/protocols.io.36wgq3bmxlk5/v1

Manuscript citation: Wells et al., PMID <u>36796362</u>

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: October 30, 2023

Last Modified: October 30, 2023

Protocol Integer ID: 90142

Keywords: fresh-frozen nuclei, nuclei extraction, snRNAseq, RNAseq, human tissue, brain tissue, multiple donors

Abstract

Isolation of nuclei from fresh-frozen brain tissue from sets of multiple (typically 2-20) human donors for analysis as a "cell village" (**Wells et al., PMID** <u>36796362</u>) in which nuclei from all donors are analyzed together. Adapted from dx.doi.org/10.17504/protocols.io.bs99nh96, Luciano G Martelotto, with modifications to optimize for human brain tissue and allow the "cell village" approach.

Materials

Supplies

- Scalpels
- Glass slides
- Dounce
- 20 µm vacuum filter
- Eppendorf tubes (1.5 mL and 5 mL)
- Eppendorf or Ranin pipette tips
- Dry ice
- Metal plate
- OCT (Optimal cutting temperature compound)
- RNAse free water
- Cell counting supplies (LUNA-FL)

Other Reagents:

- PBS
- BSA
- RNAse inhibitor (i)
- Nuclei EZ lysis buffer (NUC201-1KT)
- OptiPrep Density Gradient Medium (60% Iodixanol G60) (ab286850)

Stock solutions:

- 500 mM tricine: 8.96 g of tricine in 100 mL water
- 1 M KCI: 7.45 g of KCl in 100 mL water
- 1 M MgCl2x6H2O: 20.3 g of MgCl2x6H2O in 100 mL water

Before start

- Turn on and chill centrifuge to 4°C
- Prepare all the reagents needed on ice
- Clean glass slides with ethanol
- Gather the dounce, pestles, and scalpels on ice

Solutions to make fresh before starting experiment:

- Nuclei EZ lysis buffer -- Pour 12 mL of buffer into the dounce and reserve 1 mL for wetting down the sides of the dounce
- PBSA + i (1% BSA and 1 U/µL RNase inhibitor in 1X PBS) -- For powder BSA, add 0.01 g of BSA per 1 mL of 1X PBS
- G30: (Mix the G60 and GD together thoroughly before adding the GH)

Before Starting

1 Gather Supplies

- Scalpels
- Glass slides
- 14 mL Dounce
- 20 µm vacuum filter
- Eppendorf tubes (1.5 mL and 5 mL)
- Eppendorf or Rainin pipette tips
- Dry ice
- Metal plate
- OCT (Optimal cutting temperature compound)
- RNAse free water
- Cell counting supplies (LUNA-FL)

Other Reagents:

- PBS
- BSA
- RNAse inhibitor (i)
- Nuclei EZ lysis buffer (NUC201-1KT)
- OptiPrep Density Gradient Medium (60% Iodixanol G60) (ab286850)

Stock solutions:

- 500 mM tricine: 8.96 g of tricine in 100 mL water
- 1 M KCI: 7.45 g of KCl in 100 mL water
- 1 M MgCl2x6H2O: 20.3 g of MgCl2x6H2O in 100 mL water

1.1 **Optiprep Buffers (can be made in advance and kept on hand)**

| A | В | С | D |
|-----------------|--|----------------|--|
| GD | | GH | |
| Component | Volume | Component | Volume |
| 500 mM Tricine | 24 mL | Sucrose | 8.5 g |
| 1M KCI | 15 mL | Water | 50 mL |
| 1M MgCl2 x 6H2O | 3 mL | 500 mM Tricine | 4 mL |
| Water | 50 mL | KCI | 2.5 mL |
| 1M KOH | Adjust to pH 7.8 | 1M MgCl2x6H2O | 0.5 mL |
| | Bring volume up to 100 mL with water (keep at room temp) | 1М КОН | Adjust to pH 7.8 |
| | | | Bring volume up to 100 mL with water (keep at 4°C) |

- 1.2 Turn on and chill centrifuge to 4°C
 - Prepare all the reagents needed on ice
 - Clean glass slides with ethanol
 - Gather the dounce, pestles, and scalpels on ice
- 1.3 Solutions to make fresh before starting experiment

Nuclei EZ lysis buffer

 Pour 12 mL of buffer into the dounce and reserve 1 mL for wetting down the sides of the dounce

PBSA + i (1% BSA and 1 U/ μ L RNase inhibitor in 1X PBS)

• For powder BSA, add 0.01 g of BSA per 1 mL of 1X PBS

G30: (Mix the G60 and GD together thoroughly before adding the GH)

| _ | A | В | С | D | E | F |
|---|-------------------|--------|------------------------------|--------------------|---------------------------|-----------|
| | G30 | | PBSA + i | | Nuclei EZ Lysis Buffer | |
| _ | Component | Volume | Component | Volume | Component | Volume |
| _ | G60 (OptiPrep) | 6.0 mL | PBS (pH 7.4) (1X) | 1 mL of 10X PBS | EZ Lysis Buffer | 12,975 µL |
| _ | GD | 1.2 mL | BSA (100X) | 1 mL of 10% BSA | RNAse Inhibitor | 325 µL |
| | GH | 4.8 mL | RNase Inhibitor (200X) | 50 µL | | |
| | | | RNase free water | 7,950 μL | | |

Tissue Homogenization

1h 30m

- 2 On a glass slide **on dry ice**, shave off pieces of tissue with a scalpel until you have 50±3 mg of tissue
- 2.1 For managing tissue, use OCT to adhere tissue to glass slide

protocols.io Part of **Springer NATURE**

- 3 Place tissue pieces on a boat on dry ice. Hold the tissue on dry ice until all tissue has been cut
- 4 Transfer all tissue to dounce filled with Nuclei EZ Lysis Buffer
- 4.1 Use 1mL of the reserved Nuclei EZ Lysis Buffer to wet down the side of the dounce if tissue is stuck on the side of the dounce
- 5 Gently dounce tissue on ice 20 times with pestle "A" then 20 times with pestle "B" until the tissue is homogenized
- 6 Incubate on ice for 10 minutes

Nuclei Isolation

- 7 Wet a 20 μm vacuum filter with 1 mL of PBSA + i
- 8 Vacuum filter the homogenized tissue
- 8.1 If the filter clogs, remove any unfiltered volume with a pipette, replace the filter, and continue filtering
- 9 Transfer the homogenized tissue to (3) 5 mL tubes
- 10 Centrifuge the samples at 4°C for 5 min (500 x g)
- 11 Remove and discard the supernatant (be careful not to disturb cell pellets) and resuspend each pellet in 900 µL G30
- 11.1 Alternatively, split pellet into (9) 1.5 mL tubes (total), resuspend each pellet in 300 μL of G30, then underlay each sample with 1 mL of G30

1h 30m

1h

protocols.io | <u>https://dx.doi.org/10.17504/protocols.io.36wgq3bmxlk5/v1</u> October 30, 2023

12 Gently underlay each 5mL tube with 3 mL G30

protocols.io Part of Springer NATURE

Place the pipette tip filled with G30 at the bottom of the 5mL tube. Slowly dispense G30 underneath the cell suspension. A clear separation will form between homogenate and G30 (cloudy on top, clear on the bottom). Do not disturb the separation between the layers

- 13 Centrifuge the samples at 4°C for 10 min (4,255 x g)
- 13.1 For capturing smaller cell types, spin at 8,000 x g for 20 minutes
- 14 Remove and discard the supernatant
- 14.1 Remove the top (very viscous) layer from the tube first- this layer contains fats and must be removed before removing the rest of the supernatant

Nuclei Wash, Resuspension, and Quantification

- 15 Resuspend the pellets in 1 mL PBSA + i
- 16 Centrifuge the samples at 4°C for 5 min (500 x g) (single wash)
- 17 Remove the supernatant and resuspend the pellets in 200 1,000 μ L of PBSA + i (depending on the size of the pellet)
- 18 Count nuclei using LUNA-FL counter

Preserving Leftover Nuclei

- 19 Optional if you have excess nuclei, or if you need to pause before continuing.
- 20 Centrifuge the nuclei at 4°C for 5 min (500 x g)

15m

20m

protocols.io Part of SPRINGER NATURE

- 21 Remove the supernatant and resuspend the nuclei in CryoStor (CryoStor® CS10), mix by pipetting
- 21.1 ≤ 1k nuclei / µL of CryoStor
- 22 Transfer the nuclei to cryotubes and store them in a Mr Frosty in a -80°C freezer until frozen (1 day)
- 23 For long term storage, transfer cryotubes to regular storage in a -80°C freezer

Thawing Cryo-Stored Nuclei

24 Thaw the CryoStored nuclei

- 25 Centrifuge the nuclei at 4°C for 10 min (1200 x g)
- 26 Remove the supernatant and resuspend in 1mL PBSA + i
- 27 Centrifuge the nuclei at 4°C for 5 min (500 x g) (single wash)
- 28 Remove the supernatant and resuspend the pelleted nuclei in PBSA + i
- 28.1 Count the nuclei before proceeding to next steps

Protocol references

Wells et al., PMID <u>36796362</u>

Adapted from dx.doi.org/10.17504/protocols.io.bs99nh96, Luciano G Martelotto

