

May 16, 2020

Viral RNA extraction low-cost protocol optimized for SARS-Cov2 at AGROSAVIA

DOI

<https://dx.doi.org/10.17504/protocols.io.bggvjtw6>

Alejandro Caro-Quintero¹, Camilo Ruiz-Avila¹, Andrea Navarrete¹, Johan Bernal¹, Giuliana Daza¹, Roxana Yockteng¹

¹AGROSAVIA/Universidad Nacional de Colombia

Coronavirus Method De...

Reclone.org (The Reage...

1 more workspace



Alejandro Caro-Quintero

AGROSAVIA/Universidad Nacional de Colombia

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bggvjtw6>

External link: <http://10.13140/RG.2.2.34818.50881>



Protocol Citation: Alejandro Caro-Quintero, Camilo Ruiz-Avila, Andrea Navarrete, Johan Bernal, Giuliana Daza, Roxana Yockteng 2020. Viral RNA extraction low-cost protocol optimized for SARS-Cov2 at AGROSAVIA. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bggvjtw6>

Manuscript citation:

DOI: 10.13140/RG.2.2.34818.50881/1

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 16, 2020

Last Modified: May 16, 2020

Protocol Integer ID: 37109

Keywords: commercial viral rna extraction kit, viral rna extraction protocol, second version of an rna extraction protocol, rna extraction protocol, cost of rna extraction, rna extraction, magnetic beads from viral pathogen, viral rna, throughput nucleic acid extraction, nucleic acid extraction, viral pathogen, respiratory syndrome coronavirus, rna, novel coronavirus, cov2 at agrosavia, cov2, ct results like the commercial kit, throat swab, viral transfer media, nucleic acid, extraction, ct result, using magnetic bead, agrosavia laboratory, swab

Abstract

Here we present the second version of an RNA extraction protocol for SARs-cov2 using magnetic beads. This protocol was optimized in AGROSAVIA laboratories and is presented as an alternative to the limitations of commercial viral RNA extraction kits. The viral RNA extraction protocol was adapted from the protocol entitled "SARS-CoV-2 RNA purification from nasal / throat swabs collected in Viral Transfer Media" at <https://bomb.bio/> and which is based on methods previously reported by the same group (Oberacker et al 2019). The optimized version here achieves Ct results like the commercial kits and reduces the cost of RNA extraction per sample by 90%, that is, 2000 Colombian pesos or 50 cents per sample. We share this protocol in the hope that it can ensure the sustainability of the diagnosis in Colombia and other countries where it is difficult and expensive to import commercial kits.

REFERENCES

Oberacker, P., Stepper, P., Bond, D. M., Höhn, S., Focken, J., Meyer, V., ... & Hore, S. R. (2019). Bio-On-Magnetic-Beads (BOMB): Open platform for high-throughput nucleic acid extraction and manipulation. *PLoS biology*, 17(1), e3000107.

Yu, F., Qiu, T., Zeng, Y., Wang, Y., Zheng, S., Chen, X., & Chen, Y. (2018). Comparative evaluation of three Preprocessing Methods for extraction and detection of Influenza A Virus Nucleic Acids from sputum. *Frontiers in Medicine*, 5, 56.

He, H., Li, R., Chen, Y., Pan, P., Tong, W., Dong, X., ... & Yu, D. (2017). Integrated DNA and RNA extraction using magnetic beads from viral pathogens causing acute respiratory infections. *Scientific reports*, 7, 45199.

Corman, V.M., Landt, O., Kaiser, M., Molenkamp, R., Meijer, A., Chu, D.K., Bleicker, T., Brünink, S., Schneider, J., Schmidt, M.L. and Mulders, D.G., 2020. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *Eurosurveillance*, 25(3).

Kumar, M., Mazur, S., Ork, B.L., Postnikova, E., Hensley, L.E., Jahrling, P.B., Johnson, R. and Holbrook, M.R., 2015. Inactivation and safety testing of Middle East respiratory syndrome coronavirus. *Journal of virological methods*, 223, pp.13-18

Attachments



[English translation ...](#)

1,013KB



Guidelines

Here we present the second version of an RNA extraction protocol for SARs-cov2 using magnetic beads. This protocol was optimized in AGROSAVIA laboratories and is presented as an alternative to the limitations of commercial viral RNA extraction kits. The viral RNA extraction protocol was adapted from the protocol entitled "SARS-CoV-2 RNA purification from nasal / throat swabs collected in Viral Transfer Media" at <https://bomb.bio/> and which is based on methods previously reported by the same group (Oberacker et al 2019). The optimized version here achieves Ct results like the commercial kits and reduces the cost of RNA extraction per sample by 90%, that is, 2000 Colombian pesos or 50 cents per sample. We share this protocol in the hope that it can ensure the sustainability of the diagnosis in Colombia and other countries where it is difficult and expensive to import commercial kits.



Materials

Reagent List

	Reagent
	Guanidine Thiocyanate (GITC)
	Sarkosyl
	Magnetic beads GE Healthcare Sera-Mag Carboxylate-Modified
	Isopropanol
	EDTA
	Ethanol
	Proteinase K (20 mg/mL)

Preparation of reagents

3.1 Paramagnetic beads

- Take 1 ml of Sera-Mag pearls, be sure to mix well before you start.
- Wash the beads 3 times in 1 ml of TE buffer (this is to remove sodium azide from the solution), use the magnet to separate the beads from the washes, remove the TE.
- Re-suspend in 50 ml of TE buffer.

3.2 GITC Lysis Buffer

- 35.46 g of GITC.
- 2.5 mL of 1 M Tris HCl pH 7.6-8.0 stock solution
- 1 g of Sarkosyl.
- 2 ml of 0.5 M EDTA stock solution.
- Fill to 50 ml with sterile distilled water.

Reagent prices

The price of the reagents was quoted with national suppliers and the average prices for the commercial houses are presented.

Table 1. Reagents and prices

	Reagent	Quantity	Approximate Price * (pesos)/USD	Number of samples	Price per sample * (pesos)/USD
	Guanidine Thiocyanate (GITC)	500 g	\$2 000 000/500	7000	\$286/0.07
	Sarkosyl	25 g	\$ 1 470 000/367.50	12500	\$117/0.02
	Magnetic beads GE Healthcare Sera-Mag Carboxylate-Modified	15 mL	\$3 000 000/750	18750	\$160/0.04
	Isopropanol	500 mL	\$300 000/75	1250	\$240/0.06
	EDTA	250 g	\$ 250 000/62.5	1000	\$250/0.06
	Ethanol	2500 mL	\$200 000/50	6250	\$ 32/0.008
	Proteinase K (20 mg/mL)	6 mL	300 000/75	500	\$600/0.15
				Total	\$1685/0.45 (pesos/USD)

* price in Colombian pesos



Before start

The adjusted protocol presented here is susceptible to other improvements such as those reported for the detection of respiratory viruses (He et. al 2017). 2. Similar protocols and how to build the magnetic racks can be found on the page of the original Authors of the method (<https://bomb.bio/>). 3. We are not responsible for any accident, infection, or legal implication resulting from the use of protocols adjusted here or hosted at BOMB.bio.



1 Pre-processing of mucous samples.

- 1.1 Dilute  100 μL of the sample in  300 μL of Saline Solution [ 0.8 Mass / % volume], vortex mix  00:00:10 . Use only  100 μL of this dilution to start protocol. In step 2.3 add  540 μL isopropanol after incubation.

2 Sample Processing (in 1.5 ml tubes)

- 2.1 To  100 μL of GITC buffer add  12.5 μL of Proteinase K [ 20 mg/mL],  200 μL of the sample and mix by vortex for  00:00:10 .
- 2.2 Incubate the sample for  00:10:00 at  60 °C , vortexing every  00:02:00 .
- 2.3 Place the tube in a rack and add  270 μL of Isopropanol and  40 μL of magnetic beads (see preparation below). Mix by inverting 8 times.
- 2.4 Place the tube (s) in the magnetic rack for  00:10:00 without moving.
- 2.5 Without moving the tube from the magnetic rack, withdraw approximately  610 μL of the solution
- 2.6 Without moving the tube from the magnetic rack, add  150 μL of Isopropanol and wait 30 seconds and remove the isopropanol from the tube.
- 2.7 Without moving the tube from the magnetic rack, add  200 μL of [ 70 % (v/v) Ethanol, and remove from the tube.
- 2.8 Without moving the tube from the magnetic rack, repeat previous washing 2.7 and extract the Ethanol to the maximum.



- 2.9 Without moving the tube from the magnetic rack, allow drying for approximately  00:05:00 .
- 2.10 Add  30 μL of molecular grade water to the tube trying to discharge the liquid on top of the tube walls where the beads are located.
- 2.11 Remove the tube from the magnetic rack and vortex  00:00:10 to allow the beads to enter the solution and incubate  00:01:00 at room temperature in a normal rack
- 2.12 Place the tube back in the magnetic rack and for  00:03:00 to  00:05:00 , and transfer  20 μL of the bead-free elution to a new tube and store at  -80 °C .