

Sep 21, 2023

Virus and Prokaryote Enrichment in Coral DNA Metagenomes

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7p1y3gwz/v1>



Natascha Varona¹, Bailey Wallace¹, Cynthia Silveira¹

¹University of Miami Department of Biology



amz

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7p1y3gwz/v1>

Protocol Citation: Natascha Varona, Bailey Wallace, Cynthia Silveira 2023. Virus and Prokaryote Enrichment in Coral DNA Metagenomes. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7p1y3gwz/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 29, 2023

Last Modified: September 21, 2023

Protocol Integer ID: 87107

Keywords: coral, bacteria, virus, bacteriophage, metagenome, prokaryote enrichment in coral dna metagenomes tissue preparation, coral dna metagenomes tissue preparation, prokaryotes in dna metagenome, dna metagenome, dna extraction protocol for enrichment, coral, prokaryote enrichment, dna extraction, dna extraction protocol, prokaryote, dna

Funders Acknowledgements:

National Science Foundation GRFP to BW

Grant ID: 2023353157

Maytag Fellowship to NV

Grant ID: PG015171

University of Miami Provost Research Award to CS

Grant ID: UM PRA 2022-2547

Abstract

Tissue preparation and DNA extraction protocol for enrichment of viruses and prokaryotes in DNA metagenomes of corals.

Image Attribution

Cynthia Silveira



Materials

10X DNase Buffer:

- [M] 100 millimolar (mM) Tris-HCl (pH 7.5)
- [M] 25 millimolar (mM) MgCl_2
- [M] 5 millimolar (mM) CaCl_2
- Fresh or thawed $\sim 1\text{cm}^3$ coral sample
- Sterile mortar and pestle
- Sterile forceps
-  1.5 mL centrifuge epi tubes
-  0.02 μm filtered and autoclaved artificial seawater (ASW)
-  425-600 μm acid-washed glass beads (Sigma CAT#G9268)
-  15 mL conical centrifuge tubes
- Vortex
- Centrifuge
- Pipettes
- DNase [M] 10000 U/mL
- 100X EDTA [M] 0.5 Molarity (M)
- Sterile  10 mL Luer-LokTM syringes
-  8 μm pore size, 25mm flat filters (Whatman, Nucleopore Track-Etched Polycarbonate (CAT#110614))
- Sterile 25mm filter holders with Luer-LokTM connection (Swinnex CAT#SX0002500)
- AmiconTM Ultra-15 Centrifugal Filter Units (MilleporeSigmaTM CAT#UFC910008)
- NucleoSpin[®] Tissue Kit (MACHEREY-NAGEL Inc. REF:740952.10/.50/.250)

Before start

- Glass beads must be washed with 0.1 N HCl
- Incubator set to  56 °C &  70 °C
- Check if Buffer B5 and Proteinase K were prepared according to section 3 of NucleoSpin[®] Tissue Kit Manual.

Sample Preparation

3h

- 1 Collect a coral sample of approximately 1cm³ containing coral tissue, mucus, and skeleton. Use a fresh sample if possible. If not, flash-freeze the sample and keep it at  -80 °C until processing. Thaw the sample on ice before proceeding.
- 2 Weigh and record coral fragment weight.
- 3 Crush the coral sample using a sterile mortar and pestle until the skeleton and tissue fragments are the size of coarse sand. Do not overdo this step. Transfer the fragments, along with the mucus (coral homogenate), to a 2 mL microcentrifuge tube.
- 4 Suspend sample in  100 µL of Artificial Seawater (ASW) and  0.2 g of sterile glass beads.
- 5 Vortex at speed 3 for  00:05:00 min to dislodge tissue and mucus. 5m
* To avoid lysis of bacterial cells, do not vortex the sample at higher speeds or longer durations than described.
- 6 Pulse-centrifuge the sample (3-5 seconds at  3.200 x g).
- 7 Collect supernatant and transfer to a clean 15 mL microcentrifuge tube.
- 8 Bring sample up to  1 mL with sterile ASW.
- 9 Add  100 µL of 10X DNase Buffer and  2 µL of DNase 10,000 U/mL for a final concentration of 20 U/mL. Mix by gently inverting the tube.
- 10 Incubate at room temperature for  02:00:00 hours . 2h
- 11 Stop DNase activity by adding  10 µL of 100X EDTA.



- 12 Bring sample volume up to  5 mL with sterile ASW.
- 13 Transfer the sample to a 10 mL syringe and push through an 8.0 µm, 25 mm flat filter, collecting the flowthrough in a clean 15 mL centrifuge tube.
* If the sample clogs the filter, replacing the 8.0 µm filter may be necessary. In this case, pull back on the syringe plunger to recollect the sample in the syringe before transferring it to a new filter.
* In the case of samples with excessive mucus, you may further dilute the sample with sterile ASW prior to filtering.
- 14 Transfer flowthrough to an Amicon™ Ultra-15 Centrifugal Filter Unit.
- 15 Centrifuge at room temperature for  00:30:00 min at $\geq$  2600 x g . 30m

DNA Extraction

1h 12m

- 16 1. Add  200 µL Buffer T1 and  20 µL of Proteinase K from the NucleoSpin® Tissue Kit to the Amicon™ filter unit.
- 17 Rinse one side of the filter with the lysis buffer 3-5 times using a pipette.
- 18 Securely close the Amicon™ filter unit and incubate for  01:00:00 hour at  56 °C , being sure to place the tube in a position where the rinsed side of the filter remains covered in the Buffer T1 & Proteinase K mixture.
* Some Amicon™ filter unit lids may not seal properly. Check for leaks before placing in the incubator. Seal with parafilm if necessary. 1h
- 19 Repeat steps 17 and 18 on the other side of the filter.
- 20 Add  210 µL Buffer B3 to the Amicon™ filter unit and vortex to mix.
- 21 Incubate for  00:10:00 min at  70 °C . 10m



- 22 Add  210 μL ethanol (100%) to the AmiconTM filter unit and vortex to mix.
- 23 Proceed with steps 5-7 of the NucleoSpin[®] Tissue Kit Standard protocol for human or animal tissue and cultured cells.
- 24 On step 8 of the NucleoSpin[®] Tissue Kit, perform the following: 2m
- Place the NucleoSpin[®] Tissue Column into a 1.5 mL microcentrifuge tube (not provided in kit)
 - Add  50 μL PCR-grade water
 - Incubate at room temperature for  00:01:00
 - Centrifuge for  00:01:00 min at  11000 x g
 - Repeat steps b, c, & d once, for a final volume of 100 μL