

Dec 21, 2023

Version 4

Visium CytAssist FFPE v4 -- University of Minnesota TMCs V.4

DOI

dx.doi.org/10.17504/protocols.io.n92ldmn99l5b/v4

Laura J Niedernhofer¹, David A Bernlohr¹

¹University of Minnesota, Minneapolis, MN

Cellular Senescence Net...



Allie Abdellatif

MIBAM, University of Minnesota, Minneapolis, MN

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldmn99l5b/v4>

Protocol Citation: Laura J Niedernhofer, David A Bernlohr 2023. Visium CytAssist FFPE v4 -- University of Minnesota TMCs. [protocols.io https://dx.doi.org/10.17504/protocols.io.n92ldmn99l5b/v4](https://dx.doi.org/10.17504/protocols.io.n92ldmn99l5b/v4) Version created by **Allie Abdellatif**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 21, 2023

Last Modified: December 21, 2023

Protocol Integer ID: 92580

Keywords: visium cytassist spatial gene expression for ffpe assay, visium cytassist spatial gene expression, visium spatial gene expression, visium cytassist ffpe v4, mrna in tissue section, gene expression, ffpe assay, tissue sample, visium cytassist instrument, university of minnesota genomics center, minnesota genomics center, mrna, ffpe, tissue section

Funders Acknowledgements:

NIH

Grant ID: 5U54AG079754-02

NIH

Grant ID: 5U54AG076041-03

Abstract


The Visium CytAssist Spatial Gene Expression for FFPE assay is designed to analyze mRNA in tissue sections derived from formalin fixed & paraffin embedded (FFPE) tissue samples. The Visium CytAssist instrument requires a glass slide with intact tissue sections as input.

Tissue sections are processed with Visium Spatial Gene Expression assay performed by the University of Minnesota Imaging Center and sequenced by the University of Minnesota Genomics Center.


Troubleshooting



Deparaffinization, H&E Staining, Imaging & Decrosslinking

- 1  CG000520_Demonstrated_Protocol... 3.1MB

Library Preparation & Sequencing

- 2  CG000495_VisiumCytAssist_GeneE... 31.3MB

Note

Illumina NovaSeq 6000 with read format 28/10/10/50

FASTQ Generation

- 3 BCL data from Illumina sequencer is demultiplexed and converted into FASTQ format using bcl2fastq software version 2.20.0