

Dec 21, 2023

Version 2

Visium Direct Mount FFPE v2 -- University of Minnesota TMCs V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx3qpkg8j/v2>

Laura J Niedernhofer¹, David A Bernlohr¹

¹University of Minnesota, Minneapolis, MN

Cellular Senescence Net...



Allie Abdellatif

MIBAM, University of Minnesota, Minneapolis, MN

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx3qpkg8j/v2>

Protocol Citation: Laura J Niedernhofer, David A Bernlohr 2023. Visium Direct Mount FFPE v2 -- University of Minnesota TMCs. [protocols.io https://dx.doi.org/10.17504/protocols.io.kxygx3qpkg8j/v2](https://dx.doi.org/10.17504/protocols.io.kxygx3qpkg8j/v2) Version created by [Allie Abdellatif](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 20, 2023

Last Modified: December 21, 2023

Protocol Integer ID: 92524

Keywords: visium spatial gene expression for ffpe, visium spatial gene expression, visium direct mount ffpe v2, mrna in tissue section, visium spatial slide with intact tissue section, gene expression, tissue sample, university of minnesota genomics center, tissue section, visium spatial slide, minnesota genomics center, mrna, ffpe, intact tissue section, university of minnesota imaging center

Funders Acknowledgements:

NIH

Grant ID: 5U54AG076041-03

NIH

Grant ID: 5U54AG079754-02

Abstract

The Visium Spatial Gene Expression for FFPE is designed to measure mRNA in tissue sections derived from formalin fixed & paraffin embedded (FFPE) tissue samples and requires a Visium Spatial slide with intact tissue sections as input.

Tissue sections are processed with Visium Spatial Gene Expression assay performed by the University of Minnesota Imaging Center and sequenced by the University of Minnesota Genomics Center.



Deparaffinization, H&E Staining, Imaging & Decrosslinking

- 1  CG000409_VisiumSpatialFFPE_De... 3.1MB

Library Preparation & Sequencing

- 2  CG000407_VisiumSpatialGeneExpr... 6MB

Note

Sequencer: Illumina NovaSeq 6000 using read format 28/10/10/50

FASTQ Generation

- 3 BCL data from Illumina sequencer is demultiplexed and converted into FASTQ format using bcl2fastq software version 2.20.0