



Apr 21, 2023

Version 1

Visualisation and quantification of dendritic spines in cultured human Medium Spiny Neurons (MSNs) V.1

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge397jl47/v1>

Quyen Do^{1,2,3}

¹Oxford Parkinson's Disease Centre and Department of Physiology, Anatomy and Genetics, University of Oxford, South Park Road, Oxford OX1 3QU, United Kingdom;

²Kavli Institute for Neuroscience Discovery, University of Oxford, Dorothy Crowfoot Hodgkin Building, South Park Road, Oxford OX1 3QU, United Kingdom;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA



Cláudia C. Mendes

University of Oxford



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge397jl47/v1>



Protocol Citation: Quyen Do 2023. Visualisation and quantification of dendritic spines in cultured human Medium Spiny Neurons (MSNs). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ge397jl47/v1> Version created by **Cláudia C. Mendes**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 12, 2023

Last Modified: April 21, 2023

Protocol Integer ID: 80391

Keywords: Dendritic spine, Imaris, fluorescent, immunocytochemistry, quantification of dendritic spine, visualisation of dendritic spine, dendritic spine, cultured human medium spiny neuron, human neuron

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-020370

Abstract

This protocol describes the visualisation of dendritic spines of human neurons cultured on coverslips *in vitro* and subsequent quantification using the software Imaris.

Materials

Reagents:

- **[SlowFade™ Diamond Antifade Mountant](#)** (ThermoFisher Scientific, CAT# S36967)

Equipment:

- Olympus FluoView™ FV1000 Confocal Microscope

Software:

- Imaris v9.6.0 (Bitplane, South Windsor, CT, USA)

Visualisation of fluorescent dendritic spines in MSNs

- 1 Mount a coverslip of Medium Spiny Neurons (MSNs) immunolabelled for DARPP32 and neurobiotin onto SlowFade™ Diamond Antifade mountant (follow **Protocol: Immunocytochemistry of cultured human Medium Spiny Neurons (MSNs)**).
- 2 Place slide under under Olympus FluoView FV1000 confocal microscope with argon and solid-state laser with 488 nm and 559 nm excitation, respectively.
- 3 Use the 60x oil-immersion objective (NA = 1.40) to image and capture MSNs coexpressing Darpp32 and Neurobiotin as Z-stacks, sampling sequentially at resolution of 1024 * 1024 pixels and at 1.05 µm steps, as optimised by Nyquist sampling theorem.
- 4 Use the zoom function (3x) in the Olympus Dendritic FluoView FV1000 software to capture dendritic branches of biotinylated neurons.

Note

Approximately 2-4 images are captured per neurons, depending on image quality.

Quantification of dendritic spines

- 5 Use the 'Surface' module to threshold and segment the confocal images.
- 6 Use the 'Filament' module to automatically render the dendrites from images obtained after **step 5**.
- 7 Choose the spine detection function in the 'Filament' module to detect spines as protrusions from the previously-identified dendritic filament.
- 8 Manually exclude putative spines at branch points or disconnected dots.