


Apr 14, 2025

Visualization and Quantification of Intraepidermal Nerve Fibers (IENFs) in Mouse Epidermis via PGP9.5 Labeling and Confocal Microscopy



DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk8n95l5r/v1>

Miguel De Leon¹, Anna Weimer¹, Nicole M. Ashpole¹

¹University of Mississippi



Miguel De Leon

University of Mississippi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk8n95l5r/v1>

Protocol Citation: Miguel De Leon, Anna Weimer, Nicole M. Ashpole 2025. Visualization and Quantification of Intraepidermal Nerve Fibers (IENFs) in Mouse Epidermis via PGP9.5 Labeling and Confocal Microscopy. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.j8nlk8n95l5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 22, 2024

Last Modified: April 14, 2025

Protocol Integer ID: 103873

Keywords: IENF, Intraepidermal nerve fibers, Mice, Anti-PGP 9.5, IENF density , CIPN, Chemotherapy-induced neuropathic pain, quantification of intraepidermal nerve fiber, quantification of intraepidermal nerve fiber density, confocal microscopy intraepidermal nerve fiber, intraepidermal nerve fiber density, intraepidermal nerve fiber, assessment of nerve fiber density, early biomarker for nerve damage, unmyelinated nerve fiber, nerve fiber density, first nerve fiber, mouse epidermis, peripheral neuropathy, dermis into the epidermis, quantification of ienf, nerve damage, epidermis, thermal hyperalgesia, nociception, sensory transduction process, ienf, associated fiber loss, early biomarker, dermi, confocal microscopy for analysis, ienfd, fiber, confocal microscopy

Abstract

Intraepidermal nerve fibers (IENFs) are thin, unmyelinated nerve fibers that extend from the dermis into the epidermis. Consisting mostly of C-fibers and some A δ -fibers, these fibers are involved in sensory transduction processes, including pain (nociception), thermosensation, and itch. Due to these roles, IENFs are among the first nerve fibers to degenerate in peripheral neuropathies, contributing to conditions like mechanical allodynia and thermal hyperalgesia. As a result, IENFs serve as an early biomarker for nerve damage in conditions like diabetes, chemotherapy-induced neuropathy, and autoimmune diseases.

This protocol provides a detailed, step-by-step guide for the processing, staining, and quantification of intraepidermal nerve fiber density (IENFD) in mice. Using an immunohistochemical (IHC) approach, IENFs will be labeled with PGP 9.5, a pan-neuronal marker, and subsequently imaged via confocal microscopy for analysis.

The expected results include clear visualization and quantification of IENFs, allowing for the assessment of nerve fiber density in different experimental conditions. This method enables the detection of neuropathy-associated fiber loss and provides insight into potential neuroprotective or regenerative effects of treatments.

Image Attribution

Images were created in BioRender

Materials

Chemicals			
Chemical Name	Brand	Product #	Lot #
0.1 MPBS (pH 7.4)	Gibco	10010-031	2363510
BSA Powder	Sigma	A2153	SLCH0396
Triton x-100	Fisher Scientific	BP151	110012
Tissue Plus OCT Compound	Fisher Healthcare	4585	4573
PFA 16% w/v soln	Thermo Scientific	043368.9M	N16K030
30% Sucrose	Sigma Aldrich	S9378	308120
Invitrogen Prolong Gold + DAPI	ThermoFisher Scientific	P36931	2841592
Primary Antibody (Affinity Purified Rabbit Anti-human PGP 9.5)	Cedarlane	CL7756	2456258
Secondary Antibody (Invitrogen Alexa-Fluor 488 Goat Anti-Mouse)	ThermoFisher Scientific	A11001	1939600
DAPI	Thermo Scientific	62248	VK3126282
Sodium Citrate Dihydrate	Sigma Aldrich	W302600	MKCF6319
TWEEN 20	Sigma Aldrich	P1379	SLBR9241V
Other Chemicals			
Deionized Water	Ethanol	Dry Ice	

Materials			
Removal of the Hind Paw			
Scissors	1.5 mL MicroCentrifuge Tube	4°C Refrigerator	
Preparation of the Tissue for Slicing			
Ice Bucket	Hammer	Metal Pan	
Dissection of the Foot Pads			
Tweezers	No. 3 Scalpel	No. 10 Lancets	OCT Mold Boxes
Metal Pan (2+ in. deep)	Ice Bucket	- 80°C Freezer	
Slicing the OCT Molds			
Leica CM3050s Cryostat	Ice Bucket	12-well Plate	OCT Stands
Single-Edge Razor	Thick-Combed Paintbrush	Tweezers	
Staining the Paw Sections			
VWR p-20 Micropipette with Corresponding Tips	Gilson p-1000 Micropipette with Corresponding Tips	Microwave	Biomega Hotplate Stirrer
ThermoScientific Heat Block	Fisher Brand Plate Rotator	Aluminum Foil	4°C Refrigerator
Various Glassware for Preparation of Solutions	VWR Model 164AC Balance		
Mounting the Sections on Slides			
10 cm Petri Dish	Microscope Slides	Glass Cover Slips	Gilson p-1000 Micropipette with Corresponding Tips
10/0 Fine-Tipped Paintbrush	VWR p-20 Micropipette with Corresponding Tips		

Troubleshooting

Before start

We, the researchers, utilized a FisherBrand Plate rotator (speed set to 12) in any step that required washing or placement of solution (e.g., reagents, solutions, antibodies). While it is not a required feature of this protocol or other IENF density protocols, we highly recommend it to ensure full and even saturation of the tissues.

Additionally, all wash steps were completed with 0.1 M PBS (pH 7.4). Fiji, a distribution of ImageJ, was used to process confocal images for the quantification of IENFs.

Removal of the Hind Paw

- 1 Before removing the hind paw of the mouse, ensure that you have properly euthanized the mouse, according to IACUC-approved guidelines. Once death has been ensured, use a sharp pair of scissors and remove the entire hind paw at the base of the ankle bone.

Note

When removing hind paws, I tend to alternate between the left and right hind paws of subsequent mice.

- 2 To fix the tissue, place the hind paw in a 1.5 mL microcentrifuge tube (or any 1.5 mL tissue tube) containing a cooled 4% PFA solution and put into a 4°C refrigerator overnight. 

Note

Ensure you use enough PFA to fully submerge the hind paw.

- 3 To sucrose imbed the tissue, remove the hind paw from the PFA solution and wash it in 0.1 M PBS before placing it in a 1.5 mL tissue tube containing 30% sucrose solution. Place tubes in a 4°C refrigerator and wait until tissues have completely sunk to the bottom of the tissue tube before beginning the staining process. Once this occurs, you are ready to prepare the hind paw for slicing. 

Note

The sinking of the tissues generally takes 2-3 days.

Preparation for Tissue Slicing

- 4 Remove the sucrose-embedded tissues from the refrigerator and place tissue tubes in an ice bucket. You are now going to prepare your work bench to flash freeze your hind paws in OCT.
- 5 Acquire a bucket of dry ice and use a hammer (or any tool capable of crushing) to smash the dry ice into fine, snow-like particles.



- 6 Put half of this dry ice into a metal pan that is at least 2 inches deep and leave the remainder in the ice bucket.
- 7 The metal pan will serve as the flash-freezing reservoir for these tissues. In order to activate this reservoir, you must add ethanol to the metal pan containing the dry ice.

Safety information


Be careful not to pour the ethanol too fast for it might cause the dry ice to jump out of the pan and burn you.

- 8 You will now need to remove the hind paws from the sucrose solution and wash them in 0.1 M PBS.
- 9 The following steps will walk you through dissection of the foot pads from the paws. This part of the protocol has been slightly modified from another protocol cited in my references.



Dissection of the Foot Pads

10



Overview of paw dissection. You will position the paw face-up towards the ceiling (A) and use a scalpel or razor to remove the digits (B). You will then use the scalpel on the digit-less foot pad (C) to carefully cut the top layer off (D). When you have removed the top layer (E), put it into an OCT mold (F) and commence the flash-freezing process. *Figure made on Biorender.com*

- 11 Once the paw has been washed, remove it from the PBS and put it on a hard surface. You will need to use a pair of tweezers to hold the paw firmly in place (i.e., by the base of the ankle bone). Using a No. 3 scalpel with an attached No. 10 lancet, remove the toes of the paw at their respective joints, ensuring that you leave the paw pads intact, as shown in parts A-C of the figure .
- 12 Starting from the heel of the paw, gently use the scapel to remove the entire plantar surface of the paw, as shown in parts D and E in the figure. Ensure you are gentle with the tissue to mitigate potential damage. ⚠
- 13 Now that you have cut your tissue, you will need to place them in OCT molds. Ensure that all of your molds are oriented in the same direction to prevent any errors in slicing later. Gently place the plantar surfaces of the paw in a layer of OCT, with the external side facing up and the internal side facing down, as shown in part F of the figure. Additionally, ensure that the tissue is oriented toward the side of the mold of which you plan to begin your slicing. Add another layer of OCT on top of the tissue to fully encase it. ⚠



- 13.1 Using a fine pair of tweezers, move any bubbles that present themselves in the mold away from the tissue sample. You can also use these at this step to fine-tune the position of the tissue to make sure that it is in the middle of the mold and is parallel to the bottom of the mold.
- 14 Set the OCT molds containing tissues in the pan containing the ethanol and dry ice. You will notice that the OCT surrounding the tissue will begin to change from clear to white, starting from the bottom.
- 15 When the mold has turned completely white, remove it from the pan and place it in the ice bucket with the remaining crushed dry ice. Follow steps 10-14 for the remaining paws.
- 16 Deposit your OCT molds in a minus 80°C-safe box and store them in the freezer until you are ready to slice and stain.



Slicing the OCT Molds

- 17 Prior to slicing, ensure that:
 - 1) Your cryostat is turned on and set to -20°C.
 - 2) You place the OCT stands in a bucket of crushed dry ice to cool in preparation for attachment of the molds.
 - 3) You prepare a 12-well plate containing ~1000 μ L of PBS in each well.
- 18 Retrieve your molds from the -80°C freezer. You will now mount your OCT molds to OCT stands. Using a single-edge razor, slice the mold casing off of the OCT block to remove it. For more efficient slicing, also remove any excess OCT surrounding the mold.
- 19 Add a dollop of OCT to the bottom of the mold and press it firmly to the cooled OCT stand. Once the mold has adhered to the stand, encase the bottom of the mold with additional OCT to ensure that it does not fall off during the slicing process.
- 20 Once all molds have properly adhered to their stands, place them in the cryostat for 30 minutes. This will allow the molds to acclimate to the cryostat's temperature.
- 21 Mount the first OCT stand in the cryostat. Beginning at 100 μ m, slice the OCT block until you begin to see the top of the plantar surface. Remove shavings with a thick-combed paint brush.
- 22 Change slice size to 30 μ m and slice 15-30 slices, gently picking up the slices with fine pointed tweezers and placing them in their respective well of the 12-well plate.



30m



23 Repeat steps 16-22 until you have sliced all of your paws.

Staining the Paw Sections

50m

24 Now that you have all of your slices in their respective wells, you can begin the staining process. Start by rinsing the tissue in PBS 3x for 5 minutes each. This will ensure that there is no remaining OCT residue.

15m



25 In the event that you do not start the staining procedure right away, remove the PBS from the wells after all slices have been collected and add 500 μ L of cryopreservative solution to each well, then store your plate in a 4°C refrigerator.



26 **Permeabilization Step:**

Incubate tissues with 0.3% Triton x-100 solution (~500 μ L/well) for 20 minutes, then remove the solution from the well.

20m

Note

No washing is required after this step.

27 **Antigen Retrieval:**

- 1) Boil the citrate buffer in a microwave or on a hot plate.
- 2) Add 500 μ L of the warmed citrate buffer to each well.
- 3) Wrap your 12-well plate in aluminum foil and place it on a heating block set at 95°C for 20 minutes.

20m

28 Remove the plate from the heating block and allow it to cool to room temperature (~ 20 minutes).

20m

29 Once the plate has had time to cool, remove the citrate buffer from each well and rinse the tissues 3 times for 5 minutes each with PBS.

15m



30 **Blocking Step:**

Add 500 μ L of 5% BSA into each well and incubate for 2 hours.

2h

Primary Antibody Incubation

**31 Primary Antibody:**

- 1) Prepare the primary antibody solution using a dilution factor of 1:500 of the PGP 9.5 antibody in 1% BSA solution in 0.3% Triton x-100 solution.
- 2) Remove the 5% BSA from the wells and pipette 500 μ L of the primary antibody solution in each well and allow to incubate overnight in a 4°C refrigerator

**Secondary Antibody Incubation**

2h 20m

32 Secondary Antibody:

- 1) Prepare the secondary antibody solution using a dilution factor of 1:1000 of the Alexa Fluor 488 antibody in 1% BSA solution in 0.3% Triton x-100 solution.
- 2) Remove the primary antibody solution and rinse tissues 3 times for 5 minutes each with PBS.
- 3) Pipette 500 μ L of the secondary antibody solution in each well and allow them to incubate for 2 hours.

2h

**33 DAPI Staining:**

- 1) Prepare the DAPI solution using a dilution factor of 1:1000 of DAPI in PBS.
- 2) Rinse tissues 3 times for 5 minutes each and pipette 500 μ L of DAPI solution in each well for 5 minutes. Tissues will be ready to mount after this.

20m

**Mounting the Sections on Slides**

34 Fill a 10 cm petri dish with PBS.

35 Submerge a microscope slide 1/2 to 1/3 of the way into the dish with the dry side leaning on the rim of the dish.

36 Using a P-1000, transfer the slices from the well that you are beginning with and gently deposit them on top of the submerged slide.

Note

You can cut the end off of the P-1000 pipette tip to avoid tearing the slices.

37 Using a 10/0 fine-tipped paintbrush, gently move each slice onto the slide and adjust their position in a manner that is best for your personal microscope setup.

- 38 When you have oriented all slices onto the slide, remove the slide from the PBS and lay it out of the light to dry.

Note

We have been mounting a total of 8 slices per paw. Two total slides with four slices per slide.

- 39 Repeat steps 36-39 until all slices are mounted onto the slides.

- 40 When slices have completely dried on the slides, you will now mount coverslips with ProLong Gold antifade reagent with DAPI.

Choose the side of the slide that you will begin from (e.g., left-most or right-most). Add a drop of ProLong gold (~10 μ L) to the side of your first tissue.

- 41 Acquire your cover slip and position it at a 45° angle adjacent to where you put the ProLong gold.

- 42 Very gently move your cover slip down towards the slide and allow the slip to spread the ProLong gold evenly across all slices. 

Note

Be careful not to drop the cover slip too fast for it might cause bubbles. Additionally, if you don't allow the slices to properly dry and if you 1) add too much ProLong and/or 2) drop the coverslip too fast, this can cause the slices to move out of position and introduce bubbles to the slide.

- 43 Allow the slides to dry for 24 hours before you begin to image. 

Confocal Microscopy Imaging of IENFs

44 **Microscope Setup**

- Turn on the Leica SP8X confocal microscope and allow the system to initialize.
- Set the Fluo Turret to Scan-BF mode for imaging.
- Ensure that the HC PL APO CS2 40x/1.10 water objective lens is selected and properly immersed in water.



44.1 **Image Acquisition Settings**

- Set the image resolution to 512 × 512 pixels.
- Adjust the scan speed to 400 Hz.
- Configure averaging and accumulation settings:
 - **Line averaging:** 2
 - **Line accumulation:** 1
 - **Frame averaging:** 1
 - **Frame accumulation:** 1

44.2 **Z-Stack Acquisition**

- Enable Z-stack acquisition mode.
- Set the Z-wide mode with a Z-step size of 0.32 μm for optimal sectioning

44.3 **Laser Configuration**

Configure the laser settings as follows:

- 405 nm diode laser
- White Light Laser (WLL) for fluorophore excitation

44.4 **Imaging Execution**

- Adjust exposure and gain settings as necessary for optimal signal detection.
- Capture images with the specified settings, ensuring proper focus and clarity.
- Save images in the desired format for subsequent IENFD quantification.

Image Processing in ImageJ Software

45 **Basement Membrane Length**

- Open ImageJ
- Insert 2D Tiff files

45.1 **Process all Tiff files**

- Split Channels: Image > Color > Split Channels
- Invert Tiff: Edit > Invert
- Close out of all channels but the (blue) wavelength channel



45.2 Set Measurement

1. Analyze > Set Scale
 - Check Global Scale box
 - Click to Remove Scale

Note

The global removal of scale is done to ensure that no previous scale is being applied to the images.

1. Select straight line tool
2. Trace over scale bar on image or if you know the total x or y length of the image in microns you can trace the entire length of the image
3. Push the M key on your keyboard to get the length of pixels.
4. Set Scale:
 - distance in pixels: The length of the previously measured scale bar or x/y distance of image
 - Known distance: Length of the scale bar in microns
 - Pixel aspect ratio: 1.0
 - Unit of length: um
 - Ensure Global is checked
 - Press "Ok"

46 Tracing the Basement Membrane

Now that you only have the blue (DAPI) wavelength image visible, do the following:

1. Select the Freehand line
2. Trace the bottom of the basement membrane, as described below
3. Press "M" to record the length
4. Close out of each image and repeated 1-3 for each subsequent image

1. Identifying the BM

- The basement membrane (BM) is the boundary separating the epidermis and dermis.
- It appears as a thin, distinct interface between the densely packed keratinocytes of the epidermis and the less dense, more fibrous dermal layer.

2. Locating the Correct Tracing Position

- To accurately trace the BM, focus on the bottom edge of the BM, which is closer to the dermal layer.
- Look for a visible separation between: The tightly packed keratinocytes in the lower epidermis. The sparser, more fibrous appearance of the underlying dermis.

3. Step-by-Step Tracing Method

1. Zoom in to enhance visibility of the BM structure.
2. Identify a consistent reference point where the epidermis and dermis separate.
3. Follow the lower edge of this boundary, ensuring the traced line remains just above the dermis but does not extend into epidermal layers.
4. Maintain a smooth, continuous tracing, avoiding jagged or inconsistent movements.
5. If needed, adjust contrast to enhance BM visibility.




Image without any tracing.




Image with the BM outlined using the Freehand Line Tool

Quantification of IENFs

47 Preparing the DAPI Channel for BM Tracing

1. Open the image in ImageJ/Fiji or appropriate imaging software.
2. Select the blue (DAPI) wavelength channel (which stains nuclei and helps visualize the epidermal-dermal boundary).
3. Adjust the contrast and brightness to enhance the visualization of the keratinocyte nuclei, ensuring the BM is clearly distinguishable.
4. Locate the epidermal-dermal junction by identifying the transition between the densely packed basal keratinocytes and the sparser dermal layer.




5. Manually trace the BM along this transition using the freehand or segmented line tool, ensuring a continuous and smooth path.
6. Save the traced BM as a separate overlay or ROI (Region of Interest) for use in the next step.

Overlaying the BM Trace onto the Green Wavelength Channel

1. Switch to the green channel, which contains the PGP 9.5-labeled IENFs.
2. Invert the green channel for better fiber visibility:
 - In ImageJ/Fiji, go to Image → Adjust → Invert LUT.
3. Adjust brightness and contrast to optimize nerve fiber visualization:
 - Increase brightness to enhance faint fibers.
 - Adjust contrast to reduce background noise while maintaining fiber integrity.

Applying the BM Trace to the Green Wavelength Channel

1. Overlay the BM trace onto the green channel:
 - If using ImageJ/Fiji, use ROI Manager to load the previously traced BM.
 - Apply the BM trace to the green channel using "Add Selection" or "Overlay" functions.
1. Ensure that the BM trace correctly aligns with the epidermal-dermal junction in the green channel.
2. Save the final composite image for IENFD quantification.



Inverted green channel with the BM (outlined in yellow) overlaid on top.

Scoring and Quantification of IENFs

48 **Image Processing and Preparation for Scoring**

After processing all images and overlaying the basement membrane (BM), the finalized images were saved and compiled into a PowerPoint presentation for scorers, who were blinded to the treatment groups, to generate their counts. The scoring process followed standardized guidelines adapted from The European Federation of Neurological Societies (EFNS) to ensure consistency and accuracy in IENF quantification.


Counting Criteria for IENFs

Scoring was performed based on the following guidelines:


1. Nerve fibers were counted as they crossed the BM into the epidermis.
2. Any fiber that branches at the BM and enters the epidermis (stratum spinosum) was counted as a separate unit.
3. If a nerve fiber split below the BM before entering the epidermis, each branch was counted as an independent fiber.
4. Nerve fragments that crossed the BM and extended beyond a single layer of basal keratinocytes were included in the count.
5. Fibers that approached the BM but did not cross into the basal keratinocyte layer were excluded from the count.
6. In regions where the BM shifted across different focal planes, only epidermal fibers with a continuous connection to dermal axons were counted.
7. Epidermal axons that lacked a direct connection to dermal axons were excluded from the analysis.

Analyzing IENFD

- 49 Once you have the number of IENFs and the BM length (convert from microns to mm), you can calculate IENFD. Simply divide the **# of fibers/BM length(mm)**.



Left: Image not processed in ImageJ. Right: Yellow line represents the basement membrane. Black arrows point to IENFs that are crossing the basement membrane.



Protocol references

- (1) Krukowski, K.; Ma, J.; Golonzhka, O.; Laumet, G. O.; Gutti, T.; Van Duzer, J.H.; Mazitschek, R.; Jarpe, M. B.; Heijnen, C. J.; Kavelaars, A. HDAC6 Inhibition Effectively Reverses Chemotherapy-Induced Peripheral Neuropathy. *Pain* 2017, 158 (6), 1126–1137. <https://doi.org/10.1097/j.pain.0000000000000893>.
- (2) Ma, J.; Trinh, R. T.; Mahant, I. D.; Peng, B.; Matthias, P.; Heijnen, C. J.; Kavelaars, A. Cell-Specific Role of Histone Deacetylase 6 in Chemotherapy-Induced Mechanical Allodynia and Loss of Intraepidermal Nerve Fibers. *Pain* 2019, 160 (12), 2877–2890. <https://doi.org/10.1097/j.pain.0000000000001667>.
- (3) L. Feldman, E. U Michigan - Intra-Epidermal Fiber Density V2. August 5, 2019. <https://doi.org/10.17504/protocols.io.563g9gn>.
- (4) Lin, X.; Xu, Z.; Carey, L.; Romero, J.; Makriyannis, A.; Hillard, C. J.; Ruggiero, E.; Dockum, M.; Houk, G.; Mackie, K.; Albrecht, P. J.; Rice, F. L.; Hohmann, A. G. A Peripheral CB2 Cannabinoid Receptor Mechanism Suppresses Chemotherapy-Induced Peripheral Neuropathy: Evidence from a CB2 Reporter Mouse. *Pain* 2022, 163 (5), 834–851. <https://doi.org/10.1097/j.pain.0000000000002502>.
- (5) Wai, V.; Roberts, L.; Michaud, J.; Bent, L. R.; Clark, A. L. The Anatomical Distribution of Mechanoreceptors in Mouse Hind Paw Skin and the Influence of Integrin A1 β 1 on Meissner-Like Corpuscle Density in the Footpads. *Front. Neuroanat.* 2021, 15, 628711. <https://doi.org/10.3389/fnana.2021.628711>.
- (6) Wai, V.; Vickaryous, M.; Bent, L.; Clark, A. L. Serial Histological Sectioning of the Intact Mouse Hind Paw for Quantification of the Anatomical Distribution of Meissner-like Corpuscles across the Skin. *J. Histol. Histopathol.* 2020, 7 (1), 1. <https://doi.org/10.7243/2055-091X-7-1>.
- (7) Beiswenger, K. K.; Calcutt, N. A.; Mizisin, A. P. Epidermal Nerve Fiber Quantification in the Assessment of Diabetic Neuropathy. *Acta Histochem.* 2008, 110 (5), 351–362. <https://doi.org/10.1016/j.acthis.2007.12.004>.

Acknowledgements

All imaging data were captured and processed using LAS X software version 3.5.7.23225 provided by the core grant support of the University of Mississippi Glycoscience Center of Research Excellence NIH P20GM103460.