

Jul 28, 2022

Version 2

Visualization of a low concentration and molecular weight DNA (DNA gel stain) V.2

DOI

dx.doi.org/10.17504/protocols.io.5jyl89x46v2w/v2

Nadine Mowoh¹, Stephane Fadanka¹, Mugar Minette Shalo¹


¹Beneficial Bio, Mboalab

Beneficial Bio

Low-cost, high-quality ...



Nadine Mowoh



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl89x46v2w/v2>

Protocol Citation: Nadine Mowoh, Stephane Fadanka, Mugar Minette Shalo 2022. Visualization of a low concentration and molecular weight DNA (DNA gel stain). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl89x46v2w/v2> Version created by **Nadine Mowoh**



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 27, 2022

Last Modified: July 28, 2022

Protocol Integer ID: 67728

Keywords: Quality control test of DNA gel stain, low concentration dna with agarose gel electrophoresi, agarose gel electrophoresi, dna gel stain, based dna gel stain, gel electrophoresi, benbio dna gel stain, use of benbio dna gel, agarose gel, benbio dna gel, low concentration dna, concentration of agarose, sizes of the dna fragment, repeated agarobiose, dna fragment, agarose, molecular weight dna, dna, biomolecule, seaweed genera gelidium, pcr amplification of dna, most gel, gel, charged biomolecule, biomolecules such as rna, low molecular weight, visualization of low molecular weight, protein, molecular weight, subunit

Abstract

After PCR amplification of DNA, agarose gel electrophoresis is run to separate the DNA based on their size. The agarose gel consists of microscopic pores that act as a molecular sieve which separates molecules based on their charge, size and shape. Agarose gel electrophoresis can also be used to separate other charged biomolecules such as RNA and proteins. Agarose is isolated from the seaweed genera *Gelidium* and *Gracilaria* and consists of repeated agarobiose (L- and D-galactose) subunits. The concentration of agarose in a gel depends on the sizes of the DNA fragments to be separated, with most gels ranging between 0.5%-2%.

This protocol describes the use of BenBio DNA gel stain for the visualization of low molecular weight and low concentration DNA with agarose gel electrophoresis.

The BenBio DNA gel stain is therefore, a good alternative to EtBr-based DNA gel stains.

Guidelines

This procedure can be performed by laboratory staff that have been trained and have theoretical and practical skills in good laboratory practices. It can also be performed by molecular biology students or students of related fields under the supervision of a laboratory staff.



Materials

Reagents

- 10000x TO-DMSO gel stain stock (test sample)
- DNA amplicon (*in this case 50 or 100bp DNA ladder*)
- 10x TBE buffer
- EtBr based gel stain (positive control or standard)

Materials and equipment

- Glass beaker
- Micropipette and tips
- Micropipette and tips
- Microwave
- Gel casting tray
- Well comb
- UV transilluminator
- Voltage source (Electrophoresis unit)

Troubleshooting

Safety warnings

- ! ▪ If using UV transilluminator to visualize, care must be taken in order not to expose yourself to UV rays.
- A safer alternative to use is the blue light transilluminator.
- Care should also be taken when handling EtBr as it is a know mutagen.

Before start

Ensure all materials, equipment and chemicals to be used for this experiment are all in place before starting.

Visualization of a low concentration and molecular weight DNA

1

Note

The test for visualization of a low concentration and molecular weight DNA are applicable for DNA gel stains, Electrophoresis buffers and DNA loading dye. The visualization tests could be modified slightly to suit a particular product.

The protocol described below describes the QC steps in confirming a product's ability to allow visualization of low molecular weight DNA and low concentration DNA.

2 To confirm visualization of a low molecular weight DNA

- Follow the steps in preparing a 2% agarose gel as described [in this protocol](#),

Note

We use the TO-DMSO DNA gel stain as "test DNA gel stain" and an EtBr based DNA gel stain as "positive control or standard DNA gel stain" (*meaning 2 gels will be prepared and run simultaneously, where possible.*

If this is not possible, run one gel to finish and run the next.

- Load the gel (made from the test gel stain) by pipetting 3 to 5 uL of sample into the wells (*in this case a loading dye would not be used since the amplicon/ sample is a DNA ladder that is made with the tracking dye in it.*
- Load the second gel (*made from the control gel stain*) by pipetting 3 to 5 uL of sample into the wells (*in this case a loading dye would not be used since the amplicon/ sample is a DNA ladder that is made with the tracking dye in it).*
- Connect the electrophoresis units to a power pack and run the 2 gels simultaneously under the same conditions (*when using TBE running buffer, we run at 80 to 100 Volts for 15 to 20 minutes*) to finish.
- Transfer the gel unto a UV transilluminator and visualize the resulting band separation to confirm separation of the ladder into visible distinct bands (from the lowest to highest base pair bands).



To confirm visualization of a low concentration DNA

- Make 1:5 serial dilutions of the 100bp DNA ladder or any DNA amplicon and use the different dilutions as amplicons to load a gel well.

Note

- 2 gels would be prepared, one with the test gel stain and another with the control gel stain.
 - These gels could be run simultaneously or one after the other.
-
- Repeat the steps described above in loading a gel using the different dilutions of the DNA.
 - Connect the electrophoretic unit to the power source and allow the gel to run through (when using TBE buffer, we run at 80 to 100 Volts for 15 to 20 minutes).
 - Transfer the gel unto a UV transilluminator and visualize to identify visible, sharp bands which determine the highest dilution of the DNA that is detectable by the gel stain.