


Jan 11, 2024

Visualization of LRRK2 filaments in 293T cells

 [The Journal of Biological Chemistry](#)

DOI

dx.doi.org/10.17504/protocols.io.8epv5xpz4g1b/v1



Eva Karasmanis¹, Kathryn S Hatch¹

¹University of California, San Diego



Eva Karasmanis

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5xpz4g1b/v1>

Protocol Citation: Eva Karasmanis, Kathryn S Hatch 2024. Visualization of LRRK2 filaments in 293T cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5xpz4g1b/v1>

Manuscript citation:

doi: <https://doi.org/10.1101/2023.11.14.567123>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 11, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 93405

Keywords: ASAPCRN, visualization of Irrk2 filament, Irrk2 filaments in the presence, Irrk2 filament, cell, percentage of cell, Irrk2, visualization

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000519

Abstract

Visualization of LRRK2 filaments in 293T cells

GOAL: Express GFP-LRRK2 with or without DARPin E11 and quantify the percentage of cells with LRRK2 filaments in the presence and absence of MLi-2 in 293T cells.

Image Attribution

Eva Karasmanis

Troubleshooting



Day 1: fibronectin coating and cell plating













1

Visualization of LRRK2 filaments in 293T cells

GOAL: Express GFP-LRRK2 with or without DARPin E11 and quantify the percentage of cells with LRRK2 filaments in the presence and absence of MLI-2 in 293T cells.

Constructs needed:

- 1) CMV-8xHISDaprin C12-FLAG
- 2) CMV- His-daprin E11-FLAG
- 3) CMV-GFP-LRRK2

	DMSO	MLi-2	DMSO	MLi-2
Replica te 1				
Replica te 2				
Replica te 3				

2 Fibronectin Coating (Sigma cat# F0895, 0.1% solution, 1 mg/ml):

- 2.1 Make 0.01 µg/µL solution of fibronectin, stock @ 1 mg/mL For 6× 35 mm dishes (12 ml fibronectin working stock)- 0.12 mL fibronectin + 11.88 mL 1X PBS
- 2.2 Lay one 22 mm x 22 mm glass coverslip into each 35 mm dish or 6 well.
- 2.3 Add 2 ml of 10 ug/mL fibronectin per 35 mm dish or 6 well.



2.4 Incubate at 37 °C 5% CO₂ for 00:45:00

45m

2.5 Wash with PBS and let dry for 00:45:00 in the tissue culture hood (no UV)

45m

3 Plate cells onto Fibronectin coated dishes

3.1 Dissociate cells, count and plate 6 well plate with 200K cells /well. For transfection, plate in antibiotic-free media (DMEM+10% FBS)

Day 2 Transfect cells with GFP-LRRK2 and Darpins:

2d 3h 2m

4 Transfect 800 ng LRRK2 and 400ng Darpin with 5 µL PEI /well.
- cells should be 50-60% confluent

In a sterile tube dilute 800 ng GFP-LRRK2 and 400 ng DARPIn E11 plasmid DNA in Optimem (150 µL /well)

Per Well:

a) 150 µL Optimem (prewarmed)+ 5 µL PEI

b) 150 µL Optimem (prewarmed) + 800 ng GFP LRRK2 and 400 ng DARPIn E11 DNA.

or

c) 150 µL Optimem (prewarmed) + 800 ng GFP LRRK2

Add PEI to diluted DNA - (Per reaction: 5 µL of 1 µg/µL stock). Mix immediately.

4.1 In a sterile tube dilute 800 ng GFP-LRRK2 and 400 ng DARPIn E11 plasmid DNA in Optimem (150uL/well)

Per Well:

a) 150 µL Optimem (prewarmed)+ 5 µL PEI

b) 150 µL Optimem (prewarmed) + 800 ng GFP LRRK2 and 400 ng DARPIn E11 DNA.

or



c) 150 μ L Optimem (prewarmed) + 800 ng GFP LRRK2

Add PEI to diluted DNA - (Per reaction: 5 μ L of 1 μ g/ μ L stock). Mix immediately.

4.2 Incubate at Room temperature for 00:15:00

15m

4.3 Add DNA/PEI mixture to cells dropwise.

4.4 Swirl the plate to distribute Incubate at 37 $^{\circ}$ C 5% CO₂ for 48:00:00

2d

Day 4 MLi2 or DMSO treatment

2d 3h 2m

5 Treat Cells with MLi2 or DMSO (5 μ L per well of a 1 mM stock) for 02:00:00 ,
 37 $^{\circ}$ C 5% CO₂

2h

Fixing and Staining of Cells

2d 3h 2m

6 **Fixing:** Prewarm freshly made 3% PFA, 4% sucrose in 1X in PBS 1X. You will need 1 mL per well.

6.1 Aspirate media

6.2 Immediately add prewarmed fixation buffer (3% sucrose, 4 % (v/v) PFA in PBS 1X).

7 Incubate 00:12:00 Room temperature

12m

8 2X rinse with PBS 1x

9 2X Wash with PBS1x for 00:05:00 at Room temperature

5m



- 10 2X Wash with PBS 1x for 00:05:00 at Room temperature . 5m
- 11 **Quenching:** 2X rinse and 2x 00:10:00 of 0.4% NH₄Cl [75 millimolar (mM)] in PBS 1X 10m
- 12 **Blocking/Permeabilizing:** incubate cells blocking/permeabilizing buffer (2% BSA + 0.1 % triton X-100 in PBS 1X) for 00:20:00 Room temperature 30m
- During blocking make and spin the antibody (ab) mix. Each coverslip needs ~40-50 uL uL of ab mix. Antibodies are diluted in Blocking buffer (2% BSA in PBS 1X) and spun at 4 °C 00:10:00 m 50000 x g, 4°C to clear aggregates. Remove spun ab and place in new tube. Mix to get even concentration. HERE: we used 1:200 rabbit polyclonal anti-FLAG DYKDDDDK tag – (ptg labs Cat no : 20543-1-AP)
- 13 2X Rinse with Wash 2X (00:05:00 , Room temperature) with blocking buffer 5m
- 14 Add primary antibody mix : 3h 5m
Add 40uL on parafilm and flip coverslip on the antibody. Incubate 03:00:00
 Room temperature or 4 °C Overnight . If incubating overnight, make a humidity chamber before placing in fridge.
- HERE: we used 1:200 rabbit polyclonal anti-FLAG DYKDDDDK tag – (ptg labs Cat no : 20543-1-AP)
- 15 Rinse 2X with blocking buffer
- 16 Wash 2X 00:05:00 Room temperature with blocking buffer 5m
- 17 Add secondary mix in PBS (1:200 Goat a-rabbit Alexa568; SIGMA A-11011 +1:5000 DAPI**) 00:45:00 Room temperature 45m
** DAPI can alternatively be added as a separate step at 1:1000 dilution for 15 min
- 18 5X Rinse with PBS 1X



- 19 Mount in Fluorsave hard media (Millipore 345789)
- 20 Let dry for at least an hour. Store in Fridge 4oC if not imaging immediately. Check coverslip is set with tweezers before imaging.

Imaging and analysis

- 21 Blind your mounted slides before imaging to prevent bias during aquisition.
- 22 Find areas with transfected cells. Acquire Z stacks by determining top and bottom with a 0.3 um step size. (about 20-25 z stacks) ·
- 23 **Analysis in Fiji:** Go through each image, make max projections, and mark each transfected cells with a ROI (region of interest).
- 24 In an excel sheet keep track of each cell you mark and score as 0 if no LRRK2 filaments are present or 1 if some are. ·
- 25 Include at least 50 cells/ sample. ·
- 26 Calculate % cells with filaments (number of transfected cells with filaments /total number of transfected cells *100) ·
- 27 Unblind
- 28 Transfer values to prism to generate graph and statistics.