

Aug 08, 2022

Washing Protocol for Intact Proteoform MALDI on Human Kidney

DOI

dx.doi.org/10.17504/protocols.io.dm6gpb7qdlzp/v1

Kevin J Zemaitis¹, Dusan Velickovic¹, Ljiljana.PasaTolic¹

¹Pacific Northwest National Laboratory

Human BioMolecular Atl...

PNNL-TTD



Kevin J Zemaitis

Pacific Northwest National Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpb7qdlzp/v1>

Protocol Citation: Kevin J Zemaitis, Dusan Velickovic, Ljiljana.PasaTolic 2022. Washing Protocol for Intact Proteoform MALDI on Human Kidney. [protocols.io https://dx.doi.org/10.17504/protocols.io.dm6gpb7qdlzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gpb7qdlzp/v1)

Manuscript citation:

DOI: 10.1021/acs.analchem.2c01034

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 07, 2022

Last Modified: August 08, 2022

Protocol Integer ID: 68337

Keywords: washing protocol for intact proteoform maldi, intact proteoform maldi, human kidney scope, sample washing protocol, human kidney, washing protocol, several timed washing step, kidney, spectroglyph ep, maldi

Funders Acknowledgements:

National Institutes of Health (NIH) Common Fund, Human Biomolecular Atlas Program (HuBMAP)

Grant ID: UG3CA256959-01

Abstract

Scope:

A detailed protocol entailing the sample washing protocols developed for human kidney, this includes several timed washing steps for fixing the tissue, removing small molecule interferents, and desalting the tissue surface optimized for use with a Spectrogyph EP-MALDI-2 source mounted on a custom Thermo Scientific UHMR Q Exactive HF Orbitrap.

Expected Outcomes:

Human kidney slides ready for pre-extraction and/or matrix application.

Materials

Chemicals:

- Ethanol (200 Proof)
- Nanopure water (LC-MS Grade)
- Chloroform (LC-MS Grade)
- Acetic acid (Glacial, LC-MS Grade)
- Trifluoroacetic acid (LC-MS Grade)

Equipment:

- Regulated nitrogen gas supply with nozzle
- Coplin jars


Troubleshooting



Safety warnings



⚠ All steps of this protocol working with solvents should be performed within a fume hood as to minimize exposure to fumes from volatile organic solvents.

Before start

Prepare  100 mL of all necessary solvents for the tissue washing steps and clean Coplin jars prior to use.







Preparation

- 1 While the tissue is within the vacuum desiccator, pour all necessary solvents and solutions into the Coplin jar. These have a volume of approximately  75 mL and to reduce the variability in extraction fill the jars consistently.
- 2 Prior to submerging the tissue section in solvent and solutions, ensure that the timers are set beforehand. While three timers are not necessary, it vastly improves the process.
- 3 If any metadata is present on the slide written in permanent marker remove it at this time. Take a photo and note within a notebook prior to removal, this is necessary as this will contaminate the Coplin jars within all washes. 






Tissue fixation wash

3m

- 4 Submerge the tissue section within a Coplin jar filled with a solution of **70% ethanol** for  00:00:30 
- 5 Immediately move the tissue section with tweezers to then next coplin jar containing **100% ethanol** for  00:00:30 

Lipid depletion wash


3m

- 6 After completion of the tissue fixation steps, immediately move the tissue section with tweezers to then next Coplin jar containing **Carnoy's solution** for  00:02:00 
- 6.1 **Carnoy's solution** is a 6:3:1 volumetric solution of ethanol: chloroform: glacial acetic acid, while other solutions could be put within plastic staining jars a glass Coplin jar is highly recommended for this and other polar aprotic solvents.
- 6.2 Note: dipping, agitation, and any movement of the slide within the Coplin jar will dramatically change outcomes of all washes. Tissue is prone to detach from the surface within this step. 
- 7 Following exposure to Carnoy's solution, submerge the slide with **100% ethanol** for  00:00:30 

Desalting wash

3m




8 After completion of the lipid depletion washes, submerge the tissue sections within **nanopure water with 0.2% trifluoroacetic (TFA) acid** for  00:00:15

15s

8.1 Note: dipping, agitation, and any movement of the slide within the Coplin jar will dramatically change outcomes of all washes. Tissue is prone to detach from the surface within this step.




9 After exposure to water, submerge the slide within **100% ethanol** for  00:00:30

30s

9.1 Ensure this is a new solution of **100% ethanol** not used within previous steps.

Drying tissue

3m

10 After the completion of all steps, dry the tissue for  00:00:30 under a diffuse stream of nitrogen, take care to ensure that liquid does not pool within the surface. The tissue section is now prepared for the next steps within the tissue preparation protocol.

30s

11 Discard all solvents and solutions properly, and best practice is to use fresh solvents and solutions for subsequent washes.