

Jun 26, 2020

Water Filtration with Peristaltic Pump or Syringe & DNA Preservation

DOI

dx.doi.org/10.17504/protocols.io.bhxyj7pw

Annette S Engel¹, Audrey Paterson¹, Hannah Rigoni¹

¹University of Tennessee, Knoxville

Engel Lab at Tennessee



Annette S Engel

University of Tennessee, Knoxville



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bhxyj7pw

Protocol Citation: Annette S Engel, Audrey Paterson, Hannah Rigoni 2020. Water Filtration with Peristaltic Pump or Syringe & DNA Preservation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bhxyj7pw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our group and it is working.

Created: June 26, 2020

Last Modified: July 10, 2020

Protocol Integer ID: 38616

Keywords: field sampling, filtration, water, eDNA,



Abstract

Cells, particulates, and fine sediment (i.e., clay) in water are retained on filters during filtration. Filtered water can be collected for geochemical analyses and nucleic acids (referred to as environmental DNA or eDNA) can be extracted from the filters. This protocol describes water filtration with a 0.22 mm Sterivex (Millipore) filter cartridge by using a peristaltic pump or manually with a syringe, and the preservation of filtered material collected on the filter membrane using CTAB.

Attachments



UTK Engel SOP Water ...

484KB

Image Attribution

AS Engel

