



Apr 02, 2023

Western blot - alpha-synuclein

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l69pqdlqe/v1>

Michael J Hurley^{1,2}

¹Department of Clinical and Movement Neurosciences, UCL Queen Square Institute of Neurology, London, UK;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland, USA.



Michael J Hurley

UCL

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l69pqdlqe/v1>

Protocol Citation: Michael J Hurley 2023. Western blot - alpha-synuclein. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bp2l69pqdlqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: February 27, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 77696

Keywords: ASAPCRN, alpha-synuclein, western blot, synuclein in protein, synuclein this protocol, synuclein, peroxidase, protein

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000420

Abstract

This protocol describes how to detect alpha-synuclein in protein derived from STC-1 cells by western blot using DAB/peroxidase or ECL to visualise the bands.

Materials

RIPA buffer (50 mM Tris-HCl, pH 8.0, with 150 mM sodium chloride, 1.0% Igepal CA-630 (NP-40), 0.5% sodium deoxycholate, and 0.1% sodium dodecyl sulphate)

Protease inhibitors (0.8 μ M Aprotinin, 40 μ M Bestatin, 140 μ M E-64 at, Leupeptin at 20 μ M and Pepstatin A at 15 μ M, 1mM phenylmethanesulfonyl fluoride)

Towbin transfer buffer: 25 mM Tris, 192 mM glycine, pH 8.3, with 20% methanol (vol/vol).

To prepare 1 L of buffer:

- 800 mL distilled H₂O
- 200 mL methanol
- 3.03 g Tris base
- 14.4 g Glycine



Sample preparation

1d 4h 16m

- 1 *Sample preparation using either protein pellet from an extraction (step 1.1) **or** direct lysis (step 1.2)*
- 1.1 Dissolve protein pellets from the TRI Reagent® extraction in 2% SDS containing 8 M urea. 30m 30s
Vortex 00:00:30 twice, then orbital shaker 260 rpm 00:30:00
- 1.2 Dissolve cells in RIPA buffer containing protease inhibitors (see Materials) On ice 30m
 00:30:00
- 2 Mix samples with 4x NuPAGE™ LDS sample buffer containing 50 mM dithiothreitol 15m 20s
Vortex 00:00:20 twice, then orbital shaker 260 rpm 00:15:00
- 3 Heat samples 70 °C 00:10:00 10m
- 4 Centrifuge samples 00:00:30 to collect condensation from tube lid. Vortex 00:00:20 then centrifuge again 00:01:00 1m 50s


Electrophoresis and transfer

1d 4h 16m

- 5 Separate samples (~ 10 ug total protein) by polyacrylamide gel electrophoresis using precast Bolt™ or NuPAGE™ 12%, Bis-Tris, 1.0 mm, Mini Protein Gels at 180 V 45m
 00:45:00 or until dye front reaches the bottom of the gel. Run with pre-stained or biotinylated size markers
- 6 Wet PVDF membrane with methanol and then equilibrate in Towbin buffer (see materials)
- 7 Soak transfer sandwich components (4 sheets of filter paper and 5 blotting pads) in Towbin transfer buffer and assemble in the transfer cassette in the following order:
Starting with the cathode plate
2 x blotting pads
2 x filter paper

gel
PVDF
2 x filter paper
3 x pad (use 3 because there is less chance of the gel slipping)
Use a roller to remove any air bubbles

anode



2X Blotting pads
Filter paper
Membrane
Gel
Filter paper
2X Blotting pads

cathode

8 Place cassette in transfer tank and transfer protein onto 0.2 mm PVDF (polyvinylidene difluoride) membranes (30 V, ⌚ 01:00:00 in Towbin buffer

1h

9 Rinse membranes with PBS ⌚ 00:00:30

30s

Immunodetection of protein bands

1d 1h 5m 30s

10 Fix proteins to membrane for ⌚ 00:20:00 with 4% paraformaldehyde

20m

11 Wash with PBST (PBS containing 0.01% (v/v) Tween™-20) 4 x ⌚ 00:05:00

5m

12 Block non-specific binding sites with block solution (PBST containing 2% (w/v) BSA and 0.005 % (w/v) thiomersal) for ⌚ 01:00:00

1h



- 13 Incubate membranes in block solution containing rabbit monoclonal antibody against α -synuclein (1:1250) (ab212184) and β -actin (1:5000) (ab241153) ⌚ 20:00:00 20h
🌡 Room temperature
- 14 Wash with PBST 4 x ⌚ 00:05:00 5m
- 15 *Incubate membrane with secondary antibody in block solution (step 15.1 **or** step 15.2)*
- 15.1 Incubate membrane with biotinylated goat anti-rabbit IgG (1:1000) (Stratech) ⌚ 01:15:00 and go to **step 16** 1h 15m
- 15.2 **OR** incubate membrane with goat anti-rabbit conjugated to peroxidase (1:25,000) (Stratech) ⌚ 01:15:00 and go to **step 18** 1h 15m
- 16 Wash with PBST 4 x ⌚ 00:05:00 5m
- 17 incubated with peroxidase conjugated streptavidin (1:1000) (Roche) ⌚ 00:45:00 45m
- 18 Wash with PBST 4 x ⌚ 00:05:00 5m
- 19 Wash with PBS 2 x ⌚ 00:05:00 5m
- 20 *Visualise bands with DAB if biotinylated secondary and streptavidin-peroxidase were used (step 20.1) **or** use ECL plus (Amersham) if anti-rabbit peroxidase secondary was used (step 20.4)*
- 20.1 **For DAB detection** incubate in PBS containing 0.05 % (w/v) 3,3'-diaminobenzidine, 0.015 % (v/v) H_2O_2 and 0.05 % (w/v) nickel ammonium sulphate for $\leq$ ⌚ 00:00:30 30s
- 20.2 Rinse in running tap water > ⌚ 00:05:00 5m



20.3 Dry membrane 60°C 1 hour

20.4 **For chemiluminescent detection** use an Amersham ECL plus kit according to the manufacturer's instructions

21 Visualize ECL or digitize DAB stained membrane using a ChemiDoc™ MP imaging system (Biorad)

22 ECL reagents can be washed off the membrane and the membrane then processed for DAB staining if desired