

Jun 13, 2023

Version 1

Western Blot Analysis V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkwp25l5r/v1>

michela.deleidi^{1,2}, Bianca Marchetti^{3,4}, Federico Bertoli^{1,2}, Carmela Giachino³

¹Mitochondria and Inflammation in Neurodegenerative Diseases, DZNE, Tübingen-Germany;

²Hertie Institute for Clinical Brain Research, University of Tübingen;

³Neuropharmacology Laboratory, Oasi Research Institute-IRCCS, Troina, Italy;

⁴Biomedical and Biotechnological Sciences, Pharmacology Section, University of Catania-Italy



carmela giachino

institut imagine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkwp25l5r/v1>

Protocol Citation: michela.deleidi , Bianca Marchetti, Federico Bertoli, Carmela Giachino 2023. Western Blot Analysis. protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlkwp25l5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 31, 2023

Last Modified: June 13, 2023

Protocol Integer ID: 79825

Keywords: Western Blot Analysis, western blot analysis western blotting, western blot analysis, western blot, antibodies with fluorescent, antibody, using antibody, technique for the immunodetection, immunodetection, chemiluminescent detection, protein, fluorescent

Funders Acknowledgements:

ASAP-Aligning Science Across Parkinson's

Grant ID: ASAP-000420

Abstract

Western Blotting is a technique for the immunodetection of proteins using antibodies with fluorescent or chemiluminescent detection.

Attachments



[k4dwb48if.docx](#)

14KB

Materials

Materials

- 0.33 M sucrose
- 8 mM Hepes, pH 7.4
- Laemli SDS boiling buffer (Sigma)
- 9-12% SDS-polyacrilamide gel
- polyvinylidene difluoride membrane (Amersham Biosciences, Piscataway, NJ)
- 5% non-fat dry milk
- TBST



Western Blot Analysis

4h 40m

- 1 Prepare protein extracts as previously described.
- 2 Homogenize tissues in lysis buffer ([M] 0.33 Mass Percent sucrose/ [M] 8 millimolar (mM) Hepes,  7.4 and protease inhibitors) and quantify them using the BCA protein determination method (Bio-Rad, Hercules, CA).
- 3 Dilute protein samples to equivalent volumes containing  20 μ g of protein and boil in an equal volume of Laemli SDS boiling buffer (Sigma) for  00:10:00 . 10m
- 4 Load samples into a 9-12% SDS-polyacrilamide gel and separate it by electrophoresis for  03:00:00 at 100 V. 3h
- 5 Transfer the proteins to polyvinylidene difluoride membrane (Amersham Biosciences, Piscataway, NJ) for  01:30:00 at 300 mA. 1h 30m
- 6 After blocking of nonspecific binding with 5% non-fat dry milk in TBST, probe the membranes with primary antibodies and process.
- 7 Perform densitometric analysis using ImageQuantity One. 
- 8 Normalize data to β -actin, normalize values of phosphorylated GSK-3 β (pTyr²¹⁶ GSK-3 β); phosphorylated α -syn (pSer¹²⁹ α -syn) and phosphorylated tau (pSer³⁹⁶ tau) to total GSK-3 β , α -syn, and tau, respectively, before statistical analysis of variance and express values as percent changes (%) of WT controls.
- 9 Dashed lines (in white) indicate discontinuous bands (nonsequential lanes) taken from the same blot, at the same molecular weight (mass – kDa) in order to better represent the mean signal from all values (5-6 individual blots/genotype/treatment) for that particular group. Corresponding control bands (loading controls) match experimental bands.