

Nov 08, 2022

Western Blot

DOI

dx.doi.org/10.17504/protocols.io.q26g7y428gwz/v1

Elena Coccia¹, gustavo.parfitt Parfitt¹

¹Icahn School of Medicine at Mount Sinai

Ahfeldt Lab



Elena Coccia

Icahn School of Medicine at Mount Sinai

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7y428gwz/v1>

Protocol Citation: Elena Coccia, gustavo.parfitt Parfitt 2022. Western Blot. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7y428gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 27, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 71931

Keywords: ASAPCRN, western blot, protocol

Abstract

Western Blot Protocol

Troubleshooting



Cell lysis preparation

50m

- 1 Prepare the appropriate amount of RIPA buffer with Protease and Phosphatases Inhibitor Cocktail (Roche) and maintain On ice
- 2 Collect cells/organoids in PBS On ice in an eppendorf and spin down 500 rpm, 4°C, 00:05:00 5m
- 3 Add appropriate amount of lysis buffers (1:50 v/v proportion between cell pellet and lysis buffer) to the pellet and pipette up and down
- 3.1 For organoid samples use an homogeniser
- 4 Incubate samples for 00:30:00 On ice . Spinning every 00:10:00 40m
- 5 Centrifuge at 15000 x g for 00:10:00 at 4 °C and collect the supernatant in a new tube. 10m
- 6 If SDS fraction is desired, add 2 Mass / % volume SDS to RIPA buffer and resuspend pellet from step 5.
- 6.1 Sonicate the sample a 3 times for 00:00:10 to ensure fragmentation of genomic DNA released during lysis and consequently reduce viscosity 10s
- 7 Quantify protein concentration.

SDS-PAGE

- 8 Mix 20-30 µg of protein with Laemmli's loading buffer



- 9 Denature the proteins by incubation at 95 °C for 00:05:00 , spin briefly before loading 5m
- 10 Load pre-cast gel into Western Blot apparatus and fill with Running Buffer (BioRad).
- 11 Load samples and protein ladder into gels, Run the gel at 120V for 00:45:00 to ensure protein separation. 45m

Protein Immunodetection

- 12 Transfer proteins into a PVDF membrane (previously activated with methanol and hydrated in ddH₂O)
- 13 Remove membrane from transfer and place into a box with blocking buffer: 5% BSA in TBS-T (20mM Tris-HCl, 150mM NaCl pH8, 0.1% Tween20). Block for 01:00:00 1h
- 13.1 If alpha-synuclein is to be blotted, fix the membrane with 4% PFA in PBS for 00:30:00 prior to blocking 30m
- 14 Once blocked, sequentially probe the membrane for antibody staining and detection. Antibody dilution might need optimization.
- 14.1 Prepare primary antibody in 5% BSA in TBS-T and left incubating Overnight at 4 °C 30m
- 14.2 Wash membranes 3x 00:10:00 in TBS-T and then incubated for 01:00:00 at Room temperature with the appropriate HRP-tagged secondary antibody. 1h 10m
- 14.3 Membrane was washed 3x 00:10:00 in TBS-T 10m
- 15 Mix EZ-ECL solution 1:1 (v/v) and incubate on the membrane for 00:01:00 1m



16 Image membranes using a Licor Odyssey Chemiluminescence Imager.