

Feb 03, 2025

Western blot

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk94x5v5r/v1>

Maria Jose Perez J.¹

¹Institut Imagine

Team Deleidi



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk94x5v5r/v1>

Protocol Citation: Maria Jose Perez J. 2025. Western blot. protocols.io
<https://dx.doi.org/10.17504/protocols.io.j8nlk94x5v5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 26, 2024

Last Modified: February 03, 2025

Protocol Integer ID: 112825

Keywords: western blot

Funders Acknowledgements:

ASAP

Abstract

Western blot

For cell experiments

- 1 Extract proteins using 0.5% NP40–PBS protein extraction buffer containing protease and phosphatase inhibitors (Roche, Switzerland) on ice.
- 2 Centrifuge the suspension at 14,000 rpm for 15 minutes at 4°C

For mouse tissue experiments

- 3 Extract and clean tissue with PBS.
- 4 Dissect the cortex quickly to prevent protein degradation.
- 5 Add ice-cold homogenization buffer (50 mM Tris-HCl pH 7.4, 150 mM NaCl, 1% Triton X-100, 0.1% SDS, 1 mM EDTA, and protease and phosphatase inhibitors) at approximately 1 mL per gram of tissue.
- 6 Homogenize the tissue manually with a pestle for 10–20 passes on ice.
- 7 Centrifuge the homogenate at 14,000 rpm for 20 minutes at 4°C.

After centrifugation for both cell and mouse tissue experiments

- 8 Measure protein concentration in the supernatant using the PierceTM BCA Protein Assay Kit (Cat. number 23225, Thermo Fisher Scientific).
- 9 Load 20–50 µg of protein lysate onto a polyacrylamide gel (SurePAGETM, Bis-Tris, 10×8, 4%–12%, GenScript).
- 10 Separate proteins via electrophoresis and transfer them onto a 0.45 µm polyvinylidene fluoride (PVDF) membrane (Millipore).



- 11 Block the membrane with either 5% milk powder or 5% bovine serum albumin (BSA) in TBS containing 0.1% Tween 20 (TBST).
- 12 Incubate the membrane overnight at 4°C with primary antibodies diluted in the chosen blocking solution.
- 13 Incubate the membrane with horseradish peroxidase (HRP)-conjugated secondary antibodies (Sigma–Aldrich) for 1 hour at room temperature.
- 14 Visualize the proteins using ECL Western Blotting Detection Reagent (Pierce, Thermo Fisher Scientific) in one of the following imaging systems:
 - Odyssey DLx Imager (LI-COR),
 - Stella imaging system (Raytest), or
 - ChemiDoc MP imaging system (Bio-Rad).
- 15 Perform densitometric analysis of protein expression using ImageJ software.