



Jan 09, 2025

Western Blot - Frozen cell pellets

 [Cell Genomics](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1d562vr2/v1>

Anita Adami¹, Laura Castilla-Vallmanya¹

¹Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



Laura Castilla-Vallmanya

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1d562vr2/v1>

External link: <https://doi.org/10.1016/j.xgen.2025.100979>

Protocol Citation: Anita Adami, Laura Castilla-Vallmanya 2025. Western Blot - Frozen cell pellets. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1d562vr2/v1>

Manuscript citation:

Adami A, Garza R, Gerdes P, Johansson PA, Dorazehi F, Koutounidou S, Castilla-Vallmanya L, Atacho DA, Sharma Y, Johansson JG, Tam O, Kirkeby A, Barker RA, Hammell MG, Douse CH, Jakobsson J (2025) LINE-1 retrotransposons mediate cis-acting transcriptional control in human pluripotent stem cells and regulate early brain development. Cell Genomics 5(10). doi: [10.1016/j.xgen.2025.100979](https://doi.org/10.1016/j.xgen.2025.100979)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 09, 2025

Last Modified: January 09, 2025

Protocol Integer ID: 117993

Keywords: ASAPCRN, western blot for protein detection, protein from frozen pellet, frozen cell pellets this protocol, frozen cell pellet, frozen pellet, protein detection, western blot, protein

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-025170

Abstract

This protocol describes how to extract protein from frozen pellets and how to perform a western blot for protein detection.



Protein Extraction

1h

- 1 Thaw the cell pellets On ice . 10m
- 2 Add RIPA buffer supplemented with protease inhibitors, vortex well and lyse the cells on ice for 00:30:00 . 30m
(For a 200,000 cells pellet use 40 μ L of lysis buffer, scale up or down accordingly)
- 3 Centrifuge at 12.000 rpm, 4°C for 00:20:00 20m
- 4 Carefully transfer the supernatant to a new tube. Store the protein lysate at -80 °C or continue with the steps in the Western Blot section.

Western Blot

6h 25m

- 5 **Prepare gel samples:** 5m
 - Mix the appropriate amount of sample with 4x loading dye (LDS), 10x reducing agent and water.
 - Boil the samples for 00:05:00 at 95 °C .
 - Spin briefly.
- 6 **Run the gel:** 45m
 - Load the boiled samples in a pre-made 4-12% Tris-glycine SDS-PAGE gel.
 - Load a well with a prestained ladder.
 - Run at 200V in 1x MES SDS running buffet for 00:45:00 at Room temperature .
- 7 **Transfer to membrane:** 10m
 - Transfer the proteins to a PVDF membrane using the Transblot-Turbo Transfer system (BioRad) following manufacturer's instructions.
 - Check the ladder has transferred well to the membrane.
- 8 **Blocking:** 2h
 - Cut the membrane if the size of the proteins to be detected is different enough.



- Place the membrane in a cuvette making sure the protein side is facing up.
- Briefly rinse the membrane with Tris-buffer saline with 0,1% Tween (TBST).
- Block with 5% skimmed milk-TBST for 01:00:00 at Room temperature or Overnight at 4 °C .

9 Antibody incubation:

3h 25m

- Dilute primary antibody in 5% skimmed milk-TBST.
- Incubate on a shaker Overnight at 4 °C . For abundant targets, incubate for 01:00:00 at Room temperature .
- Briefly rinse the membrane with TBST.
- Wash three times with TBST on a shaker for 00:10:00 at Room temperature .
- Dilute species-appropriate HRP-conjugated secondary antibody in 5% skimmed milk-TBST.
- Incubate on a shaker for 01:00:00 at Room temperature .
- From this point, protect the membrane from light.
- Briefly rinse the membrane with TBST.
- Wash three times with TBST on a shaker for 00:10:00 at Room temperature .
- Wash once with TBS on a shaker for 00:05:00 at Room temperature .
- NOTE: If using an HRP-conjugated primary antibody, wash thrice and image (skip the secondary antibody step).

10 Detection:

- Prepare the detection reagent following manufacturer's instructions.
- Incubate on a shaker for following manufacturer's instructions.
- Image with chemiluminescence using an appropriate imager.

11 Stripping:

- Use a commercial stripping buffer as per instruction.
- Briefly rinse the membrane with TBST.
- Repeat blocking, antibody incubation and detection of the new target (usually the loading control).



NOTE: A large proportion of the protein will be lost during the stripping step. Always blot the most abundant target last.