



Apr 23, 2024

Western Blot in Mouse Brain Tissue for detecting pRab proteins

DOI

dx.doi.org/10.17504/protocols.io.261ge545jg47/v1

madalynn.erb Erb¹

¹Van Andel Research Institute



madalynn.erb Erb

Van Andel Research Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge545jg47/v1>

Protocol Citation: madalynn.erb Erb 2024. Western Blot in Mouse Brain Tissue for detecting pRab proteins. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ge545jg47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 23, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 98678

Keywords: ASAPCRN, prab proteins in mouse brain tissue, detecting prab protein, detecting prab proteins protocol, prab proteins protocol, western blot in mouse brain tissue, mouse brain tissue, western blot

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000592

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Protocol for detecting pRab proteins in mouse brain tissue



Materials

TBS (20L)

🧪 400 mL : 1M Tris pH 7.5

🧪 600 mL : 5M NaCl

🧪 19 L : H₂O

TBST (20L)

🧪 400 mL : 1M Tris pH 7.5

🧪 600 mL : 5M NaCl

🧪 19 L : H₂O

🧪 20 mL : Tween 20

10X Transfer Buffer

- 🧪 700 mL mH₂O
- 🧪 30.4 g Tris base
- 🧪 144 g Glycine -

Fill to 🧪 1000 mL with H₂O

Store at 🌡 Room temperature

1X Transfer Buffer:

- 🧪 100 mL 10X Transfer Buffer
- 🧪 200 mL Methanol
- 🧪 700 mL mH₂O

Store at 4°C

10X Running Buffer :

Dissolve the following components in 500 mL H₂O:

- 🧪 30.0 g Tris base
- 🧪 144.0 g glycine
- 🧪 10.0 g SDS (Sodium dodecyl sulfate)

Fill to 🧪 1000 mL with H₂O



The pH of the buffer should be 8.3 and no pH adjustment is required. Store the running buffer at room temperature and dilute to 1X before use.

5X Lammeli sample buffer

10% SDS

50% glycerol

25% 2-Mercaptoethanol


0.31M tris pH6.8

0.01% bromophenol blue

Troubleshooting










Lysate preparation

- 1 Freeze freshly dissected brain tissue on dry ice and store at  -80 °C
- 2 Make lysis buffer - keep on ice

	A	B	C	D
	Stock solutions	for 10mL	for 15mL	For 20mL
	Tris HCl 1M pH7.5	500uL	0.75 mL	1.0 mL
	EGTA 0.1M	100uL	0.15 mL	0.2 mL
	Triton X 100	100uL	0.15 mL	0.2 mL
	PIC I 100X	100uL	0.15 mL	0.2 mL
	PIC III 20X	500uL	0.75 mL	1.0 mL
	1.35M Sucrose (5X)	2.0mL	3.0 mL	4.0 mL
	Complete EDTA-free protease inhibitor cocktail (Sigma-Aldrich Cat #11836170001)	1 tablet	1 tablet	2 tablets
	water		8.55	
	Total			15.0



- 3 Remove tissue from  -80 °C and keep on dry ice before starting experiment
- 4 Weigh each tissue sample and record mass
- 5 Transfer each samples to clean 5mL test tube for homogenization
- 6 Place samples on ice and add appropriate amounts of lysis buffer based on sample mass
 - 6.1 For regional brain dissections add  1000 µL of chilled lysis buffer for every  100 mg of brain tissue
 - 6.2 For organ tissue add  1000 µL of chilled lysis buffer for every  100 mg of brain tissue
 - 6.3 For hemi brains or whole brains add  600 µL of chilled lysis buffer for every  100 mg of brain tissue
- 7 Allow samples to thaw in chilled lysis buffer for about 5 minutes
 - 7.1 Check samples with a clean pipette tip. If tissue has thawed, move onto the next step.
- 8 Use a tissue homogenizer to lyse the tissue
 - 8.1 10 seconds homogenization at maximum speed
 - 8.2 tissue should be kept on ice during and after homogenization



- 9 Allow tissue to rest on ice for 00:10:00 - 00:15:00 after homogenization 25m
- 9.1 Bubbles should disappear during this time
- 9.2 After 00:10:00 - 00:15:00 check lysate. If there are tissue chunks repeat homogenization step. 25m
- 10 Transfer the samples to Eppendorf tubes and centrifuge at 4 °C during 00:30:00 at 21.000 rcf 30m
- 11 Remove samples from centrifuge and keep on ice
- 12 Transfer supernatant to new Eppendorf tube
- 12.1 Transfer 15 µL of supernatant to new Eppendorf tube for Bicinchoninic acid (BCA) assay to measure protein concentration
- 12.2 Store supernatant at -80 °C until ready to run western blot
- 12.3 Save the pellet and store at -80 °C - pellet contains triton X insoluble proteins
- 13 Perform BCA analysis to measure protein concentrations

Western Blot


3h 55m

- 14 Dilute tissue lysates 1:5 in 5X Lammeli sample buffer
- 10% SDS
50% glycerol
25% 2-Mercaptoethanol
0.31M tris pH6.8



0.01% bromophenol blue


For  50 mL :

 7.75 mL 2M Tris pH6.8

 25 mL glycerol



 5 g SDS

 12.5 mL 2-Mercaptoethanol

Fill with milliQ H₂O to  50 mL (approx.  4.75 mL).

Add bromophenol blue

First mix Tris with 4ml H₂O and add SDS. Let it mix for about an hour. Then add BME and glycerol and continue to mix for approximately one more hour (SDS will eventually go into solution). QS to 50ml with additional H₂O. Add bromophenol blue, aliquot and store at -20C.


15 Denature the proteins at  70 °C for  00:10:00 in lysis buffer / Lammeli sample buffer mix

10m

16 Load  80 µg of protein into each well of SDS page gel

16.1 For Rab proteins use 12.5% acrylamide gels
For LRRK2 use 7.5% acrylamide gels

17 Run gels at 120V

17.1 For Rab proteins run gels  01:00:00

2h 45m


For LRRK2 run gels  01:45:00

18 Transfer gels at 20V  Overnight at  Room temperature

1h

Blot membranes

2h 35m

19 Wash membranes in H₂O for  00:02:00 then 5 min in TBST

2m

- 20 Wash membranes in TBST for ⌚ 00:05:00 5m
- 21 Incubate membrane in Ponceau S solution (Sigma P7170) for ⌚ 00:05:00 5m
- 22 Rinse membrane in H₂O and image
- 23 Wash 2 times in TBST for ⌚ 00:05:00 each wash 5m
- 24 Incubate membrane in 5% milk in TBST (BioRad 1706404XTU) for ⌚ 01:00:00 at 1h
🌡 Room temperature
- 25 Wash 2 times in TBST for ⌚ 00:02:00 each wash 2m
- 26 Incubate ⌚ Overnight at 🌡 4 °C with primary antibody diluted in 5% milk in TBST 1h

26.1

	A	B	C	D
	Target	Species	Manufacturer	Dilution
	pRab10	rabbit	Abcam ; ab230261	1:500
	Total Rab10	rabbit	Cell Signaling ; 8127S	1:1000
	pRab8a	rabbit	Abcam ; ab230260	1:500
	pRab12	rabbit	Abcam ; ab256487	1:500
	Total Rab12	rabbit	Protein Tech ; 18843-1-AP	1:1000
	GAPDH	mouse	Protein Tech ; 60004-1-Ig	1:5000
	Actin	mouse	Milipore; MAB1501	1:2000
	LRRK2	rabbit	Abcam ; ab133474	1:1000
	p935 LRRK2	rabbit	Abcam ; ab133450	1:1000
	p1292 LRRK2	rabbit	Abcam ; ab203181	1:1000



Primary antibodies for Rab proteins and LRRK2

27 Wash 3 times in TBST for  00:05:00

5m

28 Dilute HRP-conjugated secondary antibodies 1:10,000 in 5% milk in TBST


28.1

	A	B	C
	goat anti rat	112-035-175	Jackson ImmunoResearch
	goat anti mouse	115-035-174	Jackson ImmunoResearch
	mouse anti rabbit	211-032-171	Jackson ImmunoResearch
	goat anti guinea pig	ab97155	Abcam

Secondary Antibodies

29 Wash membrane 3 times in TBST for  00:05:00 and once in TBS

5m

30 Wash membrane in TBS for  00:05:00

5m

31 Image using ECL reagents (Amersham)

31.1 use  500 μ L of solution A and  500 μ L solution B for each membrane

31.2 Develop for at least  00:01:00 before imaging

1m

