

Sep 28, 2021

Western blot protocol for detecting ATP10B in mouse/rat brain

DOI

dx.doi.org/10.17504/protocols.io.byhfpt3n

María Sanchiz Calvo¹, Eduard Bentea¹, Veerle Baekelandt¹

¹Katholieke Universiteit Leuven



Eduard Bentea

KU Leuven

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.byhfpt3n

Protocol Citation: María Sanchiz Calvo, Eduard Bentea, Veerle Baekelandt 2021. Western blot protocol for detecting ATP10B in mouse/rat brain. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.byhfpt3n>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 24, 2021

Last Modified: May 31, 2024

Protocol Integer ID: 53511

Keywords: ATP10B, western blot, mouse, rat, brain, ASAPCRN

Abstract

Protocol for detection of ATP10B in rat and mouse brain tissue by Western blotting



Materials

Recipe for 100 mL 1x RIPA buffer :

- o 5 mL Tris pH 7.4
- o 1 mL Triton x100
- o 1 g Deoxycholate
- o 3 mL 5M NaCl
- o 200 µL 0.5M EDTA (pH 8.0)
- o 1 mL 10% SDS

Add water until 100 mL (MilliQ water)

Add protease and phosphatase inhibitors (Pi + PPI) right before use

Reagents used during Western blotting :

- o 6X SDS buffer
- o PageRuler™ Plus Prestained Protein Ladder
- o 3-8% NuPAGE™ Tris-Acetate gel
- o NuPAGE™ Tris-Acetate SDS Running Buffer
- o Trans-Blot® Turbo™ PVDF Membrane
- o 5% milk powder dissolved in PBS-T 0.1% (PBS + 0.1% Triton-X100)
- o ATP10B HPA034574 (Sigma) primary antibody
- o Goat anti-rabbit HRP secondary antibody (Dako)
- o ECL Prime chemiluminescence kit (GE Healthcare)

Recipe for 50 mL 6x SDS buffer (160 mM TrisHCl pH 6.8, 2% SDS, 200 mM DTT, 40% glycerol, bromophenol blue):

- o 5 mL 1.6M TrisHCl pH 6.8
- o 1 g SDS
- o 1.54 g DTT

Add water until 25 mL (MilliQ water) and dissolve

- o 20 mL glycerol

Add water until 50 mL (MilliQ water)

Add spatula tip of bromophenol blue



Before start

Perform protein extraction from snap-frozen brain tissue:

- weigh tissue
- add RIPA buffer (see Materials) : 10 X of the weight (40 mg = 400 μ L)
- homogenize samples using sample homogenizer
- sonicate samples at 4 degrees C, 3 times 15 seconds (keep the samples on ice between each sonication)
- centrifuge samples at 6000 g for 10 minutes at 4 degrees C
- collect supernatant and measure protein concentration
- aliquot protein extracts and store at -20 degrees C



Day 1

4h 5m

- 1 Prepare samples for Western blotting : 30m
 - 30 μg proteins in 12 μL of volume (adjust with milliQ water)
 - 2.4 μL of 6X SDS buffer + 10% of β -mercaptoethanol
- 2 Vortex samples and spin down
- 3 Boil the samples for 00:10:00 at 98 $^{\circ}\text{C}$ 10m
- 4 Load samples, together with 7 μL of mass marker (PageRuler™ Plus Prestained Protein Ladder) on a 3-8% NuPAGE™ Tris-Acetate gel. Use NuPAGE™ Tris-Acetate SDS Running Buffer for migration. 30m
- 5 Run the gel at 80V for 00:10:00 10m
- 6 Run the gel at 150V for another 00:45:00 45m
- 7 Transfer the proteins on PVDF membrane (Trans-Blot® Turbo™ PVDF Membrane) using Trans-Blot Turbo Transfer System (Bio-Rad), using the pre-programmed protocol STANDARD SD: 00:30:00 , up to 1.0 A, 25 V. 30m
- 8 Block membranes for 01:00:00 using 5% milk dissolved in PBS-T 0,1% at Room temperature 1h
- 9 Incubate membranes Overnight in ATP10B HPA034574 (Sigma) diluted 1/500 in 5% milk PBS-T 0,1% at 4 $^{\circ}\text{C}$ 30m

Day 2

3h 10m

- 10 Wash membranes with PBS-T 0,1% for 10 minutes, 5 times at Room temperature : 50m
 - PBS-T 00:10:00
 - PBS-T 00:10:00
 - PBS-T 00:10:00



- PBS-T 00:10:00

- PBS-T 00:10:00

11 Incubate membranes for 02:00:00 in secondary antibody solution Goat anti-rabbit HRP (Dako) diluted 1/10000 in 5% milk PBS-T 0,1% at Room temperature

2h

12 Wash membranes with PBS-T 0,1% for 10 minutes, 5 times at Room temperature :

50m

- PBS-T 00:10:00

- PBS-T 00:10:00

- PBS-T 00:10:00

- PBS-T 00:10:00

- PBS-T 00:10:00

13 Develop membranes using ECL Prime (GE Healthcare) for at least 00:10:00

10m