



Oct 14, 2019

Western Blot (semi-dry)

DOI

dx.doi.org/10.17504/protocols.io.764hrqw

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.764hrqw>

Protocol Citation: Igem Dusseldorf 2019. Western Blot (semi-dry). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.764hrqw>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 14, 2019



Last Modified: October 14, 2019

Protocol Integer ID: 28604

Keywords: equilibrate gel in cathode buffer, equilibrate gel, wash blot in tb, eliminating air bubble, membrane on top put gel, wash blot, minutes wash membrane with tbst, developer to membrane, possible to gel size, ml 10x anode buffer, ml per membrane, place on gel, blotting duration, minutes wash membrane, block membrane with tbst, expensive antibody, air bubble, cathode buffer, antibody, ad 1l h2o, secondary antibody, top put gel, gel size, membrane drain, blotting paper, run blot, place on anode put pcd, blot, anode buffer, bubble, anode put pcd, membrane activate for ca, gel, primary antibody, membrane, immunodetection, membrane activate

Abstract

- Semi-dry cathode buffer
 - 100 ml 10x conc. Cathode buffer
 - 200 ml Methanol
 - ad 1L H₂O
- Semi-dry anode buffer
 - 100 ml 10x anode buffer
 - 200 ml Methanol
 - ad 1L H₂O
- TBST
 - 100 mM Tris/HCl pH=8, 150 mM NaCl, 0,1% Tween
- TBST + 5% milk powder
- Primary / secondary antibodies

Western Blot

Equilibrate gel in Cathode buffer for 30 minutes (no shaking)

PCDF-Membrane

Activate for ca. 20 seconds in MeOH, incubate for 2 minutes in dH₂O (no shaking)

Equilibrate for 5 minutes in Anode buffer

Wet two Blotting papers (thickness 1 mm each) in Anode buffer and place on Anode

Put PCDF-membrane on top

Put Gel on top

Wet two blotting papers (thickness 1 mm each) in Cathode buffer and place on gel.

Cut membranes and papers as exactly as possible to gel size (Careful! Gel gains some size while equilibrating) and place each component exactly on top of each other. After each step consider eliminating air bubbles by rolling them out with a glass rod with buffer.

Run blot at 1,5 mA/cm² for 90 minutes (depending on the size of the protein, the blotting duration will have to be adjusted)

Immunodetection

Block membrane with TBST + 5% milk powder for 45 minutes

Wash membrane with TBST and primary antibody for at least 1h at RT or 4°C over night. Anti-gfp antibody will be diluted 1:5000 in TBST, use as little volume as possible to preserve the expensive antibody.

Wash for 10 minutes in TBST. Repeat 3 times.

Add secondary antibody and incubate at least 45 minutes at RT or at 4°C over night. Swivel gently.

Wash for 10 minutes in TBST. Repeat 3 times.



Mix Bio-Rad Clarity Western 1:1 (ca. 1 ml per membrane). Let membrane drain, evenly add developer to membrane (with 1ml tip). Spread between film (remove bubbles), incubate for 1 minutes and detect.

Blots; Chemi or Chemi-high Sensitivity, Signal Accumulation mode, transmission light for ladder detection

Don't forget to take a picture of membrane for merge!

Stripping

Wash blot in TBS

immerse in stripping buffer

incubate at RT for 5-15 minutes (Under the hood!)

remove and wash again in TBS

Troubleshooting

