

Aug 09, 2023

Western blotting using the BioRad Criterion Blotter system

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkowy1v5r/v1>

Louise Uoselis¹

¹Lazarou Lab, WEHI



Louise Uoselis

WEHI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkowy1v5r/v1>

Protocol Citation: Louise Uoselis 2023. Western blotting using the BioRad Criterion Blotter system. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlkowy1v5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 09, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 86225

Keywords: ASAPCRN, biorad criterion blotter system, western blot analysis, biorad criterion blotter system protocol, page western blot analysis, western blotting, sd

Abstract

Protocol for performing an SDS-PAGE Western blot analysis using the BioRad Criterion Blotter system.



Day 1

1h 47m

- 1 Load samples onto a 4-12% Bis-Tris (NuPage) gel with MOPS buffer (ThermoFisher). If using SDS loading dye (with 5% SDS final concentration), all lanes must contain the same volume of sample. Ensure all lanes are filled with 1x SDS loading dye at the same volumes as the sample/marker lanes. If this is not done, samples will spread according to the total sample volume, and the wells will not run evenly.
- 2 Run gels according to the following gradient: 100 V (⌚ 00:10:00), 150 V (⌚ 00:10:00), 200 V (until the dye front has reached the desired distance – typically ~40 min).
- 3 Cut 1x piece of PVDF (0.45 um), and 2x pieces of Whatman paper per gel, according to the gel dimensions.
- 4 Place the PVDF membrane into methanol for ~ 5 seconds, remove the methanol and wash the membrane three times in distilled water.
- 5 Place the membrane in cold (⌚ 4 °C) 20% methanol Bis-Tris transfer buffer (ThermoFisher)
- 6 Remove the gel from the cassette and place in ice cold 20% methanol Bis-Tris transfer buffer (ThermoFisher)
- 7 Add a frozen ice pack to the Criterion Blotter transfer tank, and fill the tank with cold transfer buffer.
- 8 Soaking each Whatman paper piece in transfer buffer for ~ 5 sec before use, assemble the following transfer stack: 1x Whatman, SDS-PAGE gel, PVDF membrane, 1x Whatman, rolling each layer to prevent bubbles accumulating between pieces
- 9 Place the transfer stack between sponges and secure them in the cassette in the correct orientation, and place the cassette in the transfer tank. If transferring 1 gel, ensure an empty cassette containing only sponges is placed in the second cassette holder in the tank.
- 10 Run the transfer at 100 V for ⌚ 01:00:00 .

20m

1h



- 11 Remove the membrane from the transfer stack and place into destain (50% methanol/7% acetic acid) for 00:02:00 rocking at Room temperature . 2m
- 12 Wash the membrane 3x in 0.05% v/v Tween/PBS for 00:02:00 rocking at Room temperature . 2m
- 13 Incubate the membrane in 3% w/v milk powder in 0.05% v/v Tween/PBS for 00:20:00 rocking at Room temperature . 20m
- 14 Wash the membrane 3x in 0.05% v/v Tween/PBS for 00:03:00 rocking at Room temperature . 3m
- 15 Place the membrane in the desired primary antibodies made up in 3% w/v BSA in 0.05% v/v Tween/PBS and incubate rocking, either overnight at 4 °C or for 01:00:00 at Room temperature . 1h

Day 2

1h 14m

- 16 Wash the membrane 3x in 0.05% v/v Tween/PBS for 00:03:00 rocking at Room temperature . 3m
- 17 Place the membrane in the desired secondary antibody in 3% w/v milk powder in 0.05% v/v Tween/PBS for 01:00:00 rocking at Room temperature . 1h
- 18 Wash the membrane 2x in 0.05% v/v Tween/PBS for 00:03:00 rocking at Room temperature . 3m
- 19 Wash the membrane once in PBS for 00:03:00 rocking at Room temperature . 3m
- 20 Incubate the membrane in ECL developing reagent for 00:05:00 , prior to developing the blot using a ChemiDoc Imaging System. 5m

