

Feb 04, 2020

Version 1

Whole Blood Processing (SepMate) V.1

DOI

dx.doi.org/10.17504/protocols.io.ba87ihzn



Yang Sun¹, David Lee¹, Jimmie Ye¹

¹UCSF

Human Cell Atlas Metho...



Alyssa Ward

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ba87ihzn>

Protocol Citation: Yang Sun, David Lee, Jimmie Ye 2020. Whole Blood Processing (SepMate). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ba87ihzn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 12, 2020

Last Modified: February 04, 2020

Protocol Integer ID: 31743

Keywords: sepmate tube, whole blood processing, whole blood, sepmate, pbmc, donor, lymphoprep

Abstract

Purpose: To purify PBMCs from 20 or 50mL whole blood (per donor). SepMate tubes will be used with Lymphoprep to ease processing time and effort.

Attachments



PBMC isolation from ...

17KB



PBMC isolation from ...

17KB



Materials

MATERIALS

✕ 50ml Conical Tubes, Green Cap, 25/Bag **Bio Basic Inc. Catalog #CT788-G.SIZE.1PK**

✕ DMSO **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1435**

✕ SepMate™-50 (IVD) 100 Tubes **STEMCELL Technologies Inc. Catalog #85450**

✕ Lymphoprep™ 250 mL **STEMCELL Technologies Inc. Catalog #7801**

✕ FBS **Invitrogen - Thermo Fisher**

✕ Protein LoBind 1.5mL microcentrifuge tubes **Eppendorf Catalog #0030108116**

✕ EasySep™ Buffer **STEMCELL Technologies Inc. Catalog #20144**

✕ Falcon™ 15mL Conical Centrifuge Tubes **Fisher Scientific Catalog #14-959-53A**

✕ Cryovial 2.0 ml round base internal thread screw cap writing area sterile **greiner bio-one Catalog #121277**

✕ DPBS, no calcium, no magnesium **Thermo Fisher Catalog #14190136**

✕ Mr. Frosty Freezing Container, 2mL tubes, Nalgene Mr. Frosty Freezing Container for 1-2mL cryogenic tubes, PC, clear w/ blue lid, 1/Cs. **Thermo Fisher Catalog #5100-0001**

- SepMate-50 (STEMCELL)
- Lymphoprep (STEMCELL)
- EasySep Buffer (STEMCELL)
- DPBS (Ca/Mg free; Fisher)
- 15mL conical tubes
- 50mL conical tubes
- 1.5mL low bind tubes
- ACK lysis buffer
- Freezing media A (100% FBS)
- Freezing media B (20% DMSO in FBS)
- 2.0mL Cryovials
- Freezing container

Troubleshooting

Safety warnings

⚠ Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.



Before start

Lymphoprep should be at 🌡 Room temperature (~ 🌡 20 °C) before use. It is important to keep it at 🌡 Room temperature for the density to remain 1.077g/ml.



- 1 Purify 20mL or 50mL whole blood using the following protocol:

STEP CASE

50mL whole blood

21 steps

To purify PBMCs from 50mL whole blood (per donor).

For each donor, 50mL blood will be collected in 5 BD heparin tubes (green top).

- 2 Centrifuge blood samples at  1200 x g, Room temperature, 00:10:00 swinging-out rotor and brake off within  06:00:00 of collection. 

- 3 Carefully transfer  3 mL plasma from each tube in a 15mL tube avoiding the buffy coat, use remaining blood for PBMC isolation. 

Note

Note: To minimize PBMC loss during plasma separation, leave the clear layer level at least 1cm above the RBC.

- 4 Centrifuge plasma at  16000 x g, 4°C, 00:10:00 in fixed-angle rotor (if a high-speed centrifuge for 15mL tubes is not available, aliquot plasma in multiple 1.5mL tubes and centrifuge at  16000 x g). Transfer plasma to new tubes avoiding the pellet. 

- 5 Aliquot  1 mL plasma in multiple 1.5mL low bind tubes, freeze at  -20 °C on the day of collection and move plasma samples to  -80 °C after  24:00:00 . 

PBMC isolation

- 6 Prepare SepMate tubes with Lymphoprep

- 6.1 Prepare 3 SepMate tubes per donor and add  15 mL Lymphoprep below the insert by pipetting into the center hole. 

**Note**

Take care to minimize any air bubbles below the plastic divider.

Video from SepMate manufacturers (~2min)

<https://www.youtube.com/embed/K9lIBxjgLt4>

7

Dilute whole blood 1:2 with DPBS (Ca/Mg free).

**Note**

For 105mL total:  35 mL residual blood and  70 mL DPBS .

8

Slowly pipette  35 mL diluted blood down the side of the tube.

**Note**

Some mixing may occur between the medium and blood, but take care not to mix under the divider.

Unlike conventional gradient separation, do not tilt the tube when adding diluted blood.

9

Spin at  1200 x g, 20°C, 00:15:00 with brakes ON.

**Note**

Prepare three destination 50mL tubes (~30mL per Sepmate tube will be yielded.)



- 10 When the spin finishes, pour off the top layer from the SepMate tubes into your prepared destination tubes.

Note

Do not invert the tubes for more than 2 seconds.

- 11 Spin at  450 x g, Room temperature, 00:10:00 with brake ON. 

Note

First wash at higher RCF because gradient medium is mixed with the cells.

- 12 Pour off supernatant and combine pellets in three tubes to one, add EasySep buffer to 50mL. 

- 13 Spin at  350 x g, Room temperature, 00:10:00 with brakes ON. 

Note

If yield is critical, ensure that the supernatant is clear; opaque supernatant indicates incomplete centrifugation of cells. If necessary, perform the spin at  450 x g .

- 14 Pour off supernatant, wash in  50 mL EasySep buffer for a second wash. 

- 15 Spin  350 x g, Room temperature, 00:10:00 with brakes on. 

- 16 Pour off supernatant and resuspend in  20 mL EasySep buffer . 



- 17 Perform cell count.
- 18 Spin at  350 x g, 00:10:00 and resuspend in appropriate volume of freezing media A to get 10e7 cells/mL. 
- 19 Prepare cryovials with  500 µL freezing media B . 
- 20 Add  500 µL cell suspension from previous step to cryovials and gently mix by pipetting. 
- 21 Place in freezing container for 1 day at  -80 °C and move to LN2 tank for long term storage on the following day.