



Jan 08, 2020

Version 1

🌐 Whole-body clearing of beetles by successive treatment of hydrogen peroxide and CUBIC reagents V.1

DOI

dx.doi.org/10.17504/protocols.io.ba33igqn



Monto Kuroda¹, Shinya Kuroda²

¹Keika High school, I, 5-6-6 Hakusan, Bunkyo-ku, 112-8612, Tokyo, Japan;

²Department of Biological Sciences, Graduate School of Science, University of Tokyo, 7-3-1 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan

Optical Clearing of Tissue



Shinya Kuroda

University of Tokyo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ba33igqn>

Protocol Citation: Monto Kuroda, Shinya Kuroda 2020. Whole-body clearing of beetles by successive treatment of hydrogen peroxide and CUBIC reagents. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ba33igqn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our group and it is working.

Created: January 08, 2020

Last Modified: January 08, 2020

Protocol Integer ID: 31579

Keywords: Whole-body clearing, hydrogen peroxide, CUBIC, Trypoxylus dichotomus, Rhomborrhina japonica, Dorcus titanus typhon, body clearing method for large beetle, body clearing of beetle, body clearing of the large beetle, body clearing method, body clearing, mammalian method because of pigment, cubic reagents internal tissues of multicellular organism, combined method of hydrogen peroxide, insects such as beetle, hydrogen peroxide treatment, mammalian method, successive treatment of hydrogen peroxide, beetle, large beetle, body imaging cocktail, internal tissue, insect, advanced clear, multicellular organism, unobstructed brain, pigment, exoskeleton, tissue, such as melanin

Abstract

Internal tissues of multicellular organisms cannot directly be seen because they contain pigments. For this reason, whole-body clearing methods have been developed and applied to mammals such as mice. Insects such as beetles, however, cannot be cleared by the mammalian method because of pigments such as melanin in their exoskeletons. In this study, we tried to develop a whole-body clearing method for large beetles. We first bleached the exoskeleton using a hydrogen peroxide treatment, and applied the the advanced Clear, Unobstructed Brain/Body Imaging Cocktails and Computational analysis (CUBIC) reagents to make the internal tissues transparent. The combined method of hydrogen peroxide and the advanced CUBIC allowed us to successfully perform whole-body clearing of the large beetles.

Attachments



[Monto Kuroda.pdf](#)

408KB



Materials

10% Formaldehyde Neutral Buffer Solution (Nacalai Tesque, cat. no. 37152-51)

4% Paraformaldehyde Phosphate Buffer Solution (FUJIFILM Wako Pure Chemical Co. cat. no. 163-20145)

Ethanol (FUJIFILM Wako Pure Chemical Co. cat. no. 057-00451)

Hydrogen Peroxide (H₂O₂) (FUJIFILM Wako Pure Chemical Co. cat. no. 081-04215)

Urea (FUJIFILM Wako Pure Chemical Co. cat. no. 219-00175)

N,N,N',N'-tetrakis (2-hydroxypropyl) ethylenediamine (Quadrol) (Tokyo Chemical Industry, cat. no. T0781)

Triton X-100 (Sigma-Aldrich, cat. no. X100-500ML)

Nitrilotriethanol (FUJIFILM Wako Pure Chemical Co. cat. no. 145-05605)

Sucrose (FUJIFILM Wako Pure Chemical Co. cat. no. 196-00015)

Troubleshooting

